CLEANING AND STERILISATION MANUAL

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

Plus / Flex

Copyright

© Mectron S.p.A. 2024. All rights reserved. No part of this document can be reproduced in any form without the written consent of the copyright owner.
Images are for demonstration purposes only.

TABLE OF CONTENTS

1	Introduction	1
2	Disassembly of Parts	2
3	Cleaning	6
3.1	Cleaning the Device Body and the Foot Pedal	6
3.2	Cleaning the PIEZOSURGERY plus/flex Components	8
3.2.1	Pre-Cleaning	8
3.2.2	Manual Cleaning	8
3.2.3	Automatic Cleaning	14
4	Cleaning Check	18
5	Drying and Lubrication of Components	18
6	Sterilisation	20

EN

PAGE INTENTIONALLY LEFT BLANK

1 INTRODUCTION

Read this manual carefully before proceeding to install, use, maintain or carry out any other interventions on the device. Make sure it is available at all times in the reprocessing rooms of the device and its components.

IMPORTANT: to avoid any injury to people or damage to property, please pay careful attention to all the safety precautions included in this manual. Depending on the level of risk involved, the safety requirements are classified as follows:

 **WARNING:** always related to personal injury.

 **CAUTION:** relating to possible damage to objects.

A **NOTE** identifies special information that may allow for simpler system maintenance or clarify or emphasise important instructions.

The main aim of this manual is to protect patients, surgeons and hospital personnel from possible infections caused by contaminated instruments. It is therefore essential to pay careful attention when performing each step of the reprocessing procedure. These instructions provide detailed information to allow hospital personnel to safely manage all operations and carry out effective cleaning of the:

- PIEZOSURGERY plus/flex device;
- Foot pedal;

and effective cleaning and sterilisation of:

- Handpiece;
- Torque wrench.

The hospital or site that physically carries out the procedures is responsible for cleaning and sterilising the instruments.

The procedures and parameters described in this manual have been validated by laboratory tests carried out by an independent entity. These tests have confirmed that the sterility level (SAL) of 10^{-6} is only achieved if the instructions described in this manual are followed perfectly.

All stages of sterilisation must be performed by the operator in compliance with the current revision standards: UNI EN ISO 17665-1, UNI EN ISO 556-1 and ANSI/AAMI ST.46.

 **CAUTION:** The instructions reported in the following are intended for qualified personnel able to treat medical devices for their re-use. These instructions have been validated by Mectron. The site (the hospital or facility that physically performs the cleaning and sterilisation procedures of instruments) remains responsible for ensuring that reprocessing is carried out in such a way as to achieve the desired results, by using the instructions, the equipment and the materials described here, and qualified personnel. The process must be regularly monitored and validated. Any deviation applied to the reprocessing procedures with respect to that specified in these instructions must be validated and the effectiveness and potential adverse consequences assessed.

NOTE: The handpieces and the torque wrenches are supplied non-sterile and must be cleaned and sterilised at first use and at each re-use. For the handling of reusable items, always refer to the relevant user manuals, which are supplied with the item.

 **WARNING:** The procedures described in this manual do not apply to sterile, single-use components. These components cannot be reused. Single-use parts are supplied in sterile packaging and are not intended to be re-used. The "cleaning and sterilisation manual" must be used together with the "use and maintenance manual" supplied with the PIEZOSURGERY plus/flex unit.

LIMITATIONS TO REPEATED CLEANING AND STERILISATION CYCLES:

The products are designed to tolerate a high number of sterilisation cycles. The materials used have been selected specifically for this purpose. Regardless, each time the products are reprocessed for use, the thermal and mechanical stresses to which they are subjected causes the materials to age. Always replace products showing signs of wear or premature deterioration, regardless of the number of sterilisation cycles to which they can still be theoretically subjected.

⚠ WARNING: Cleaning and sterilisation of new or repaired tools. All new and repaired device components are supplied in non-sterile conditions. At first use, and after each treatment, they must be cleaned and sterilised in strict compliance with the instructions given in this manual.

⚠ WARNING: Infection control. Special care must be taken in handling contaminated inserts as they are a cause of a potential danger of infection. Hospital personnel must wear suitable personal protection equipment.

⚠ WARNING: Handle sharp and pointed instruments carefully.

2 DISASSEMBLY OF PARTS


Before proceeding with the cleaning procedures described in *Chapter 3 on page 6*, disconnect all the device components.

» PROCEDURE

⚠ WARNING: Disconnect the device from the power supply using the general switch on the device console and disconnect it from the electricity supply. Disconnect all leads from the device (handpiece cord and foot pedal).

⚠ WARNING: Supplementary equipotential cable. If present, disconnect the supplementary equipotential cable before cleaning or sterilising.

1


Disconnect the device foot pedal: press the foot pedal connector release sleeve and pull the connector backwards;

⚠ CAUTION: Do not attempt to unscrew or rotate the connector when disconnecting: the connector may become damaged.

⚠ CAUTION: When disconnecting the foot pedal cable, always hold the cable connector only. Never pull the cable itself.

2


Disconnect the handpiece from the device;

⚠ CAUTION: Grip the connector only when carrying out this operation. Do not carry out this manoeuvre using the cord. Do not attempt to unscrew the connector when disconnecting the handpiece. The connector may become damaged.

⚠ CAUTION: The handpiece and its cord cannot be separated.

3


Remove the irrigation kit as follows:

- Remove the spike from the irrigation bag;


NOTE: Close the irrigation kit clamps before disconnecting it from the irrigation bag.

4


- Open the peristaltic pump cover completely and remove the tube;
- Close the peristaltic pump cover;

5


- Remove the tube and the clips from the handpiece cord;

6

- Disconnect the tube from the nozzle on the back of the handpiece;
- Dispose of the irrigation kit in accordance with the current regulations governing hospital waste;


7

Remove the insert from the handpiece using the torque wrench supplied with the device;

Dispose of the single use insert in accordance with current regulations governing hospital waste. For the handling of reusable items, always refer to the relevant user manuals, which are supplied with the item;

8


⚠ CAUTION: Only use the Mectron torque wrench to tighten and remove the insert on/from the handpiece. Do not use any other tools, such as pliers, scissors, etc.

⚠ WARNING: Sterile single-use.
Sterile single-use objects must only be used in one surgical procedure and on one patient. Single-use objects must not be reused. Separate and dispose of each single-use object in accordance with current regulations governing hospital waste

Unscrew the insert cover or metallic front cone of the handpiece by twisting in an anticlockwise direction;
Dispose of the insert cover in accordance with current national regulations governing hospital waste;

⚠ CAUTION: The metallic cone should not be disposed of; rather it should be reprocessed with the other reusable components.


9

FOR INSERTS WITH EXTENSION

Remove the extension from the handpiece using the torque wrench supplied;
Dispose of the extension in accordance with current national regulations governing hospital waste;

10


Connect the protection cap to the end part of the cord, at the connector cable gland;
Insert the cap into the handpiece cord connector.

⚠ CAUTION: Always install the protection cap before cleaning.

⚠ CAUTION: Before starting cleaning, always visually inspect the protection cap for any damage. Replace if any anomalies are seen.

11

3 CLEANING

Device console and foot pedal:

The only cleaning method applicable is the manual method

Components

(Handpiece and torque wrench):

The operator needs to choose from the following cleaning procedures:

- Manual cleaning;
- Automatic cleaning.

The procedures applicable to the PIEZOSURGERY plus/flex components described in the following manual have been validated by an independent body.

3.1 Cleaning the Device Body and the Foot Pedal

» PREPARATION

Check that all the following components have been removed/disconnected from the device (see *Chapter 2 on page 2*):

- electric power supply lead;
- irrigation kit;
- handpiece;
- irrigation bag support rod;
- supplementary equipotential cable connector (if present);
- foot pedal.

 **WARNING:** Always turn the device off by means of the switch and disconnect it from the electric grid before cleaning or disinfecting.

 **WARNING:** The device console is not protected against the penetration of liquids. Do not spray liquids directly onto the surface of the device.

 **CAUTION:** Do not immerse the device console or foot pedal in liquids and/or solutions of other types.

 **CAUTION:** Do not sterilise the device console or the foot pedal.

 **CAUTION:** Do not use machinery/ devices to clean the device console or foot pedal.

 **CAUTION:** Do not rinse the electric contacts of the foot pedal under running water.

» REQUIRED MATERIALS

- Clean, soft lint-free cloths;
- Detergent (pH 6-9) and, if necessary, a non-aggressive neutral pH (pH 7) disinfectant solution.

 **CAUTION:** For disinfection, if required, it is recommended to use water-based disinfectant solutions, with neutral pH (pH7). Alcohol-based disinfectant solutions and hydrogen peroxide are contraindicated as they can discolour and/or damage plastic materials. This also applies to chemicals such as acetone and alcohol.

» **PROCEDURE**


Clean the device and foot pedal surfaces with a clean, soft lint-free cloth, dampened with a detergent solution (pH 6-9), following the manufacturer's instructions;

For disinfection, if required, it is recommended to use a non-aggressive neutral pH (pH7) solution prepared according to the manufacturer's instructions;


Dry the device body using a clean, non-abrasive, soft lint-free cloth;

Dry the foot pedal using a clean, non-abrasive, soft lint-free cloth.


1

2

3

3.2 Cleaning the PIEZOSURGERY plus/flex Components

The Operator can choose between the following cleaning procedures:

- Manual cleaning;
- Automatic cleaning.

The procedures applicable to the PIEZOSURGERY plus/flex components have been validated by an independent body.

⚠ WARNING: The automated cleaning must not be used as the sole cleaning method for surgical instruments. The automated cleaning cycle must always be preceded by a manual cleaning procedure.

3.2.1 Pre-Cleaning

The "pre-cleaning" cycle is optional and preparatory to the chosen cleaning method (manual or automatic) and is recommended in the following cases:

- if the components to be cleaned are particularly contaminated and/or dirty;
- if a considerable time has passed since last use and pre-cleaning.

The person in charge of reprocessing activities will assess the need for a possible pre-cleaning cycle according to the conditions of the components, tools and specific requirements. For directions concerning the pre-cleaning of the various components, follow the manual cleaning instructions (*Chapter 3.2.2 on page 8*).

3.2.2 Manual Cleaning

» PREPARATION

The process applies to the following components:

- Handpiece;
- Metallic front cone;
- Torque wrench.

After having disconnected the components, following the instructions provided in *Chapter 2 on page 2* proceed with the cleaning procedure described in the following.

⚠ CAUTION: Before cleaning the handpiece, remove the insert.

⚠ CAUTION: Always install the protection cap before cleaning.

⚠ CAUTION: Before proceeding with the cleaning, always visually inspect the protection cap for any damage. Replace it if any anomalies are detected.

⚠ CAUTION: Do not immerse the handpiece in ultrasonic tank.

» MATERIALS REQUIRED

- Enzymatic detergent (pH 6-9);
- Tap water;
- Demineralized water;
- Suitably-sized container;
- Clean soft bristled, nylon brush;
- Syringe;
- Ultrasonic tank;
- Compressed air.

⚠ WARNING: Run the procedure immediately at the end of the intervention. Do not allow contaminated instruments to dry before reprocessing.

⚠ CAUTION: Cleaning instruments.

Use appropriate cleaning instruments as specified in this manual. Never use brushes, metallic sponges, abrasive instruments or sharp objects to remove residues during manual cleaning. These may ruin the product surface and permanently and definitively damage the instruments. Only use brushes with soft, nylon bristles.

⚠ CAUTION: If you intend to disinfect the components, it is recommended to use water-based disinfectants with neutral pH (pH7). Alcohol-based disinfectant solutions and hydrogen peroxide are contraindicated, because they can fade the colour and/or damage the plastic materials. This is also valid for chemical products such as acetone and alcohol. Always rinse with sterile water to maintain disinfection.

⚠ WARNING: Ensure the continued cleanliness and efficiency of the ultrasonic tank before use.

⚠ CAUTION: When cleaning, do not use organic solvents such as acetone or isopropyl alcohol to clean the instruments. Saline cleaning agents containing aldehydes, mercury, active chlorine, chloride, bromine or bromide, iodine or iodide or fixing or corrosive proteins must never be used. Alcohol based disinfectant and hydrogen peroxide are contraindicated as they can discolour and/or damage plastic materials. Do not place or soak instruments in a rinsing solution. Do not use solvents, lubricants or other products with chemical characteristics that differ from those specified. Use of these products can cause the instruments to malfunction or the release of undesired substances at the operating site.

» PROCEDURE

Prepare an enzymatic detergent solution with pH 6-9, according to the manufacturer's instructions^{a)}.


⚠ CAUTION: Once used, dispose of the enzymatic detergent correctly, do not recycle it.

a) Process validated by independent body with enzymatic detergent Enzymec.


Place the handpiece, the metallic front cone and the torque wrench, in a container. Add enough enzymatic solution to completely cover the components.

⚠ CAUTION: Do not immerse the handpiece if warm.

1


2

Leave the handpiece, the metallic front cone and torque wrench soaking in the enzymatic detergent for at least 10 minutes.

This will reduce the organic residues on the instruments.


3

For the handpiece only:


Using a disposable syringe (min. 20 ml, without needle), inject the enzymatic solution at least 3 times into the cavities and ducting of the handpiece to effectively remove any residue from the surfaces of the internal ducts.

4

Gently clean the surface of the handpiece the metallic front cone and torque wrench using the enzymatic detergent solution and a soft-nylon-bristle brush; brush the cavities with a soft-nylon-bristle pipe cleaner to remove all traces of visible soil in the internal and external parts. Brush each component for at least 2 minutes. Pay particular attention to the following areas:

5

- handpiece thread;
- stem and terminal part;
- external and internal parts of metal front cone.

Rinse the components under running water and brush the device with a clean, soft-nylon-bristle brush for at least 1 minute. Use a clean brush to brush the internal cavities and grooves under running water. Also rinse the handpiece cord.

⚠ WARNING: The following parts must be brushed carefully:

Handpiece

- Threaded part onto which the inserts are screwed;
- Threaded part onto which the metallic front cone is screwed;
- Parts adjacent to threaded parts.

Metallic front cone

- The internal threaded part.




Torque wrench

- External metallic ring;
- Internal channels and cavities.

NOTE: Rinse the cavities, grooves and fissures and critical parts thoroughly.

⚠ CAUTION: To prevent the accumulation of enzymatic detergent residues, rinse the instrument surfaces thoroughly.

6

For the torque wrench only:




Prepare a fresh solution of enzymatic detergent with pH 6-9, according to the manufacturer's instructions (use preferably warm water).

Place the torque wrench in a stainless steel instrument tray basket;
Place the basket in the ultrasonic tank;
Fully submerge the wrench in the enzymatic detergent solution and sonicate for at least 10 minutes.


Rinse the torque wrench, under running water and brush the device with a clean, soft-nylon-bristle brush for at least 1 minute. Use a clean brush to brush the internal cavities and grooves under running water.

Rinse/submerge the handpiece, the torque wrench and the metal cone, with demineralised water for at least 1 minute.

7

8

For the handpiece only:


Rinse the internal duct, injecting it 3 times with demineralised water using a 20 ml syringe (without needle) to remove any enzymatic detergent residue.

NOTE: Rinse the channel within the handpiece by injecting the liquid from the back, taking care to use a waste water container to avoid sprays. Leave the water to drain through the internal channel. Repeat the procedure until the rinsing water is free from residues.




Blow the internal ducts and all surfaces of handpiece, wrench and cone with compressed air until there is no more visible water residue.

Pay more attention to the inner channel by blowing from the back of the handpiece

9

10

EN

3.2.3 Automatic Cleaning

» PREPARATION

The process applies to the following components:

- Handpiece;
- Metallic front cone;
- Torque wrench.

Before proceeding with automated cleaning a manual pre-cleaning procedure of the disassembled components has to be performed (see *Chapter 3.2.2 on page 8*). Use washer-disinfectors compliant with ISO 15883 with baskets able to hold small-sized

products and with rinse connectors.

Properly place the product in a suitable basket, connect the lumens to the rinse connectors and start automatic cleaning.

⚠ CAUTION: Before cleaning in the washer-disinfector, evaluate the soil level and, if necessary, to avoid repeating the entire cleaning cycle, carry out a pre-wash according to the instructions in *Chapter 3.2.1 on page 8*.

» REQUIRED MATERIALS

- Alkaline detergent ¹⁾;
- Neutral detergent ²⁾
- Neutralising liquid ³⁾;
- Washer-disinfector ⁴⁾;
- Metal basket;
- Adapters for washer-disinfectors;
- Brush with soft nylon bristles.

1) For example, Neodisher® FA (0.2 % v/v).

2) For example, Renu Klenz® (Steris®).

3) For example, Neodisher® Z (0.1 % v/v).

4) For example, Miele or AT-OS washer-disinfector.

NOTE: For best and long-lasting results, Mectron recommends using automatic washer-disinfectors compliant with ISO 15883.


Always follow the instructions provided by the manufacturer of the washer-disinfector. Carefully follow the instructions provided by the producer of the disinfectant solution.

» **PROCEDURE**

Place the handpiece, the metallic front cone and torque wrench in a metallic tray.

⚠ CAUTION: Do not kink the handpiece cord; this may damage it.

1


2

Connect the specific adapter (supplied as optional) to the back of the handpiece and then to the appropriate water-jet connections of the washer/disinfector machine.

NOTE: Make sure that the components to be cleaned are appropriately blocked in the basket and cannot move during washing. Any blows could damage them. Position the instruments in a way that the water can flow through all the surfaces, even internal.

⚠ WARNING: Avoid overloading the thermal disinfector as this could compromise cleaning effectiveness.

It is possible to choose between two different cleaning cycles:

- neutral pH detergent cycle
- alkaline detergent cycle.

NEUTRAL PH DETERGENT CYCLE


Sequence and parameters applicable to the cycle:

- 3 mins. Pre-wash with cold water;
- 5 mins. Washing with 0,5% neutral pH detergent at 45°C (113°F);
- 2 mins. Rinsing with tap water;
- 2 mins. Rinsing with cold demineralised water;
- 5 mins. Thermodisinfection at 90°C-93°C (194°F-199.4°F) with demineralised water

Thermal disinfection has not been tested experimentally. In compliance with ISO 15883-1 thermal disinfection at a temperature of 90°C (194 °F) for 5 minutes determines an A0 value of 3000.

⚠ WARNING: The automatic cleaning cycle, as described above, must always be preceded by the manual cleaning steps.

3

ALKALINE DETERGENT CYCLE

Sequence and parameters applicable to the cycle:


- 1 min. Pre-wash with cold water;
- 5 mins. Washing with alkaline detergent at $55^{\circ}\text{C} \pm 2^{\circ}\text{C}$ ($131^{\circ}\text{F} \pm 3.6^{\circ}\text{F}$);
- 1 min. Neutralisation with suitable solution (1/3 cold e 2/3 warm);
- 1 min. Rinsing with tap water (1/3 cold e 2/3 warm);
- 5 mins. Thermodisinfection at 93°C (199.4°F) with demineralized water.

Thermal disinfection has not been tested experimentally. In compliance with ISO 15883-1 thermal disinfection at a temperature of 90°C (194°F) for 5 minutes determines an A0 value of 3000.

⚠ WARNING: The automatic cleaning cycle, as described above, must always be preceded by the manual cleaning steps.

⚠ WARNING: On completion of the cleaning cycle in the washer-disinfector, the handpiece remains at the wash temperature for a long time. Use appropriate precautions when extracting the handpiece from the washer-disinfector to prevent injury to the operator.

⚠ CAUTION: Due to its form the handpiece can roll. When not used the handpiece must always be put back on its support stand.

EN


4 CLEANING CHECK

All instruments must be inspected before they are prepared for sterilisation. A visual inspection must be performed by means of a magnifier under a suitable light source.

All instrument surfaces must be controlled to reveal the presence of soils, corrosion and wear. Special attention must be paid to the control of critical points such as lumens, threaded parts, etc.

If soils should be seen, repeat the chosen cleaning cycle.


1

5 DRYING AND LUBRICATION OF COMPONENTS


Dry the handpiece, metallic front cone and torque wrench with a clean lint-free cloth;

1

Use filtered medical compressed air to effectively dry the handpiece, metallic front cone and torque wrench;

2

⚠ WARNING: Dry the following thoroughly:

Handpiece and metallic cone:

- Internal channel;

Torque wrench:

- Cavities, grooves and fissures and critical parts.



Lubricate the wrench using a medical lubricant for surgical instruments. Spray the lubricant directly on the internal peripheral contact surfaces of the torque wrench.

After the lubrication remove the excess of lubricant with a cloth.

⚠ WARNING: Do not use oil- or silicone-based lubricants as they prevent direct contact of the surface with steam.

⚠ CAUTION: It is important to lubricate the wrench and gearing after each cleaning cycle, only using a medical lubricant for surgical instruments in order to avoid any malfunction.

3

EN

6 STERILISATION

» PREPARATION

The following products require sterilisation:

- Handpiece with cord;
- Metallic front cone;
- Torque wrench.

The packaging must ensure the sterility of the instruments until opened for use at the sterile field.

The instruments must be packed with a double medical grade steam sterilization pouch, in compliance with the ISO 11607 requirements.

» MATERIALS REQUIRED

- Single-use pouches for steam sterilisation;
- Steam steriliser.

 **CAUTION:** Sterilise only using a steam autoclave. Do not use any other sterilisation method as this may be incompatible with the materials used.

 **CAUTION:** During the sterilisation cycle, do not exceed the allowed load of the steam sterilizer.

 **CAUTION:** Ensure that the inner pouch is large enough to contain the instruments without stressing the seals or tearing the packaging.

NOTE: The handpiece must be packed with the handpiece protection cap open.

NOTE: Each instrument must be packaged individually.

» **PROCEDURE**

The sterilisation process must be carried out in a steam autoclave with a pre-vacuum cycle. Mectron guarantees SAL 10^{-6} , by setting the parameters according to only one of the two validated procedures listed below.

Procedure A:

- **Cycle type:** 3 times pre-vacuum
- **Minimum Sterilisation temperature:**
132°C (270°F) (0°C ÷ +3°C).
(range 132°C ÷ 135°C (270°F ÷ 275°F))
- **Sterilisation time:** 4 mins.
- **Drying time:** 10 mins.

or **Procedure B**


- **Cycle type:** 3 times pre-vacuum
- **Minimum Sterilisation temperature:**
134°C (273°F) (-1°C ÷ +1°C).
(range 133°C ÷ 135°C (270°F ÷ 275°F))
- **Sterilisation time:** values equal to or between 3 mins and 18 mins.
- **Drying time:** 10 mins.

All stages of sterilisation must be performed by the operator in compliance with the current revision standards: UNI EN ISO 17665-1, UNI EN ISO 556-1 and ANSI/AAMI ST:46.

⚠ CAUTION: The handpiece with cord is made of materials that withstand a maximum temperature of 135 °C (275 °F) for up to 20 minutes.

⚠ CAUTION: Do not sterilise the device console or pedal.

1



EN

Sterile instruments must be stored in a clean, dry place at a storage temperature of between 5°C and 40°C.

⚠ CAUTION:

- The handpiece and its cord cannot be separated.
- Do not sterilise the handpiece with the insert screwed into it.
- After sterilisation, the handpiece can only be used when it has cooled down to room temperature. The cooling process must not be accelerated.
- Re-screw the front cone to the handpiece before use.

2

» FOR US MARKET

In order to organize, transport, store and protect the devices, it is suggested to use protective sterilization cassettes/trays that are FDA-cleared, validated and legally marketed in the U.S.

Always follow the cassettes/trays manufacturer's instructions for proper use and check the compatibility with Mectron's sterilization parameters.


⚠ **CAUTION:** The cassettes/trays are intended as a protection system and must be used in conjunction with a legally marketed, FDA cleared sterilization wrap in order to maintain sterility of content.

The shelf life is dependent on the sterile barrier employed, storage manner, environmental conditions and handling.

A maximum shelf life for sterilized instruments should be defined by each health care facility.

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© Mectron S.p.A. 2024. Tutti i diritti riservati. Nessuna parte di questo documento può essere riprodotta, in qualsiasi forma, senza il consenso scritto del proprietario del copyright.
MECTRON S.p.A.
Piezosurgery, Inc.