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English

INSTRUCTION FOR USE

This device consists of two parts: a tip holder and a tip; and it has been designed for removing plaque or biofilm from the surfaces of implants, metal and ceramic dentures and natural teeth.

- IC1 (Fig. 1 - A) Tip (operational); can be used on ultrasound oral surgery devices produced by Mectron through the use of dedicated tip holders;
- ICP (Fig. 1 - B) Tip holder for ultrasonic oral surgery devices produced by Mectron.

⚠ WARNING: This device must be used in accordance with the correct settings described in Table 1.

The device must be used exclusively by specialized and duly trained personnel such as the Surgeon / Dentist, adult of any weight, age, height, gender and nationality, non-disabled.

NOTE: Regarding the intended patient group, patient selection criteria, indications for use, contraindications, and warnings, please refer to the use and maintenance manual supplied with your oral surgery device manufactured by Mectron.

IMPORTANT: The terms **"WARNING"** and **"CAUTION"** and **"NOTE"** identify conditions that require particular caution and must, therefore, be examined carefully. In particular:

⚠ WARNING: Identifies conditions and practices that present a risk of death or serious injury for the patient and/or user.

⚠ CAUTION: Identifies conditions and practices that may cause minor injuries and/or damage to the devices.

NOTE: Warning not associated with a hazardous or imminently dangerous condition.

IMPORTANT: Information regarding this user guide. The information contained in this document is limited to that strictly necessary for the first use of the product and must be consulted as a supplement to the "USE & MAINTENANCE MANUAL" supplied with the device. Read the "USE & MAINTENANCE MANUAL" supplied with the device carefully, paying particular attention to the "Safety Provisions" chapter, before performing any type of operation on the system and commencing any type of treatment.

⚠ WARNING: Before first use and subsequent uses. The tip holder and the tip and the torque wrench are supplied "non-sterile", therefore, before use, they must be reconditioned in accordance with the procedures described below and reported in detail in the "Cleaning and Sterilisation" section of the manual supplied with the device. After each use, the product must be cleaned and sterilised.

⚠ WARNING: If a serious accident attributable to the device occurs during correct and intended use, it is recommended to report it to the Competent Authority and to the company indicated on the product label.

1- Cleaning procedures

IMPORTANT: The cleaning procedures must be performed immediately after each use. Submerge the tip and/or instrument in demineralised water in an enzymatic cleaning solution immediately after use. Do not leave residue or blood deposits on the tips and instruments, eliminate larger impurities with a disposable cloth or paper towel.

The cleaning procedure can be performed using two different methods: manual cleaning and automatic cleaning. These two methods are alternative and only one or the other must be performed.

Limitations on repeated processes: Repeated processing has a minimal effect on these instruments. The end of their useful life is generally determined by wear and damage due to use.

⚠ CAUTION: Do not use hydrogen peroxide to clean the tips.

⚠ WARNING: To avoid injury to operators, treat the tips separately.

Particular attention must be paid while handling sharp and potentially infected instruments to avoid the risk of wounds and infections. Health personnel in charge of use and cleaning of the devices must use suitable protection and always use care when handling instruments with sharp edges and pointed tips.

1.1 - Manual Cleaning

NOTE: The tip holder, tip, and torque wrench must be processed separately.

- Do not leave the tip holder inserted into the torque wrench.
- Do not leave the tip installed on the tip holder.
- Do not place the tip holder, tip or torque wrench in the same container.
- Do not use the same toothbrush for the tip holder, tip, and torque wrench.

NOTE: Process independently validated with Enzymec enzymatic detergent and 0.8% dilution.

- Prepare a pH-neutral (6-9) enzymatic detergent solution, following the manufacturer's instructions.
- ⚠ CAUTION** Once used, dispose of the enzymatic cleaner solution properly, do not recycle.

- Place the tip holder, tip, or torque wrench in a clean container, in a horizontal position. Add enough of the prepared solution to completely cover the device to be cleaned

- Soak the device in the enzyme detergent solution for at least 10 minutes at room temperature. This process reduces the amount of blood, protein, and mucus on the instrument.

- After 10 minutes of immersion in the enzyme solution, gently brush all surfaces until any visible dirt is removed.

For the tip holder and tip: use a clean toothbrush with soft nylon bristles for the outer surfaces, a clean brush with soft nylon bristles for the inner cavities and crevices. Thoroughly clean hard-to-clean areas such as sharp edges.

For the torque wrench: Use a clean soft-bristled nylon toothbrush. Brush thoroughly, for about 20 seconds, all of the following parts of the torque wrench:

- through holes and internal channels;
- external metal ring;
- internal cavities, grooves and fissures.

- Remove the device being cleaned from the enzymatic detergent solution.

- Rinse thoroughly and brush all surfaces of the device being cleaned (see point 4) under running water for at least:

- 10 minutes for the wrench,
- 1 minute for the tip holder and the tip.

- For the tip holder and the tip add this operation: use a 20 ml disposable syringe to aspirate and inject the enzymatic detergent solution into hard-to-reach areas (through holes/cannulas). Repeat this operation three times to ensure effective removal of dirt from the internal surfaces of the through hole and inject fresh enzymatic detergent solution (see point 1).

- Place the device to be cleaned in a stainless-steel instrument tray with a mesh bottom. Place the tray in the ultrasonic cleaning tub and make sure that in the bath, at 24°C ±2°C (75,2°F ±3,6°F), the instruments are submerged in fresh enzyme solution (see step 1).

- Ultrasonicate the torque wrench for at least 20 minutes and the tip holder and tip for 10 minutes, or as instructed by the manufacturer of the enzyme cleaner and ultrasonic tank.

- Remove the device from the ultrasonic cleaning tank.

1.2- Automatic Cleaning

NOTE: The tip holder, tip, and torque wrench must be processed separately. Do not leave the tip installed on the tip holder. Do not leave the tip holder inserted into the torque wrench.

Automatic cleaning involves the use of a washer-disinfector¹⁾ and the following materials:

- Alkaline detergent, for example Neodisher Mediclean Forte;
- Neutralizing liquid, for example Neodisher Z (0.1 % v/v);
- Metal basket;
- Adapters for washer-disinfector

¹⁾For example AT-OS instrument washer-disinfector.

NOTE: Make sure that the accessories are appropriately blocked in the basket and cannot move during washing. Any blows could damage them. Position the instruments in a way that the water can flow through all the surfaces, even internal.

⚠ WARNING: Avoid overloading the thermal disinfector as this could compromise cleaning effectiveness.

- Place the accessories in a metal basket using the special adaptors (supplied as optionals).
- Set the sequence and parameters for the washer-disinfector cycle:
 - 3 min, Prewash with cold deionized water;
 - 5 min, Wash with alkaline detergent and deionized water at 55°C ±2°C (131°F ±3,6°F);
 - 2 min, Neutralization with suitable solution and deionized water at 40°C ±2°C (104°F ±3,6°F);
 - 2 min, Rinse with cold deionized water at 32°C ±2° (89,6°F ±3,6°F);
 - 5 min, Termodisinfection at 93°C (199,4°F) with deionized water.

Automatic thermal disinfection is not experimentally tested. In compliance with ISO 15883-1, Table B.1, thermal disinfection at a temperature of 90°C (194°F) for 5 min results in an A0 3000 value.

⚠ CAUTION: After cleaning

The use of alkaline detergents removes the lubricant from the torque wrench, causing malfunctions and increasing its wear. Therefore, it is important to lubricate the torque wrench after each cleaning procedure using medical-grade lubricants (see chapter "2- Lubrication"). The dilution of the alkaline detergent required for automatic cleaning is such that it does not damage the silicone parts of the torque wrench (O-ring).

- Make sure that the device is completely dry even internally before starting the sterilization cycle. Use filtered compressed air to flush out any moisture from cavities, grooves, crevices, and other hard-to-reach areas. This prevents the appearance of stains or streaks on the surface and oxidation.
- Check

➤ **CAUTION:** Do not grip the handpiece by the end and/or cord, but only by the body. Do not rotate the body while securing the tip holder. Firmly grip the body of the handpiece and turn the torque wrench only.

6.2- Removing the Tip holder

Position the torque wrench on the tip holder and screw in an anti-clockwise direction.

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Piezosurgery, Inc.