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English

INSTRUCTIONS FOR USE

This insert has been specifically designed to be used in conjunction with bone surgical devices manufactured by Mectron S.p.A. according to their intended use/applications.

The insert is a bone surgery instrument used in applications including osteotomy, osteoplasty, soft tissue elevation (e.g. periosteum), endodontic and periodontal surgery, scaling, extraction and implantology.

This insert must be used in accordance with the correct settings described in the attached leaflet "APPROPRIATE SETTINGS FOR INSERTS".

The insert must be used exclusively by specialized and duly trained personnel such as the Surgeon / Dentist, adult of any weight, age, height, gender and nationality, non-disabled.

NOTE: Regarding the intended patient group, patient selection criteria, indications for use, contraindications, and warnings, please refer to the use and maintenance manual supplied with your oral surgery device manufactured by Mectron.

IMPORTANT: The words **"WARNING"** and **"CAUTION"** contained in this document carry special meaning and must be carefully reviewed.

⚠. WARNING: Identifies conditions and practices that present a risk of death or serious injury for the patient and/or user.

⚠. CAUTION: : Identifies conditions and practices that may cause minor injuries and/or damage to the devices.

NOTE: A condition not associated with a risk or imminent danger.

IMPORTANT: Information regarding this user guide. The information contained in these operating instructions is limited to that which is essential for the first use of the insert and must be employed in conjunction with the use and maintenance manual supplied with your oral surgery device manufactured by Mectron.

It is therefore necessary to carefully read the use and maintenance manual supplied with your oral surgery device manufactured by Mectron, before starting the treatment, set up the entire system and check each function, paying special attention to the "Safety Precautions" paragraphs.

⚠. WARNING: Before first use and subsequent uses. The insert is not sterile upon receipt and must be treated, before use, according to the procedures described below and/or reported in the "Cleaning and Sterilization" manual supplied with your oral surgery device manufactured by Mectron.

Once used, and prior reuse, reprocess it again thoroughly following the procedures described below or contained in the "Cleaning and Sterilization" manual supplied with your Mectron oral surgery device.

⚠. WARNING: If an adverse event and/or serious accident attributable to the device occurs during correct and intended use, it is recommended to report it to the Competent Authority and to the manufacturer indicated on the product label.

1 - Cleaning Procedures

IMPORTANT: Cleaning processes must immediately be carried out after each use. Do not allow contaminated inserts to dry prior reprocessing.

The cleaning procedure can be performed using two different methods: manual cleaning and automatic cleaning. These two methods are alternatives.

Limitations on repeated processes: Repeated processing has a minimal effect on these instruments. The end of their useful life is generally determined by wear and damage due to use.

⚠. CAUTION: Do not use hydrogen peroxide to clean the inserts.

⚠. WARNING: To avoid injury to operators, treat inserts separately.

Special caution is necessary due to the WARNING of injury and infection from handling soiled cutting tools. Healthcare personnel involved in device operation and cleaning must be appropriately protected. Handle sharp-edged and pointed instruments with care.

1.1 - Manual Cleaning

NOTE: Process validated by independent body with Enzymec enzymatic detergent.

1 Prepare a neutral pH (6-9) enzymatic detergent solution according to the manufacturer's specifications. (preferably use warm water).

⚠. CAUTION: Once used, dispose of the Enzymec detergent solution properly, do not recycle.

2 Place the insert in a clean container, add enough solution to the container so as to cover the insert. Soak the insert in the enzymatic detergent solution for at least 10 minutes.

3 After immersion in the enzymatic solution, gently brush the surface of the insert with a soft-bristled nylon toothbrush until all visible soil has been removed. Use a clean nylon soft-bristled toothbrush for the outer surfaces, a clean nylon soft-bristled brush for the inner cavities and crevices. Thoroughly clean difficult areas such as sharp edges, and particularly the interstices between the cutting cusps.

4 Rinse the insert under running water and brush the device with a clean nylon soft-bristled toothbrush for at least 1 minute. Repeat this procedure until all visible soils has been removed.

5 Using a syringe (min. 20 ml), inject the enzymatic detergent solution three times into the insert cavity to effectively remove residues from the inner surface.

6 Place the insert in the basket in the ultrasonic clean tank at room temperature; fully immerse the insert in the fresh enzymatic detergent solution (refer to point 1) and sonicate for at least 10 minutes.

7 Rinse the insert under running tap water and brush the device with a clean nylon soft-bristled toothbrush for at least 1 minute.

8 Rinse the lumen by injecting demineralized water three times using a 20 ml syringe and then place the insert in a clean container covered with demineralized water for at least 1 minute.

9 Blow out the internal channels and all surfaces with compressed air until the insert appears completely dried.

1.2 - Automatic Cleaning

Automatic cleaning involves the use of a washer-disinfector and the following materials:

- Alkaline detergent ¹⁾;
- Neutralizing liquid ²⁾;
- Metal basket;
- Washer-disinfector ³⁾

- Adapters for washer-disinfector

^[1] For example Neodisher Mediclean Forte.

^[2] For example Neodisher Z (0.1 % v/v).

^[3] For example AT-OS instrument washer-disinfector.

NOTE: Make sure that the accessories are properly locked in the basket and cannot move during washing. Any shock could damage them. Place the tools so that water can flow through all surfaces, including internal ones.

⚠. WARNING: Avoid overloading the washer-disinfector,as this could compromise cleaning effectiveness.

Use ISO 15883 compliant washers-disinfectors with baskets suitable for holding small products and with flush connectors. Place the product correctly in a suitable basket, connect the lumens to the rinse fittings and start the automatic cleaning.

1 Place the insert in a metal basket. Attach the appropriate adapter for the inserts (supplied as an option).

2 Set the sequence and parameters for the washer-disinfector cycle:

- 3 min, Pre-wash with cold water;
- 5 mins, Wash with alkaline detergent at 55°C ±2°C (131°F ±35.6°F);
- 2 min, Neutralization with suitable solution;
- 2 min, Rinse with cold water;
- 5 min, Thermal disinfection at 93°C (199.4°F) with demineralized water.

Automatic thermal disinfection is not experimentally tested. In compliance with ISO 15883-1, Table B.1, thermal disinfection at a temperature of 90°C (194°F) for 5 min results in an A0 3000 value.

⚠. CAUTION: Make sure the insert is completely dry even internally before starting the sterilization cycle. Use filtered compressed air to remove any moisture from the through-hole. This prevents the appearance of stains or halos on the surface and the internal oxidation of the insert.

⚠. WARNING: Check before sterilization. The inserts must be inspected prior to sterilization. A visual inspection conducted in good lighting conditions and generally without the use of magnification means is sufficient. All parts of the inserts should be checked for visible soil, damage and/or corrosion. All the inserts should be carefully inspected visually for damage and signs of wear. Particular attention should be paid to:

- Places where soil can get trapped, such as:
 - holes
 - threaded parts,
 - recessed features (holes, cannulations),
 - sides of insert teeth, if any,
 - cutting edges sharpness.
 - check edges that may be deformed or largely affected areas. The edges must be continuous.

If the surfaces are not visually clean, repeat the cleaning and inspect again. Discard damaged tools.

1.3 - Packaging

The insert must be sterilized using a medical-grade pouch suitable for autoclave sterilization in compliance with the requirements of ISO 11607-1.

⚠. CAUTION: Each insert must be packed individually.

⚠. CAUTION: Make sure the pouch is large enough to hold the instrument without tensioning the seals and without tearing the package.

1.4 - Sterilization

Sterilize only by pre-vacuum steam sterilization in autoclave.

⚠. CAUTION: If it is necessary to autoclave several instruments in a single cycle, do not exceed the maximum allowable load.

⚠. CAUTION: Risk of contamination. Do not use gravity autoclaves to sterilize inserts. The operating cycle of gravity autoclaves does not guarantee the suitable sterilisation of the internal channel, of the cavities and hard-to-reach points.

⚠. CAUTION: Sterilize using only a steam autoclave. Do not use any other sterilization process as it may be incompatible with the materials used to produce the instruments (dry heat, irradiation, ethylene oxide, gas, low-temperature plasma, etc.).

⚠. CAUTION: Do not use the following substances/systems to sterilize instruments: hydrogen peroxide, peracetic acid, Formaldehyde, Gluteraldehyde, or other equivalent solutions/systems. Perform sterilization immediately after cleaning and disinfection. For the sterilization of products with moist heat, comply with the instructions contained in the ISO 17665-1 standard and the legal regulations in force in your country. Follow the directions, instructions, and warnings of the steam sterilizer manufacturer.

⚠. WARNING: The sterilization process validated by Mectron S.p.A., in a steam autoclave, guarantees a SAL 10-6 by setting the parameters listed below:

- 3 times Pre-vacuum.
- Sterilization temperature 132 °C (range 0 °C ÷ +3 °C).
- Sterilization time 4 minutes.
- Minimum drying time 10 minutes.

⚠. CAUTION: Do not exceed the temperature of 135 °C.

2 - Important information regarding the safety of inserts

⚠. WARNING: Breakage and wear of the insert.

In rare cases, the high frequency oscillations and the wear may lead to insert breakage. Do not bend, reshape, and re-sharpen inserts in any way. Bending, forcing, and sharpening can cause the insert to fracture during use. Any insert modified in this way should definitely not be used. If an insert breaks during use, extreme care must be exercised to ensure that all insert fragments are retrieved and removed from the patient. In oral surgery it is necessary to instruct the patient to breathe through his/her nose during the treatment, or use a dental dam, in order to avoid ingestion of broken off fragment insert fragments.

⚠. WARNING:

The insert is consumable. Thoroughly inspect the insert before and during each treatment for damage and/or excessive wear (see table "INSERTS WEARING OUT EXAMPLES"). Do not use an insert if it is scratched or rusty. The insert may break during use. If you notice any damage or loss in performance, replace the insert with a new one.

NOTE: If an insert is deformed, the device may stop to operate, and a message may appear on the display.

- Before treatment, make sure that the chosen insert is firmly attached to the handpiece. Properly tighten the insert on the handpiece with the correct torque wrench supplied with the device.
- Do not change the shape of the insert in any way. Bending or forcing the insert may cause it to fracture. Never use bent inserts.

replace the insert(s) only with original Mectron spare parts. The use of inserts other than the original ones will affect the warranty of the device. Never use inserts other than the original ones and never use inserts with incompatible devices as they will damage the device and may cause injury to the operators or patient. The use of inserts other than the original ones will damage the threaded pin of the handpiece with the consequence that the attachment of the inserts, even the original ones, will no longer be secure. The use of inserts on non-compatible devices will damage the threaded pin of the handpiece with the consequence that the attachment of the inserts, even the original ones, will no longer be secure.

- It is advisable to avoid applying excessive force or prolonged contact of the insert to soft tissues to prevent thermal damage and/or injury.
- Allow the ultrasonic vibrations to work, do not exert excessive pressure on the inserts during use. Apply a gentle pressure to the insert for best efficiency
- Applying excessive pressure could cause the insert to fracture which in turn can lead to injury.
- The insert must be kept always moving. If the insert is blocked, thermal damage to the treated part can occur. It is recommended to use a continuous motion to minimize contact between the tip and the part. It is recommended to use high levels of irrigation as the power level increases.
- Mectron inserts vibrate with a longitudinal oscillation, with back-and-forth movement.
- Check the threaded parts of the insert and those on the handpiece. These parts must be cleaned thoroughly.
- Nitrided inserts should be replaced as soon as wear on the titanium nitride coating becomes visible.
- The use of a worn insert will lower the yield in terms of cutting power.

⚠ WARNING: Hospital waste
Treat worn or damaged inserts as hospital waste.

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Piezosurgery, Inc.