Italiano

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English

INSTRUCTIONS FOR USE

This insert is designed specifically to be used in conjunction with all Piezosurgery® devices manufactured by Mectron S.p.A. according to their intended use / applications. Make references to the user manual booklet supplied with the device - paragraph "Intended use". This insert must be used in accordance with the proper settings described in "Table 1".

The insert must be used exclusively by specialized and duly trained personnel, such as the Surgeon/Dentist, adult of any weight, age, height, gender and nationality, non-disabled.

With regard to the intended patient group, patient selection criteria, indications for use, contraindications and warnings, refer to the use and maintenance manual supplied with the Piezosurgery® device, produced by Mectron.

IMPORTANT: The words "Warning", "Caution" and "Note" contained in this document carry special meaning and must be carefully reviewed.

⚠ **WARNING:** Implies that death or serious injury could occur if the advice is not followed.

⚠ **CAUTION:** Implies that minor injury or device damage could occur if the advice is not followed.

NOTE: Implies advice that is not related to harm.

IMPORTANT: Information regarding this user guide. This user guide provides only basic essential information regarding the first use of the insert tip and must be employed in strict coordination with the "Instructions for use and maintenance" booklet provided with your device. Therefore, before starting the treatment read carefully the "Instructions for use and maintenance" booklet, set up the entire system and verify each function. Pay special attention to the paragraph "Safety Requirements".

⚠ WARNING: First and every use

The insert is not sterile upon receipt and must be prepared prior the use by applying the procedures described below and/or at paragraph "Cleaning and Sterilization" of the user manual booklet provided with the unit. Once used, and prior reuse, the insert must be reprocess thoroughly by following the processes described below and/ or at paragraph "Cleaning and Sterilization" of the user manual booklet provided with the unit.

⚠ **WARNING:** If an adverse event and/or serious accident attributable to the device occurs during correct and intended use, it is recommended to report it to the Competent Authority and to the manufacturer indicated on the product label.

1- Intended use

Extraction of teeth, fractured roots, root fragments.

The EXL inserts series is composed of 3 different types according to curvature, conformation and sizing of the stem and the apical operating part:

- EXL1 – "Hook" curvature;
- EXL2 – Robust "S" curvature and short stem;
- EXL3 – Slim "S" curvature and long stem;

The shape of the stem allows performing different extraction procedures according to the position of the tooth/ root and the anatomical condition.

Once the insert is placed in the correct position, the sharp edges of the operative part engage the tooth/root to be extracted while the opposite smooth and chamfered side offers support to exert the lever effort which in combination with the micro-vibrations transmitted to the tooth/root will cause rapid ligament dislocation which will facilitate extraction.

⚠ **WARNING:** For a correct extraction technique, it is recommended to carefully inspect the tooth morphology through RX / TC scan, to determine its position/size and do not to cause damage.

2- Operation

1. On the device, set the parameters indicated in Tab. 1.

⚠ **CAUTION:** The recommended settings for power and irrigation reported in Tab. 1 allow minimizing the thermal stress of the anatomical support part.

2. The insert must engage the tooth at the root level and have the point of support on the adjacent tooth or on the interdental septum in the same way as a conventional manual lever.

3. The "hook" lever can be used for dislocation when applied as a wedge between the roots. Subsequently engaging the roots by hooking them they can be extracted by pulling them along their axis.

4. Start the device, always with the insert not in contact with the parts, engage the tooth / root and apply the lever effort alternately in the mesio-distal direction in order to facilitate dislocation and therefore extraction.

5. The lever effort must be balanced to allow the insert to always vibrate, despite the effort, and to transmit the vibrations to the tooth/root.

6. In the specific case of extraction of roots/root fragments, the operating procedure within the alveolus follows the same precautionary criteria in applying the effort.

7. In the specific case of extraction of the 3rd molar, it is recommended to first perform the total or partial coronotomy and any rhizectomy. The support side of the insert will lie on the distal surface of the 2nd molar or between the roots of the 3rd molar.

⚠ **CAUTION:** The normal working conditions, in contact with the tooth/root surface, will inevitably cause the deterioration of the insert.

⚠ **WARNING:** During extraction, do not force the insertion of the operative part of the insert into anatomical spaces that are not large enough to accommodate it.

⚠ **WARNING:** During extraction, do not use as point of support of the insert the interdental septa or thin alveolar walls, because they are more sensitive to lever efforts and overheating.

⚠ **WARNING:** The irrigation flow is directed towards the point of support of the lever and must be high enough to cool it and avoid damage to the tissues caused by excess temperature.

⚠ **WARNING:** During extraction, do not keep the insert in the same position for a long time, but move it alternatively in opposite directions to avoid overheating of the parts and to facilitate dislocation.

⚠ **WARNING:** During extraction, do not apply lever efforts in the vestibul-palatal direction in the presence of thin alveolar walls that are more sensitive to this type of effort.

⚠ **WARNING:** During extraction, do not exceed in the lever effort so as not to block the insert vibration and its transmission to the tooth/root because it could cause overheating of the parts in contact.

3- Main warnings for cleaning procedures

IMPORTANT: The cleaning procedures must be performed immediately after each use. Do not let the contaminated inserts dry before the cleaning treatment.

Once used, the insert must be reprocessed thoroughly by following the processes described at paragraph "Cleaning and Sterilization" of the user manual provided with the unit.

4- Sterilization

Sterilise only by pre-vacuum steam autoclaving.

The steam sterilisation parameters have been validated by Mectron S.p.A. to guarantee a sterility level (Sterility Assurance Level - SAL) of 10⁻⁴.

⚠ **CAUTION:** If multiple instruments need to be sterilised in the autoclave in a single cycle, do not exceed the maximum allowed load.

⚠ **CAUTION: Risk of contamination.** Do not use gravity autoclaves to sterilise the inserts. The operating cycle of gravity autoclaves does not guarantee the suitable sterilisation of the internal channel, of the cavities and hard-to-reach points.

⚠ **CAUTION:** Use only a pre-vacuum autoclave to sterilise the inserts and the torque wrench. Do not use other methods of sterilisation insofar as they may be incompatible with the materials used to produce the instruments.

Do not use the following methods of sterilisation: ethylene oxide sterilisation, hot air sterilisation, flash autoclaving, STERRAD sterilisation, sterilisation with STERIS or similar sterilisation systems. **Do not use the following substances/systems to sterilise the instruments:** oxygenated water, peracetic acid, Formaldehyde, Glutaraldehyde based systems or other equivalent solutions/systems.

4.1- Sterilisation - Parameters

Minimum sterilisation parameters to guarantee a sterility level (Sterility Assurance Level - SAL) of 10⁻⁴:

⚠ **WARNING:** these sterilisation parameters, validated by an independent laboratory, refer exclusively to a PRE-VACUUM autoclave operating cycle.

- Type of cycle: PRE-VACUUM - (3 pre-vacuum stages - 60mbar).
- Minimum temperature: 132 °C + 3 °C (270 °F + 5 °F).
- Pressure: 2,9 bar (42 psi).
- Minimum exposure time (in package): 4 minutes.
- Minimum drying time: 20 minutes.

⚠ **CAUTION:** Do not exceed a temperature of 135 °C.

All stages of sterilisation must be performed by the operator in compliance with the current revision standards: UNI EN ISO 17665-1, UNI EN ISO 556-1 and ANSI/AAMI ST-46.

5- Important information relative to safety of inserts

⚠ **WARNING: Breakage and wear of insert.**

In rare cases, high frequency oscillations and wear may lead to insert breakage.

Do not bend, reshape and re-sharpen inserts in any way. Bending, forcing and sharpening can cause the insert to fracture during use. Any such inserts should definitely not be used. Should an insert fracture during use, extreme care must be exercised to ensure that all insert fragments of the insert are retrieved and removed from the patient. In oral surgery it is necessary to instruct the patient to breathe through his/her nose during the treatment, or use a dental dam, in order to avoid ingestion of broken off fragment insert fragments.

⚠ **WARNING:** the insert is consumable. Carefully inspect the insert before and during each treatment for any damage and/or for excessive wear. Do not use insert that is scratched or rusted. In use, the insert may break. If a damage or drop in performance is noted replace with a new insert.

NOTE: Should an insert be deformed, the device could stop functioning and a massage could be displayed on the device.

- Before treatment, make certain the appropriate insert is securely attached to the handpiece. Use the torque wrench provided with the Piezosurgery® device for securing the insert to the handpiece.
- Replace the insert/s exclusively with original Piezosurgery® spare parts. Use of non-authentic Piezosurgery® inserts will void the warranty. Never use non-authentic Piezosurgery® inserts as they will damage your Piezosurgery® device and may cause injury to you or your patient. The use of non-authentic Piezosurgery® insert/s damages the threaded pin of the handpiece, with consequent poor fastening of original insert on it.
- It is advisable to avoid applying excessive force or prolonged contact of the insert on soft tissues to avoid thermal damage and / or injury.
- Applying excessive pressure could cause the insert to fracture which in turn can cause injury.
- The insert must always be kept in motion. If the insert is blocked, thermal damage of the treated part can be generated.
- Check the treated parts of the insert and the ones of the handpiece. They must be cleaned.

⚠ **WARNING: Hospital waste**

Treat the insert/s, when worn or broken as hospital waste.

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