mectron

medical technology

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USE AND MAINTENANCE MANUAL

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IT

EN

FR

MICROPIEZO

CE
0051

mectron

medical technology

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IT

MICROPIEZO

Copyright

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MICROPIEZO

Copyright

© Mectron S.p.A. 2024. All rights reserved. No part of this document may be reproduced, in any form whatsoever, without the prior written consent of the copyright holder.
Images are for demonstration purposes only.

CONTENTS

1	Introduction	1
1.1	Intended Use	1
1.2	Description of the Device	2
1.2.1	Eligible Patient Group	2
1.2.2	Patient Selection Criteria	2
1.2.3	Indications for Use	2
1.2.4	Users	2
1.3	Disclaimer	3
1.4	Safety Requirements	4
1.5	Symbols	6
2	Identification Data	8
2.1	Device ID Plate	8
2.2	Handpiece Identification Data	8
2.3	Inserts Identification Data	9
3	Delivery	9
3.1	List of Components	9
4	Installation	12
4.1	First Installation	12
4.2	Safety Requirements in the Installation Phase	12
4.3	Connecting Accessories	13
5	Use	16
5.1	Power On/Off	16
5.2	Description of the Keyboard	17
5.2.1	Functions and Power Levels of the Treatment	17
5.2.1.1	Soft Mode	18
5.2.2	Water Delivery	18
5.2.3	Filling and Cleaning the Water Circuit	19
5.2.3.1	FLUSH - Short Mode	19
5.2.3.2	FLUSH - Long Mode	19
5.2.3.3	FLUSH - Interruption	19
5.2.4	Light	20
5.2.5	Wi-Fi - IoT Technology	20
5.2.5.1	Connection to a Wi-Fi Network	20
5.2.5.2	Wi-Fi Network Requirements	21
5.2.5.3	Software Update	21
5.2.6	Diagnostic System	21
5.3	Safety Requirements Before and During Use	22
5.4	Instructions for Use	25
5.5	Important Information on the Inserts	28
6	Long Flush function	29
7	Disassembling Parts for Cleaning and Sterilisation	31
8	Cleaning and Sterilisation	33
8.1	Preparation	34
8.2	Manual Cleaning with Enzymatic Detergent	35
8.3	Manual Cleaning with Enzymatic Detergent and Running Water	36
8.4	Disinfecting the Tank Irrigation Circuit	38
8.5	Cleaning the Sterilisable Accessories	43

8.5.1	Manual Cleaning	43
8.5.1.1	Scaler Handpiece	44
8.5.1.2	Inserts	47
8.5.1.3	Inserts Tightening Wrench	49
8.5.2	Automatic Cleaning	51
8.6	Cleaning Control	53
8.7	Drying and Lubrication	54
8.8	Sterilisation	55
8.8.1	Special Information	56
9	Maintenance	57
9.1	Replacing the Peristaltic Pump	57
10	Methods and Precautions for Disposal	59
11	Technical Data	59
11.1	Software	60
11.2	Electromagnetic Compatibility IEC/EN 60601-1-2	61
11.3	Guide and Manufacturer's Declaration - Electromagnetic Emissions	61
11.4	Accessible Parts of the Casing	62
11.5	Guide and the Manufacturer's Declaration - Electromagnetic Immunity	63
11.5.1	Power Connection BC Input	63
11.5.2	Points of Contact with the Patient	65
11.5.3	Parts Accessible to the Input / Output Signals	66
11.6	Specifications of the Tests for the Immunity of the Accessible Parts of the Casing to the Wireless RF Communications Device	67
11.7	Immunity to Proximity Magnetic Fields in the Frequency Range 9 kHz to 13,56 MHz	68
12	Troubleshooting	69
12.1	Diagnostic System and Symbols on the Keyboard	69
12.2	Quick Troubleshooting	73
12.3	Shipping to an Authorised MECTRON Service Centre	75
13	Warranty	76

PAGE INTENTIONALLY LEFT BLANK

1 INTRODUCTION

Carefully read this manual before proceeding with installation operations, use, maintenance or other intervention operations on the device. Always keep this manual at hand.

Important: To avoid damages to persons or property, carefully read all the "Safety Requirements" sections in the manual.

According to the degree of severity the safety requirements are classified with the following indications:

 **WARNING:** Always referred to damage to persons

 **CAUTION:** Referring to possible damage to property

The aim of this manual is to make the operators aware of the safety regulations, installation procedures, instructions for proper use and

maintenance of the device and its accessories. Use of this manual for aims other than those strictly linked to the installation, use and maintenance of the device is prohibited.

The information and illustrations in this manual were updated on the edition date shown on the last page.

MECTRON is engaged in continuously updating its products with possible changes to the components of the apparatus.

In case you encounter discrepancies between the descriptions found in this manual and the equipment in your possession you can:

- check for any available updates in the *MANUALS section of the MECTRON website*¹;
- ask clarifications to Your Dealer;
- contact MECTRON After Sales Service.

1.1 Intended Use

MICROPIEZO is a piezoelectric ultrasonic scaler.

With the appropriate inserts it is possible to perform the following treatments:

- **Scaling:** all the procedures for the removal of deposits of bacterial plaque and supragingival, subgingival and interdental stones, and removal of stains;
- **Periodontics:** periodontal therapy for scaling and root-planing/debridement, including cleaning and irrigation of the periodontal pocket;
- **Cleaning treatment of the implant surface;**
- **Endodontics:** all the treatments for preparation of the channels, irrigation, filling, condensation of gutta-percha, endodontic reprocessing and retrograde preparation;
- **Restorative and prosthesis:** preparation of the cavities and removal of caries, removal of the prosthesis and of restoration materials in excess, condensation of the amalgam, finishing of the prosthetic stump.

The Mectron piezoelectric ultrasonic scaler devices can be used on patients of any age and gender requiring a tooth polishing treatment. There are no contraindications for specific population segments.

 **WARNING:** The device must be used in a dental practice or clinic, or in professional oral hygiene and preventive care centres. Do not use the device in environments where the atmosphere is saturated with flammable gases (anaesthetic mixtures, oxygen, etc.).

 **WARNING: Qualified and specialised staff.**

The device is intended exclusively for use by a dental surgeon, dentist, or dental hygienist, with suitable medical/health training.

¹ <https://manuals.mectron.com/>

⚠ WARNING: Intended use.

Use the device exclusively for the use for which it is intended. Failure to observe this regulation can cause serious injury to the patient, the operator and damage/faults on the device.

1.2 Description of the Device

MICROPIEZO is a multi-purpose piezoelectric ultrasonic scaler designed to equip the operator with a product featuring an innovative design and unique technical characteristics, and to ensure patient comfort during treatment.

The user interface is optimised to make all functions immediately available, by integrating

them into the touch keyboard. The handpiece is equipped with a circular LED light and is autoclavable at 135°C. The device is equipped with an automatic tuning circuit that optimises power and frequency for each insert available, in order to always operate in conditions of maximum efficiency.

1.2.1 Eligible Patient Group

This medical device is designed to be used with the following patient groups:

- Children;
- Adolescents;
- Adults;
- Seniors.

This medical device can be used on patients of any age, weight, height, gender, and nationality, where applicable.

1.2.2 Patient Selection Criteria

It is not recommended to use the device in the following cases:

1. Patients with active implantable medical devices (for example: pacemakers, hearing aids, and/or other electromagnetic prostheses) without prior authorisation from their doctor;

2. Patients with clinical conditions not suitable for treatment of the sites (for example: local anaesthesia).

All piezoelectric ultrasonic scaler device models are intended for professional use only. The user is therefore the only person able to decide if and how to treat their patients.

1.2.3 Indications for Use

Use of the device is suitable for all eligible patients (see *Chapter 1.2.1 on page 2*) whose doctor has prescribed professional scaling, periodontology, endodontic cleaning, implant site surface cleaning, a restorative or prosthetic technique, as part of the intended use of the device (see *Chapter 1.1 on page 1*).

1.2.4 Users

The device must be used exclusively by specialised and properly trained personnel, specifically the dental surgeon/dentist or dental hygienist, who must be able-bodied, adult, of any weight, age, height, gender, and nationality.

1.3 Disclaimer

The manufacturer MECTRON disclaims all liability, express or implied, and cannot be held responsible for direct or indirect personal injury and/or property damage, occurring as a result of incorrect procedures linked to the use of the device and its accessories.

The manufacturer MECTRON cannot be held liable, expressly or implicitly, for any type of personal injury and/or property damage, inflicted by the user of the product and its accessories and which occurs in the following cases, by way of example and not of limitation:

- Misuse or use during procedures other than those specified in the destination of use of the product;
- The environmental conditions for preservation and storage of the device are not complying with the requirements indicated in *Chapter 11 on page 59*;
- The device is not used in accordance with all the instructions and requirements described in this manual;
- The electrical system of the places where the equipment is used do not comply with the laws in force and the related regulations;
- Assembly operations, extensions, re-adjustments, upgrades and repairs of the device are carried out by personnel not authorised by MECTRON;
- Misuse, abuse, abnormal use, negligent use, intentional misconduct or use exceeding the limits of the device indicated and allowed and/or normal wear or deterioration, ill-treatment and/or incorrect interventions;
- Any attempt to tamper with or modification of the device under every circumstance;
- Use of non original MECTRON inserts resulting in permanent damage to the thread of the handpiece with impaired functioning and risk of injury to the patient;
- Use of non original MECTRON inserts used in accordance with the settings designed and tested on the original MECTRON inserts. Correct use of the settings is only guaranteed with original MECTRON inserts;
- Shortage of stock material (handpiece, inserts, wrenches) to be used in case of a standstill due to faults or problems;
- Incorrect/omitted maintenance compared to what is stated in *Chapter 9 on page 57* of this manual;
- Breach of the requirements and the information contained in *Chapter 5.5 on page 28* of this manual;
- Breach of the requirements and the information contained in *Chapter 8 on page 33* of this manual;
- Unauthorised repairs in accordance with the indications contained in *Chapter 12.3 on page 75* this manual.

1.4 Safety Requirements

WARNING: Contraindications.

Do not use the device on patients with Pacemakers or other implantable electronic devices. This regulation also stands for the operator.

 **WARNING: Contraindications.** Do not perform scaling treatments without water spray in order to avoid an overheating of the insert that can cause damage to the tooth. Treatments without water spray can only be those performed with the "Dry Work" inserts without water passage.

 **CAUTION: Contraindications. Ultrasonic scaler.** Do not perform treatments on metal or ceramic prosthetics. Ultrasonic vibrations could lead to decementation of the dentures.

WARNING: Contraindications. Interference from other equipment.

An electrosurgical scalpel or other electrosurgical units positioned near the MICROPIEZO device may interfere with the correct operation of the device itself.

 **WARNING: Contraindications. Interference with other equipment.** Even if compliant with standard IEC 60601-1-2, MICROPIEZO may interfere with other devices in the vicinity. MICROPIEZO must not be used near or stacked with other devices. However, if this should become necessary, it is necessary to check and monitor the correct functioning of the device in that configuration.

 **WARNING: Risk of explosion.** The device cannot operate in environments where the atmosphere is saturated with flammable gases (aesthetic mixtures, oxygen, etc.).

 **CAUTION:** In case the final user, operating in their own medical room or surgery, in order to comply with mandatory requirements, must periodically inspect the equipment present in the surgery, the test procedures to apply to medical electrical equipment and medical electrical systems for the safety assessment must be carried out following the standard EN 62353 'Medical electrical equipment - Periodic inspections and tests to be carried out after repair of medical electrical equipment'. The frequency of periodic inspections in the intended conditions of use described in this "Use and Maintenance" manual is once per year or every 2000 hours of use, whichever condition is satisfied first.

WARNING: Checking the condition of the device before treatment.

Always check that there is no water under the device. Before every treatment always check the perfect operation of the device and efficiency of the accessories. If operation anomalies are observed do not carry out the treatment. If the faults concern the device, contact an Authorised MECTRON Service Centre.

 **CAUTION:** The electrical system of the premises in which the device is installed and used must comply with the rules in force and the relevant requirements of electrical safety.

 **CAUTION:** To avoid the risk of electrical shocks, this device has to be connected exclusively to a power supply with an earth safety device.

 **WARNING: Cleaning and sterilising new or repaired instruments.** All the new or repaired device accessories are not sterile. At first use and after each treatment they must be cleaned and sterilised by carefully following the instructions in *Chapter 8 on page 33*.

⚠ WARNING: Infection control.

For maximum patient and operator safety, before using all the parts and reusable accessories, make sure they have been previously cleaned and sterilised by following the instructions in *Chapter 8 on page 33*.

⚠ CAUTION: Contraindications. After having autoclaved the handpiece, inserts, torque wrench or any other sterilisable accessory, wait for them to completely cool down before reusing them.

⚠ WARNING: Insert Breakage and Wear. High frequency oscillations and wear can, in rare cases, lead to the breakage of the insert. Deformed or damaged inserts are susceptible to breakage during use. Broken or worn inserts must not ever be used. In the event of breakage check that no fragments remain in the treated part and at the same time use suction effectively to remove them.

It is necessary to instruct the patient to breathe through the nose during the treatment, or use a dental dam, so as to avoid to ingest fragments of broken inserts. Check the state of wear of the insert and its integrity before and during each use. If a drop in performance occurs, see to its replacement.

The state of wear of the most common inserts (S1, S1-S, S2, S5, P2, P4, P10) can be checked using the supplied INSERT-CARD. To use the INSERT-CARD correctly:

- Position the insert on the INSERT-CARD so that the profile corresponds to the one printed on the card. The printed profile on the card has a red line

indicating the limit of wear;

- If the insert is shorter than the limit of wear, its performance will be significantly inferior compared to that of a new insert, and should therefore be replaced.

If the layer of titanium nitride (gold-plated surface), where present, is visibly worn, the insert must be replaced. The use of a worn insert reduces its efficiency.

Diamond Inserts: diamond inserts must be replaced when the layer of titanium nitride is visibly worn and in any case after a maximum of 10 treatments.

When the nitriding is consumed, the cutting edge loses effectiveness; a possible regrounding harms the insert, therefore it is prohibited. Verify that the insert is not worn. During the operation frequently check that the insert is intact, especially in the apical part.

During the operation avoid prolonged contact with the retractor or with metal instruments in use. Do not exert excessive pressure on the inserts during use.

⚠ WARNING: Exclusively use original MECTRON inserts, accessories and spare parts.

⚠ CAUTION: No modification is allowed to this device.

⚠ WARNING: In case of an adverse event and/or accident attributable to the device during correct use and in accordance with the intended use, a report must be made to the Competent Authority and manufacturer indicated on the product label.

1.5 Symbols

Symbol	Description	Symbol	Description
	Class IIa device compliant with Regulation (EU) 2017/745. Notified body: IMQ S.p.A.		Nemko brand Compliance with UL - CSA regulations
	Medical Device		Caution
	Consult instructions for use or consult electronic instruction for use		Manufacturer
	Date of manufacture		Serial Number
	Batch Number		Catalogue number
	Unique Device Identifier		Model number
QTY.1	Quantity in package: 1		Health Industry Bar Code
	Non-sterile		Sterilizable up to a max. temperature of 135°C
	Type B applied part		Alternating current
	Direct current		Polarization of the connector.
	Power button.		Connection of the control pedal
	Class II equipment		For internal use only.
	Atmospheric pressure limitation		Humidity limitation

Symbol	Description	Symbol	Description
	Temperature limit		Water flow
	Power level		Soft Mode
IP 22	Protection degree of the case	IP 20	Protection degree of the case
	Keep dry		Fragile
Rx Only	For US market only CAUTION: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist or physician.		The device and its accessories must not be disposed of or treated as municipal solid waste

NOTE: For other symbols, refer to *Chapter 12.1* on page 71

2 IDENTIFICATION DATA


A correct description of the model and of the serial number of the device will allow the After Sales Service to provide fast and effective answers.

Always provide this information every time that you contact MECTRON After Sales Service.

2.1 Device ID Plate

Each device is equipped with an identification plate indicating the main technical characteristics and traceability data, UDI code included. The identification plate is located under the device. The complete technical specifications are provided in Chapter 11 on page 59.

NOTE: The complete list of symbols and their description is provided in Chapter 1.5 on page 6.

Additional symbols and characteristics of the device are provided on a separate plate. This identification plate is located under the device.


NOTE: The complete list of symbols and their description is provided in Chapter 1.5 on page 6.

2.2 Handpiece Identification Data

On the scaler handpiece are lasermarked the traceability data, UDI code included.


NOTE: The complete list of symbols and their description is provided in Chapter 1.5 on page 6.

2.3 Inserts Identification Data

Each insert is laser-marked with the traceability data.

Their packaging contains traceability data including the UDI code.

3 DELIVERY

3.1 List of Components

MICROPIEZO includes a basic equipment and a set of accessories that can be ordered separately, variable in relation to the configuration and customer requests (see Table on page 11).

NOTE: Both the items included in the standard supply and all accessories can be ordered separately by the client.

Refer to the MECTRON website, or your local retailer, for the list of available and compatible accessories

The packaging of the device cannot undergo strong impacts as contains electronic components, therefore the transport and the storage must be carried out with particular care.

Do not stack various boxes in order to avoid pinching the underlying packaging.

All the material sent by MECTRON was controlled at the time of dispatch.


The device is shipped appropriately protected and packaged.

Upon receipt of the device, check for any possible damage caused during transport and in case any damage and/or defects is found, complain to the transporter.

Keep the packaging in case any items need to be sent to an Authorised MECTRON Service Centre and to store the device during long periods of non-use.

⚠ WARNING: Before starting the treatment, always make sure to have a stock of material (handpiece, inserts, wrenches) to use in the event of failures or drawbacks.

EN

Basic equipment		
Ref	Item	Note
A	Device core unit ^{a)}	
B	Irrigation tank complete with cap ^{a)}	
C	Complete FS-11 Pedal ^{a)}	
D	Power supply ^{b)}	
E	Scaler Handpiece ^{a)}	MANIPOLO ABLATORE SLIM
F	Use and Maintenance Manual ^{a)}	
	Information: Online Documentation ^{a)}	

Accessories		
Ref	Item	Note
G	Torque wrench ^{a)}	K10 torque wrench
		K7 torque wrench ^{c)}
H	Inserts ^{a)}	
D1	Power adapter ^{b)}	UK plug
D2		AU plug
D3		EU plug
D4		CN plug
D5		US-JP plug
I	Tank cap ^{b)}	
J	Peristaltic pump ^{b)}	

a) Produced by Mectron.

b) Distributed by Mectron.

c) To be used only for the "Crown Preparation" kit.

4 INSTALLATION

4.1 First Installation

The device must be installed in a suitable place allowing its comfortable use.

⚠ WARNING: The place where the device is installed must meet the requirements found in the *Chapter 4.2 on page 12*.

MICROPIEZO can be purchased ready for use

or may need to be enabled by typing in an activation key.

In case your device requires an activation key, the procedures to be followed may vary from country to country.

Always refer to your dealer for details.

4.2 Safety Requirements in the Installation Phase

⚠ WARNING: Contraindications.

Interference with other equipment.

Even if compliant with standard IEC 60601-1-2, MICROPIEZO may interfere with other devices in the vicinity. MICROPIEZO must not be used near or stacked with other devices. However, if needed, it is necessary to check and monitor the correct operation of the apparatus in that condition.

⚠ WARNING: Contraindications.

Interference from other equipment.

An electrosurgical scalpel or other electrosurgical units positioned near the MICROPIEZO device may interfere with the correct operation of the device itself.

⚠ CAUTION: The electrical system of the premises in which the device is installed and used must comply with the rules in force and the relevant requirements of electrical safety.

⚠ CAUTION: Place the device in such a way that the power plug is easily accessible insofar as considered a disconnection device.

⚠ WARNING: Risk of explosion.

The device cannot operate in environments where there are saturated atmospheres of flammable gases (aesthetic mixtures, oxygen, etc.).

⚠ WARNING: Install the device in a safe place protected from impact or accidental water or liquid spray.

⚠ WARNING: Do not install the device above or near sources of heat. Arrange during installation for a suitable circulation of air around the device.

⚠ CAUTION: Do not expose the device to direct sunlight or sources of UV light.

⚠ CAUTION: The device can be transported but must be handled with care when moved. Place the pedal on the ground in such a way as to be activated only intentionally by the operator.

⚠ CAUTION: Before connecting the handpiece to its cord, check that the electrical contacts are perfectly dry, on both parts. If necessary dry them with compressed air.

⚠ CAUTION: The tank can hold liquids up to a maximum of 500 ml.


4.3 Connecting Accessories

The device is provided with a "wall" type external power supply unit and with interchangeable plug. Identify the appropriate type of plug and insert it into the power supply unit in such a way that the guides present on both parties are aligned. Slide the plug until you hear a "click".

⚠ CAUTION: If correctly assembled, the plug must not move or disconnect. If the plug is not properly secured, remove it and re-insert it. To remove the plug, press the hold button and slide the plug outwards.

⚠ CAUTION: Use only the power supply unit supplied with the device or an original part MECTRON.

1


2

In order to prevent permanent deformation of the water circuit or the manifestation of wear due to non-use, the system is supplied with the peristaltic pump not installed.


To install the peristaltic pump:

1. Remove the cover located on the back of the device;
2. Insert the peristaltic pump and make sure that the pin of the impeller and the clip (upper and lower) coincide with the respective seats.
3. Press up to the complete insertion of the clips.
4. Insert the water conduits in the relative sleeves for the whole length of the latter.
5. Close the cover.

Connect the connector of the pedal to the connection interface identified with the symbol , until you hear a 'click'. The connection interface is located in a recess formed in the bottom of the device.


3

Connect the connector of the power supply to the connection interface identified with the symbol  and located in a recess formed in the bottom of the device. Insert the connector up to the bottom.

Connect the power supply unit to the wall socket;

4

CAUTION: Position the device in a way that the power plug is easily accessible at all times; this is considered a means of isolation.

Remove the cap of the tank and fill it with the desired liquid;


5

CAUTION: The tank can hold liquids up to a maximum of 500 ml.

The tubing and all parts of the irrigation circuit are made from a material compatible with the main clinical solutions and liquids generally used for the treatments forming part of the intended use of the device, specifically:

- Water (distilled, demineralised, potable);
- Saline solution at 0.9%;
- Chlorexidine <0.3%;
- Ethyl Alcohol <15%;
- Citric Acid <5%.


CAUTION: To avoid chemical reactions that may damage the device, run a Long FLUSH cycle before using a different substance (see *Chapter 6 on page 29*).

Holding the tank in a vertical position, insert it on the machine core unit of the device until it comes into contact;

⚠ CAUTION: Do not overturn the tank. The plug is not watertight. The outflow of physiological solution or aggressive liquids may cause damage to the surfaces.


6

EN


Correctly insert the scaler handpiece on its cord by aligning the dot on the connector of the handpiece with the groove on the connector of the cord. Check that the electrical contacts of both are perfectly dry and if necessary dry them by blowing with compressed air;

7

8

Place the handpiece in the appropriate housing available on the device.

5 USE

5.1 Power On/Off

To turn on the system, after completing the procedure described in the *Chapter 4.3* on page 13, it will be necessary to press and hold the Power key for 2 seconds.

All the indicators on the keyboard will light up and the system, once the start-up procedure is complete, will emit an acoustic signal.

NOTE: During turning on, check that all the indicators illuminate correctly and that the system emits the acoustic signal. If any indicator does not light up or the system does not emit any sound, contact an Authorised MECTRON Service Centre.

At the first ignition and after an unexpected restart (for example, as a result of the interruption not provided of the power supply), the system will be configured with the following settings:

- power 0
- water 0
- **SOFT MODE** non active
- **LIGHT** active
- Wi-Fi active
- **FLUSH** non active

NOTE: For the description of the above settings refer to *Chapter 5.2* on page 17.

In all the other conditions, the system will be configured with the last settings used before proper shut-down.

⚠ CAUTION: With the pedal pressed down, the system will still boot, as provided in conditions of normal use, but will inhibit all the commands and configuration changes. This condition will be signalled by the flashing  symbol and 0 value on the irrigation scale. For information, refer to *Chapter 12.1* on page 69.

To turn the system off, press and hold again the Power key for 2 seconds.


Similarly to what occurs during start-up, the system will disable all the indicators and will emit an acoustic signal.

⚠ CAUTION: Place the device in such a way as to always have the power plug easily accessible as it is considered to be a disconnection device.

5.2 Description of the Keyboard

The system is provided with a "touch" keyboard.

The operator can configure the machine simply by touching the touch keyboard. When pressing each button the system will emit an acoustic signal as confirmation and/or execution of the command received. Depending on the option selected, the system will automatically adjust the correct working frequency.

EN

5.2.1 Functions and Power Levels of the Treatment

The power level of the treatment can be adjusted by selecting the numbers along the right arc .

It is possible to select different power levels for each type of treatment.

The system also indicates the selected treatment by the colour of the LED located at the base of the water tank.

For details regarding the available power levels, see the following table.

Treatment	Tank LED Colour	Power Level	Description
No treatment selected.	Blue	0	Null Intensity. Only irrigation will flow from the insert on the handpiece (if set with a value greater than 0).
endo	Green	1	Dedicated to endodontic treatments such as cleansing of root canals and retrograde approach.
		2	
perio - scaling	Yellow	1	Dedicated to all the procedures of supra- and sub-gingival prophylaxis, of root planing and cleaning of the implant surface.
		2	
		3	
prosthetics	Purple	.	Dedicated to conservative techniques and prosthesis.

5.2.1.1 Soft Mode

For patients with particularly sensitive teeth, the distributed power can be automatically regulated by activating the **SOFT MODE** function .

In **SOFT MODE**, the system will not attempt to ensure maximum performance during the treatment, but rather adjust and limit the power based on the opposing resistance on the insert, adopting a less invasive reaction.

The **SOFT MODE** can be enabled for all types of treatment and power levels with the exception of power level 0.

Furthermore, by using **SOFT MODE**, the operator has double the number of selectable intensity/power levels for each type of treatment.

The following table shows the available power levels in **SOFT MODE**.

Treatment	Power Level	Soft Mode
-	0	-
endo	1	Active
	1	Non Active
	2	Active
	2	Non Active
perio - scaling	1	Active
	1	Non Active
	2	Active
	2	Non Active
	3	Active
	3	Non Active
prosthetics	.	Active
	.	Non Active

Once selected, **SOFT MODE** remains active even if the power levels of a specific treatment are modified, for example, when switching from **perio-scaling 1** to **perio-scaling 2**. The **SOFT MODE** function is deactivated on

the other hand when switching from a power level of a different treatment, for example, when switching from **perio-scaling 2** to **endo 1**.

5.2.2 Water Delivery

The flow rate of the peristaltic pump can be adjusted by selecting the numbers along the left arch .

There are 7 flow rate levels:

- 0 : the operation of the pump is off: no irrigation; flows from the insert;
- from 1 to 6 : the pump flow rate goes from 10 ml/min to 33 ml/min $\pm 10\%$.

 **WARNING:** The water flow level 0 can be selected for all treatments and all power levels. For the safety of the patient, it is not recommended to disable the delivery of the water during the treatments: **perio-scaling** (with any power) e **prosthetics** with inserts that are not specifically intended for use without irrigation.

NOTE: The minimum flow rate of the peristaltic pump can vary from 6 to 12 ml/min depending on the time of use and state of wear.

5.2.3 Filling and Cleaning the Water Circuit

The **FLUSH** button allows a cleaning cycle to be run on the irrigation circuit and has two operating modes: Short and Long.

- **Short FLUSH:** the SHORT FLUSH cycle allows to get the irrigation at the tip before starting the treatment;
- **Long FLUSH:** The LONG FLUSH cycle allows a cleaning cycle to be run on

the irrigation circuit at the end of the treatment and before proceeding with the cleaning and sterilisation of the parts used.

5.2.3.1 FLUSH - Short Mode

1. Press the button **FLUSH** once to start the short mode.
The water flow rate control will illuminate and the **Flush** button will flash.
2. Press and release the pedal to start the procedure.
The system will activate the water pump with flow rate 6 for 10 seconds.
The numbers on the water scale will progressively turn off (from 6 to 0) indicating the progress of the procedure.
The **FLUSH** button will remain lit for the whole duration of the procedure (10 seconds).

5.2.3.2 FLUSH - Long Mode

1. Press the button **FLUSH** twice to start the long mode.
The water and power flow rate control will illuminate.
The button **FLUSH** will flash.
2. Press and release the pedal to start the procedure.
The system will activate the water pump with flow rate 6 for 20 seconds.
The numbers on the water scale and those on the power scale will progressively turn off indicating the progress of the procedure: first the numbers from 6 to 0 on the water scale, and then values from prosthetics to 0 on the power scale.
The **FLUSH** button will remain lit for the whole duration of the procedure (20 seconds).

5.2.3.3 FLUSH - Interruption

The cycle can be interrupted as soon as the liquid exits from the handpiece and with the following modes:

1. Press the button **FLUSH** again or, alternatively, press the pedal again or wait until the end of the timer.
The procedure will be interrupted, the keyboard will be active again and will show the last configuration used.

5.2.4 Light

The **LIGHT** button is used to deactivate and activate the light on the handpiece tip.

The **LIGHT** function is active by default.

The circular LED light lights up when pressing the pedal and stays on for the next 3 seconds after releasing the pedal.

If the **LIGHT** button is pressed, the light on the handpiece cone stays off even when the pedal is pressed.

5.2.5 Wi-Fi - IoT Technology

MICROPIEZO is designed for the use of IoT (Internet of Things) technology via the Mectron IoT platform, in the countries where this technology is active. The platform is designed to create a network of Mectron devices connected to the Internet and equipped with software that allows you to exchange data with Mectron remote servers. The use of said platform allows the device to send diagnostic information of the device itself with the aim of allowing Mectron to improve its usability, be more effective in maintenance and to improve its durability. The platform allows the updating of the device and to receive device firmware improvements.

In order to take advantage of the IoT technology, MICROPIEZO must be connected to a Wi-Fi network with Internet access and the user must use the official mobile app "My Mectron" to manage its connectivity and the authorization to exchange data (protected and private) with the Mectron cloud.

All data sent to the private MECTRON server are made available to users via the mobile App. Each data access is protected by authentication and authorization mechanism. The data collected through IoT technology will be managed, stored, and used by MECTRON with maximum transparency for purposes associated with the improved performance of MICROPIEZO and any other

investigation into methods of using the device.

MECTRON reserves the right not to activate, or, if already activated, to interrupt the use of MICROPIEZO IoT technology.

NOTE: For specific support regarding IoT (Internet of Things) functions, contact an Authorised Mectron Service Centre.

NOTE: Activation of the IoT function is subject to the availability of the My Mectron mobile application (APP) in the country of sale of the device and the reachability of the Mectron server.

It is possible to download the App at this link: <https://app.mectron.com> or by scanning the QR-Code:

If the Mectron Server cannot be accessed by the device, the IoT function will be automatically disabled. without any impact on the operation of the device. **In the absence of communication with the Mectron Server, the functions managed by IoT technology will not be available.**

5.2.5.1 Connection to a Wi-Fi Network

The **WiFi** button flashes only if the device's Wi-Fi interface is active and when searching for a valid connection with the Mectron server. Once a connection has been established with the server, the **WiFi** button will stay lit (not flashing).

NOTE: During normal operation, the Wi-Fi button turns off when the pedal is pressed. If the Wi-Fi button was lit, it will light up again when the pedal is released.

If the button is pressed, the connection will be deactivated; to reactivate it, press the button again.

5.2.5.2 Wi-Fi Network Requirements

Check that the Wi-Fi network used to configure MICROPIEZO satisfies the following requirements:

1. WPA safety protocol enabled - for network protection;
2. DHCP enabled - for automatic assignment of a local network address;
3. If there is a firewall, check that it allows outbound connections from the configured local network, at least on ports 123, 443 e 8883.

If requirements 2 and 3 are not satisfied, the

5.2.5.3 Software Update

Using the My Mectron APP it is possible to check for and eventually update the software installed in the device.

NOTE: The My Mectron APP has a guided procedure and user information. Refer to the My Mectron APP for details regarding the supported functions.

If updated software is available for the device, the user will receive notification via the My Mectron APP.

It will be the user's choice whether or not to update the device.

The device has software protection mechanisms allowing to complete the update if the process is accidentally interrupted.

An auto-repair mechanism is provided which allows you to reactivate the latest version installed if the machine's flash memory

IoT functions may be inhibited due to a lack of communication with the Mectron server.

If the Wi-Fi network used does not satisfy safety requirement 1, the network is not accepted nor presented among the usable networks.

NOTE: Failure of the IoT technology does not affect the intended functionality of the device in any way whatsoever. The main functions of the device described in this manual do not depend on IoT technology.

becomes corrupt.

If the user runs a software update, at the end of the process the device must be switched off and on again to complete the procedure.

The last phase of the update will be managed with the restart: the procedure requires that during this phase the LED of the wrench remains active for approximately two minutes. At the end of the procedure, the device will be turned off and a courtesy push notification will be sent to the mobile app.

NOTE: It is important that the user does not turn off the device during this phase to avoid triggering a reset procedure.

NOTE: The device will not be rebooted and will not undergo any updates or changes during operation.

5.2.6 Diagnostic System

MICROPIEZO is provided with a diagnostic circuit that allows the detection of malfunctions and viewing of their type on the touch keyboard by means of a symbol.

To help the user identifying the part that is non-functioning there are four symbols described in *Chapter 12.1 on page 69*.

5.3 Safety Requirements Before and During Use

⚠ WARNING: Before starting the treatment, always make sure to have a stock of material (handpiece, inserts, wrenches) to use in the event of failures or drawbacks.

⚠ WARNING: Exclusively use original MECTRON inserts, accessories and spare parts.

⚠ WARNING: Use of non-original inserts MECTRON: this implies a permanent damage to the thread of the handpiece and compromise the correct operation with the risk of causing harm to the patient.

⚠ CAUTION: Contraindications - Ultrasonic Scaler. Do not perform treatments on metal or ceramic prosthetics. Ultrasonic vibrations could lead to decementation of the dentures.

⚠ WARNING: Contraindications. Do not use MICROPIEZO on patients with Pacemakers or other implantable electronic devices. This requirement is applicable also for the Operator.

⚠ WARNING: Do not perform scaling treatments without irrigation to avoid an overheating of the insert that can cause damage to the tooth.
The treatments provided without irrigation can be exclusively those performed with "Dry Work" inserts, without passage of water.

⚠ CAUTION: In the treatments that require irrigation use only inserts with passage of liquids.

⚠ WARNING: Treatments that require irrigation. Always check the operation of the irrigation before and during use. Make sure that the liquid exits from the insert. Do not use the device if the irrigation does not work or if the pump is faulty.

⚠ WARNING: Checking the condition of the device before treatment. Always check that there is no water under the device. Before every treatment always check the perfect operation of the device and efficiency of the accessories. If operation anomalies are observed do not carry out the treatment. If the faults concern the device, contact an Authorised MECTRON Service Centre.

⚠ WARNING: Do not use the device in case of breaks or cracks on the casing.

⚠ WARNING: Infection control. First use: All parts and reusable accessories (new or returning from an Authorised MECTRON Service Centre) are delivered in non-sterile conditions and must be treated before each use, following the instructions in *Chapter 8 on page 33*.

Subsequent uses: After each treatment, clean and sterilise all parts and reusable accessories according to the instructions in *Chapter 8 on page 33*.

⚠ WARNING: Infection control. Do not leave liquids in the tank for long periods of time. The tank must be filled after having cleaned and sterilised all parts, and before a treatment. If the tank has been filled without having used the device, it must be emptied at the end of the day, and all parts and accessories must be cleaned and sterilised.

⚠ WARNING: To ensure that the handpiece cools down activate it always with irrigation circuit installed and properly filled. To fill the irrigation circuit always use the FLUSH function.

⚠ CAUTION: For correct use of the device, the pedal must be pressed and the device started with the insert not in contact with the part to be treated, so that the electronic circuit is able to recognise the best point of resonance of the insert without interference, thus allowing optimum performance.

⚠ WARNING: Before every treatment make sure that on the handpiece is inserted the insert suitable for the treatment. Only use the torque wrench MECTRON to fix the insert to the handpiece.

⚠ WARNING: The patient must not come into contact with the machine core unit or with the pedal.

⚠ WARNING: To avoid causing injury to the Operator, do not replace the insert while the handpiece is in function.

⚠ CAUTION: LONG FLUSH function. After each use it is necessary to run a cleaning cycle on the irrigation circuit with the **LONG FLUSH function** (*Chapter 6 on page 29*). If a cleaning cycle of the irrigation circuit is not carried out, salts crystallization could damage the device.

⚠ CAUTION: FLUSH function. The **FLUSH** function must be used after each treatment, before starting the cleaning and sterilisation procedures.

⚠ WARNING: During the intervention on the patient, do not perform any maintenance tasks on the system.

⚠ WARNING: Insert breakage and wear. High frequency oscillations and wear can, in rare cases, lead to the breakage of the insert. Do not bend, change shape or sharpen an insert in no way.

Folding an insert or levering on it can cause the insert to break.

Deformed or otherwise damaged inserts, are susceptible to breakage during use. These inserts should never be used.

Excessive pressure on the inserts during use can lead to breakage.

In case of breakage check that no scraps remain in the treated part and at the same time effectively suck in order to remove them. It is necessary to instruct the patient to breathe through the nose during the treatment, or use a dental dam, in order to avoid ingestions of fragments of broken inserts.

Check the state of wear of the insert and its integrity before and during each use.

If a drop in performance occurs, see to its replacement.

The state of wear of the most common inserts (S1, S1-S, S2, S5, P2, P4, P10) can be checked using the supplied INSERT-CARD. To use the INSERT-CARD correctly:

- Position the insert on the INSERT-CARD so that the profile corresponds to the one printed on the card. The printed profile on the card has a red line indicating the limit of wear;
- If the insert is shorter than the limit of wear, its performance will be significantly inferior compared to that of a new insert, and should therefore be replaced.

If the layer of titanium nitride (gold-plated surface), where present, is visibly worn, the insert must be replaced. The use of a worn insert reduces its efficiency.

Diamond Inserts: diamond inserts must be replaced when the layer of titanium nitride is visibly worn and in any case after a maximum of 10 treatments.

When the nitriding is consumed, the cutting edge loses effectiveness; a possible regrinding harms the insert, therefore it is prohibited. Verify that the insert is not worn. During the operation frequently check that the insert is intact, especially in the apical part.

During the operation avoid prolonged contact with the retractor or with metal instruments in use.

 **CAUTION: Contraindications.** After autoclaving the handpiece, the inserts, the torque wrench, or any other sterilisable accessory, wait until it has completely cooled before reuse.

 **CAUTION: The electrical contacts inside the connectors on the handpiece and cord must be dry.** Before connecting the handpiece to its cord, make sure that the electrical contacts of the connector on both sides are perfectly dry, especially after the sterilisation cycle in an autoclave. If necessary dry the contacts by blowing them with compressed air.


 **CAUTION:** For its conformation, the handpiece can roll. The handpiece when not used must always be placed on its support.

5.4 Instructions for Use

After having connected all the accessories as shown in *Chapter 4.3* on page 13 proceed as follows:


Remove the handpiece from its seat.

1

To load the irrigation circuit use the function **SHORT FLUSH** by pressing once the button **FLUSH** on the keyboard. The water flow rate control will illuminate and the **FLUSH** button will flash;

2


NOTE: With handpiece put back, the **FLUSH** cycle will not start. For information, refer to *Chapter 12.1* on page 69.

Press the pedal and release it to start the procedure.

The system will activate the water pump. The numbers on the water scale will progressively turn off (from 6 to 0) indicating the progress of the procedure. The **FLUSH** button will remain lit for the duration of the procedure;

3

NOTE: Putting the handpiece back, the **FLUSH** cycle will be interrupted and the keyboard will be active again by returning the last configuration used.

The cycle can be interrupted as soon as the liquid exits the handpiece by pressing again the FLUSH button or, alternatively, when the pedal is pressed again.

Both in the case where the procedure is interrupted and when the cycle successfully concluded, the keyboard will be active again by returning the last configuration used.

Screw in the preselected insert on the MICROPIEZO handpiece until it comes into contact;


Insert the insert inside the torque wrench MECTRON.

Firmly grip the body of the handpiece;


⚠ CAUTION: Do not grab the handpiece by the nozzle and/or by the cord. Firmly hold the body of the handpiece and rotate only the torque wrench. Do not rotate the body during tightening.

Turn the torque wrench in a clockwise direction until the notch clicks (the external body of the torque wrench turns with respect to the body of the handpiece, emitting mechanical sounds). The insert is now optimally locked;


4


5

6


7


Select the type of function, power and irrigation necessary and the light if desired, on the keyboard;

8

At the end of treatment stow the scaler handpiece on its housing.

9

5.5 Important Information on the Inserts

WARNING:

- Before using the sterilised insert, check the integrity of the sterile packaging and inspect the product to ensure there is no damage. The insert will no longer be sterile if the packaging is broken or damaged. If the packaging is damaged, the insert **MUST** be sterilised again before use.
- Before starting to work, properly tighten the insert to the handpiece using the torque wrench.
- When the layer of titanium nitride is visibly consumed the insert must be replaced. The use of an insert too worn decreases its efficiency.
- **Diamond Inserts:** diamond inserts must be replaced when the layer of titanium nitride is visibly worn and in any case after a maximum of 10 treatments.
- Do not activate the handpiece while the insert is in contact with the part to be treated so that the electronic circuit can recognise the best point of resonance of the insert allowing optimum performance.
- Check the state of wear of the insert and its integrity before and during each use. If you experience a loss of performance provide for its replacement.
- Use only original MECTRON inserts. The use of non-original inserts will void the warranty and damage the thread of the MICROPIEZO handpiece, with the risk of no longer being able to correctly screw in the original inserts for subsequent use. In addition, machine settings are tested and guaranteed for a correct operation using only original inserts MECTRON.
- Do not change in any way the shape of the insert, by bending it or filing it. This could cause it to break.
- Do not use an insert that has undergone deformation of any type.
- Do not attempt to sharpen an insert used.
- Always verify that the threaded parts of the insert and the handpiece are perfectly clean - see *Chapter 8 on page 33*.
- Excessive pressure applied to the insert can cause it to break and potentially injure the patient.
- The Mectron inserts vibrate with longitudinal oscillation, with a back and forward movement. During the treatment, always keep the instrument in a tangential direction with respect to the surface of the tooth. Move the handpiece back and forward, applying slight side pressure.
- Do not point the instrument directly onto the surface of the enamel or implant. Position the tip/operating part only in a tangential direction with respect to the surface of the tooth or implant.
- The insert must be held in constant movement. If the insert is blocked it may cause overheating of the treated part. It is recommended to use continuous motions to minimise contact between the tip and the part. Do not block them against the tissues so as not to cause their overheating. It is advisable to increase the level of irrigation as the power level increases.
- Leave the ultrasonic vibrations to work, do not exert excessive pressure on the inserts during use. Apply slight force on the insert for optimum performance.
- When the insert is used in interproximal spaces, do not block the instrument or leverage on the operating part. The inserts must be left free to vibrate.

- In endodontic root canal therapy, never operate the files outside of the root canal to prevent them from breaking. To prevent breakages, create a smooth path with a manual endo file and plan the most direct possible path of access to minimise any bending of the insert.


Make delicate movements. Examine the file often for signs of wear. If the file breaks inside the canal, avoid contact between the instrument and the broken file to stop it from being pushed further in. Do not apply pressure on the insert in an axial direction.

6 LONG FLUSH FUNCTION

The Long **FLUSH** function allows a cleaning cycle to be run on the irrigation circuit.


Fill the tank with water, preferably demineralised water. Connect the tank to the device;

1

Place the handpiece, with or without insert, above a container for the liquid that will be released during the cleaning cycle;


2

NOTE: With handpiece put back, the **FLUSH** cycle will not start. For information, refer to *Chapter 12.1* on page 69.

To enter the cleaning mode (LONG **FLUSH**) press the button **FLUSH** on the touch keyboard twice. The water and power flow rate control will illuminate. The button **FLUSH** will flash;

3


NOTE: It is possible to exit from the cleaning mode at any time by pressing again the button **FLUSH**; the keyboard will be active again by returning the last configuration used;

To start the **LONG FLUSH** cycle, press the pedal and release while the **FLUSH** button flashes;

The cycle can be stopped at any time by pressing again the button **FLUSH** or, alternatively, when the pedal is pressed again.

The keyboard will be active again by returning the last configuration used.

4


5

The system will activate the water pump with flow rate 6 for 20 seconds.

The numbers on the water scale and those on the power scale will progressively turn off in succession indicating the progress of the procedure: first the numbers from 6 to 0 on the water scale, and then the values from prosthetics to 0 on the power scale.

The **FLUSH** button will remain lit for the whole duration of the procedure (20 seconds).

Once the cycle is complete, the keyboard will be active by setting the last configuration used.

After cleaning the irrigation circuit, disassemble the individual parts (see *Chapter 7 on page 31*) and proceed to clean and sterilise them (see *Chapter 8 on page 33*).


7 DISASSEMBLING PARTS FOR CLEANING AND STERILISATION

Before carrying out the cleaning procedures described in *Chapter 8 on page 33*, disconnect all accessories and components of the MICROPIEZO.

⚠ WARNING: Switch the device off.


Before proceeding with cleaning and sterilisation operations, if the device is switched on, turn it off by pressing and holding the Power button for 2 seconds, disconnect the adapter from the wall socket and disconnect the power cable from the machine core unit.

1


2

If present, remove the insert from the handpiece using the torque wrench;


3

Disconnect the handpiece from the relative cord;

4


⚠ CAUTION: Do not try to unscrew or turn the connector when disconnecting the handpiece.
The connector can be damaged.

Unscrew the cone from the handpiece;


⚠ CAUTION: The cone has a light pipe. If the cone is unscrewed, the light pipe will no longer be held in place and may slide and be disconnected. Be careful not to lose the light pipe.

5


Remove the light pipe from the cone;

6


Disconnect the tank from the machine core unit, lifting it upwards;

7

Disconnect the pedal from the device: hold the connector of the pedal, press the release tab and pull the connector;

8

⚠ CAUTION: Do not try to unscrew or turn the connector during disconnection: this might damage the connector.

⚠ CAUTION: While disconnecting the cable from the pedal, keep always and only the connector of the cord. Never pull the cord itself.

8 CLEANING AND STERILISATION

⚠ CAUTION: To disinfect the device and/or its accessories, it is recommended the use of water-based disinfectant solutions with a neutral pH (pH7).

DO NOT USE as disinfectants:

- Alcohol-based disinfectant solutions;
- Peracetic acid, Formaldehyde, Gluteraldehyde based products or other equivalent solutions/products;
- Very alkaline products (pH > 9);
- Products containing sodium hypochlorite;
- Products containing hydrogen peroxide;
- Products containing abrasive substances;
- Very acid products (PH < 4);
- Products containing aldehyde, amine and/or phenols

- Acetone;
- Methyleneethylketone;

since they may discolour and/or damage the materials of the device and its accessories. The manufacturer disclaims all liability for any damage caused by the substances mentioned above. In case of damage caused by those substances, the Warranty will be void.

The following tables represent only mere indications of the cleaning and sterilisation activity that must be ensured by the user. For the full procedures for cleaning and sterilisation of individual parts refer to the Chapters referenced in the tables.

⚠ CAUTION: It is prohibited to use methods not covered in the tables below.

Machine Core Unit, Pedal, End Part of Cord		
Phase	Chapter	Procedure
I	8.1	Preparation
II	8.2	Manual cleaning with enzymatic detergent

Tank and Cap		
Phase	Chapter	Procedure
III	8.3	Manual cleaning with enzymatic detergent and running water
IV	8.4	Disinfecting the Tank Irrigation Circuit

Accessories (Scaler Handpiece, Scaler Cone, Inserts, Torque Wrench)							
Phase	Chapter	Procedure	Scaler Handpiece	Scaler Cone	Inserts	Torque Wrench	
V	8.5.1	Manual cleaning	Manual cleaning with Enzymatic Detergent		Immersion into enzymatic detergent		
					Ultrasonic Cleaning		
	8.5.2	Automatic cleaning	Thermal-disinfector with detergent				
VI	8.6	Check cleaning	Visual inspection with Lens 2.5x				
VII	8.7	Drying and lubrication	Drying			Drying and Lubrication	
VIII	8.8	Sterilisation	Packaging, Sterilisation and Storage				

NOTE: Repeated reconditioning has a minimal effect on the devices and their accessories. The end of the service life of the accessories is generally determined by wear or damage resulting from use. Mectron guarantees the integrity of its sterilisable scaler handpieces for up to 250 reconditioning cycles.

8.1 Preparation

1. Perform the **FLUSH** function (see *Chapter 6 on page 29*);
2. Check that all of the following accessories have been removed/disconnected from the machine core unit (see *Chapter 7 on page 31*):
 - Power supply cable;
 - Pedal;
 - Scaler handpiece;
 - Inserts.

 **CAUTION:** The cleaning and sterilisation operations described in the following chapters must be carried out at the time of first use and at all subsequent uses.

 **WARNING:** Always switch off the device by pressing and holding the Power button for 2 seconds and disconnect it from the mains before carrying out any cleaning.

 **CAUTION:** Always disconnect the insert from the handpiece before cleaning and sterilising it.

 **CAUTION:** Do not immerse the handpiece in disinfecting solutions or other liquids as these could damage it.

 **CAUTION:** Do not immerse the handpiece in the ultrasonic tank.

8.2 Manual Cleaning with Enzymatic Detergent

The following procedure must be performed on all **non-sterilisable** parts of the device, except for the tank and its cap. The parts in question are:

- Machine core unit;

- Pedal and relative connection cable to the machine;
- End part of the cord.

» REQUIRED MATERIALS


- Clean, soft, low-lint cloths;
- Cleaning solution (pH 6-9).

» CLEANING METHOD


Clean the surface of the parts in question using a clean, dry, low-lint cloth, dampened with a detergent solution (pH 6-9);

Dry the parts using a dry, non-abrasive, low-lint cloth.

1

2

⚠ CAUTION: Do not sterilise the parts in question. They may stop working and cause damage to people and/or property.

⚠ WARNING: Always switch off the device by pressing and holding the Power button for 2 seconds and disconnect it from the mains before carrying out any cleaning.

⚠ WARNING: The device and its parts that cannot be sterilised are not protected against the penetration of liquids. Do not spray liquids directly on the surface of the device or its parts that cannot be sterilised.

⚠ CAUTION: Do not use running water to clean the parts in question.

⚠ CAUTION: Do not soak these parts in liquid and/or various kinds of solutions.

8.3 Manual Cleaning with Enzymatic Detergent and Running Water

The following procedure must be carried out **on the tank and cap** of the device.

» PREPARATION

1. Disconnect the tank from the machine core unit (see *Chapter 7 on page 31*);
2. Unscrew the cap from the tank.

⚠ CAUTION: Do not sterilise the tank and cap in an autoclave. They could be damaged.

» REQUIRED MATERIALS

- Water;


- Cleaning solution (pH 6-9);
- Clean, soft, low-lint cloth;
- Demineralised water.

NOTE: Do not use water with a temperature higher than 60° C.

» CLEANING METHOD


Thoroughly rinse the inside and outside of both the tank and cap under running water;

1


Clean the external and internal surfaces of the tank and cap with a clean, soft, low-lint cloth, dampened with a mild detergent solution (pH 6-9);

2


Thoroughly rinse the inside and outside of the tank and cap under running water to eliminate all residue of the detergent solution;

3


Remove any residue on the external or internal surfaces of the tank and cap using a soft, low-lint cloth, dampened with demineralised water;

4

When cleaning operations are over, check the pieces under an appropriate light source, paying attention to any soil residue and, if necessary, repeat the cleaning cycle.

5

EN

8.4 Disinfecting the Tank Irrigation Circuit

After having manually cleaned the tank and cap with enzymatic detergent and running water (see *Chapter 8.3 on page 36*) and before carrying out the subsequent cleaning procedures, proceed to disinfect the tank irrigation circuit.

To disinfect the irrigation circuit, proceed as follows:


» REQUIRED MATERIALS

- CIDEX® OPA disinfectant solution;
- Distilled water;
- Clean, soft, low-lint cloth.

» PROCEDURE


Disconnect the tank from the machine core unit, lifting it upwards;

1


Lift the cap off the irrigation solution tank and empty it;
Rinse the irrigation solution tank with running water;

2

Fill the irrigation solution tank with 250ml of disinfectant (CIDEX® OPA);


3

Holding the tank in a vertical position, insert it on the machine core unit of the device until it comes into contact;

CAUTION: Do not overturn the tank. The plug is not watertight. The outflow of physiological solution or aggressive liquids may cause damage to the surfaces.


4

EN


Remove the handpiece connector from its seat.

5


Position the scaler handpiece connector above a container or sink to catch the liquid released during the disinfection cycle;

6

Press the **FLUSH** button on the touch keyboard twice (LONG FLUSH - see Chapter 6 on page 29). The water and power flow rate control will illuminate. The button **FLUSH** will flash;

7


NOTE: It is possible to exit from the cleaning mode at any time by pressing again the button **FLUSH**; the keyboard will be active again by returning the last configuration used;

To start the **LONG FLUSH** cycle, press the pedal and release while the **FLUSH** button flashes;

The cycle can be stopped at any time by pressing again the button **FLUSH** or, alternatively, when the pedal is pressed again.

The keyboard will be active again by returning the last configuration used.

8


The system will activate the water pump with flow rate 6 for 20 seconds.

The numbers on the water scale and those on the power scale will progressively turn off in succession indicating the progress of the procedure: first the numbers from 6 to 0 on the water scale, and then the values from prosthetics to 0 on the power scale. The **FLUSH** button will remain lit for the whole duration of the procedure (20 seconds).

Once the cycle is complete, the keyboard will be active by setting the last configuration used.

9


Repeat the **LONG FLUSH** cycle (see points 7 to 9 of the procedure) a second time;

10

At the end of the **LONG FLUSH** cycles, wait 12 minutes;

11


Repeat the **LONG FLUSH** cycle (see points 7 to 9 of the procedure) 3 times;

12

At the end of the LONG FLUSH cycles, wait 12 minutes;

13


EN

Repeat the LONG FLUSH cycle (see points 7 to 9 of the procedure) 3 times;

14


At the end of the LONG FLUSH cycles, wait 12 minutes;

15


Disconnect the irrigation solution tank from the device (see point 1);
Lift the cap off the irrigation solution tank and empty it (see point 2);
Rinse the irrigation solution tank with running water;

16

With the tank disconnected, run a LONG FLUSH cycle (see points 7 to 9) to empty the irrigation circuit of the liquid inside;


17

Fill the irrigation solution tank with 300ml of demineralised water;

NOTE: Use of distilled water is recommended for long periods of non-use. If the device is used after the disinfection procedure, use common running water.


18

Holding the tank in a vertical position, insert it on the machine core unit of the device until it comes into contact;

CAUTION: Do not overturn the tank. The plug is not watertight. The outflow of physiological solution or aggressive liquids may cause damage to the surfaces.

19


On the touch keyboard (see Chapter 5 on page 16), select:

- Power Level: 0;
- Irrigation Level: 6.

Position the scaler handpiece connector above a container or sink to catch the liquid released during the disinfection cycle (see point 6 of the procedure); Press the pedal to start the selected program until the tank is completely drained of all the water;


20

Disconnect the irrigation solution tank from the device (see point 1) and lift the cap;

Dry the valve area and tank with a soft, low-lint cloth.

21

8.5 Cleaning the Sterilisable Accessories

The sterilisable parts of the device are:

- Scaler handpiece;
- Scaler cone;
- Inserts;
- Inserts tightening wrench.

Before proceeding with operations to check cleanliness (*Chapter 8.6 on page 53*), drying and lubrication (*Chapter 8.7 on page 54*) and then sterilisation (*Chapter 8.8 on page 55*), one of the three possible cleaning methods explained in detail in the following sub-chapters must be selected.

⚠ CAUTION: The following instructions have been validated by MECTRON in order to enable the medical device to be reused. The person responsible for the maintenance, cleaning and sterilisation process, must ensure that repeated processes are effectively carried out using the equipment, materials and staff in the reprocessing structure, in order to achieve the desired result. This generally requires the validation and systematic monitoring of the process. Similarly, all deviations from the instructions provided by the person responsible for the processes should be properly evaluated for judging the effectiveness and the potential undesirable consequences.

8.5.1 Manual Cleaning

Manual cleaning can be carried out as an alternative to the automatic cleaning as described in *Chapter 8.5.2 on page 51*.

» REQUIRED MATERIALS


- Enzymatic detergent at pH 6-9;
- Water;
- Container for immersion in the enzymatic liquid;
- Ultrasonic tank;
- Clean, soft, low-lint cloths;
- Brush with soft nylon bristles;
- Syringe;
- Demineralised water.

Prepare a solution of enzymatic detergent^{a)} with pH 6-9, according to the instructions of the manufacturer;

⚠ CAUTION: Once used, dispose of the enzymatic detergent correctly, do not recycle it.

- a) Process validated by independent bodies with Mectron ENZYMEC enzymatic detergent, 0.8% v/v.


1

8.5.1.1 Scaler Handpiece

Completely immerse the scaler handpiece, cone, and light pipe in the enzymatic solution. Leave to soak for 10 minutes at $22^{\circ}\text{C} \pm 2^{\circ}\text{C}$;




2

3


Gently brush the entire surface of the scaler handpiece, cone, and light pipe for at least 20 seconds using a brush with soft nylon bristles, paying particular attention to the following areas:

- thread of the scaler handpiece;
- titanium shaft;
- external and internal parts of cone;
- external and internal parts of light pipe.


Rinse the inner channel of the scaler handpiece using a 20 ml syringe previously filled with a new enzymatic solution. Repeat twice;

4


Remove the handpiece, cone, and light pipe from the enzymatic solution and gently brush their surfaces using a brush with soft nylon bristles under hot ($40^{\circ}\text{C} \pm 5^{\circ}\text{C}$) running water, for at least 10 minutes;

5


Fill the ultrasonic tank with the enzymatic detergent solution prepared according to the manufacturer's instructions. Use water at ambient temperature ($22^{\circ}\text{C} \pm 2^{\circ}\text{C}$);

6


Place the scaler handpiece, cone, and light pipe in the ultrasonic tank immersed in the enzymatic detergent solution at $22^{\circ}\text{C} \pm 2^{\circ}\text{C}$ and run a cycle for at least 10 minutes;

7


8

Remove the handpiece, cone, and light pipe from the enzymatic solution and gently brush their surfaces using a brush with soft nylon bristles under hot ($40^{\circ}\text{C} \pm 5^{\circ}\text{C}$) running water, for at least 10 minutes;


9

Place the scaler handpiece, cone, and light pipe in the ultrasonic tank immersed in demineralised water at $40^{\circ}\text{C} \pm 5^{\circ}\text{C}$ and run a cycle for at least 10 minutes;


10

Rinse the inner channel of the scaler handpiece using a 20 ml syringe previously filled with demineralised water. Repeat five times;

11

Dry the surface of the scaler handpiece, cone, and light pipe with a clean, non-abrasive, low-lint cloth.

8.5.1.2 Inserts

Place the insert in a clean container, in a horizontal position. Add enough enzymatic detergent solution to completely cover it.

Leave the insert to soak in the enzymatic detergent solution for 10 minutes at ambient temperature ($22^{\circ}\text{C} \pm 2^{\circ}\text{C}$);

While soaking in the enzymatic solution, gently brush all the surfaces until all visible soil has been removed.

Use a clean brush with soft nylon bristles for the external surfaces, and a clean pipe cleaner with soft nylon bristles for the internal cavities and grooves.


NOTE: Thoroughly brush all the following parts of the insert for about 20 seconds:

- Through-holes and internal channels;
- Hard-to-clean areas such as sharp edges and in particular, the interstices between the cutting cusps;
- Internal cavities, grooves and fissures.


Remove the insert from the enzymatic detergent solution.

Thoroughly rinse and brush all the surfaces of the insert (see previous point) under running water for at least 10 minutes;


12

13


14

EN


Use a disposable syringe to withdraw and inject the enzymatic cleaning solution into the hard-to-reach areas (through-holes/cannulae). Repeat this operation three times to ensure the effective removal of the soil from the internal surfaces of the through-hole,;

15


Rinse the internal channel of the insert with demineralised water at a pressure of 3.8 bar for at least 10 seconds so as to eliminate any soil residue,;

16


Place the insert in the ultrasonic tank immersed in the enzymatic detergent solution at 24 °C $\pm$ 2 °C and run a cycle for at least 20 minutes;

17


Remove the insert(s) from the ultrasonic tank and rinse them in demineralised water;

18


Brush the internal and external surfaces of the insert with a clean brush with soft nylon bristles, under running water;

19


Rinse the insert(s) in demineralised water;

20

Rinse the internal channel of the insert with demineralised water at a pressure of 3.8 bar for at least 10 seconds so as to eliminate any soil residue.


21

8.5.1.3 Inserts Tightening Wrench

Place the torque wrench in a clean container, in a horizontal position. Add enough enzymatic detergent solution to completely cover the torque wrench. Leave the torque wrench to soak in the enzymatic detergent solution for 10 minutes at ambient temperature (22 °C ± 2 °C). This procedure reduces the amount of blood, protein, and mucus on the instrument(s);

22

While soaking in the enzymatic solution, gently brush all the surfaces until all visible soil has been removed.

Use a clean brush with soft nylon bristles for the external surfaces, and a clean pipe cleaner with soft nylon bristles for the internal cavities and grooves.

NOTE: Thoroughly brush all the following parts of the torque wrench for about 20 seconds:


- Through-holes and internal channels;
- Hard-to-clean areas such as sharp edges and in particular, the interstices between the cutting cusps;
- External metal ring;
- Internal cavities, grooves and fissures.

Remove the torque wrench from the enzymatic detergent solution.


Thoroughly rinse and brush all the surfaces of the torque wrench (see previous point) under running water for at least 10 minutes;

Place the torque wrench in the ultrasonic tank immersed in the enzymatic detergent solution at $24^{\circ}\text{C} \pm 2^{\circ}\text{C}$ and run a cycle for at least 20 minutes;


23

24


25

Remove the torque wrench from the ultrasonic tank and rinse it under running water so as to eliminate all detergent residue.

Brush the internal and external surfaces of the wrench with a clean brush with soft nylon bristles, under running water.

26

EN

8.5.2 Automatic Cleaning

Automatic cleaning can be carried out as an alternative to the manual cleaning described in Chapter 8.5.1 on page 43.

NOTE: For the handpiece, manual cleaning is recommended.

» REQUIRED MATERIALS

- Alkaline detergent^{b)};
- Neutralising liquid^{b)};
- Water;
- Metal basket;
- Adapters;
- Thermal-disinfector^{b)}.

b) Process validated with:

- Miele Professional Thermal-Disinfector, Model PG 8536;
- DES-VAR-TD program (Miele);
- Alkaline detergent: Neodisher® FA (0.2 % v/v);
- Neutralising liquid: Neodisher® Z (0.1 % v/v)

NOTE: Make sure that the accessories are suitably blocked in the basket and that they cannot move during the wash. Any blows could damage them. Position the instruments in a way that the water can flow through all the surfaces, even internal.

⚠ WARNING: Avoid overloading the thermal-disinfector as this could compromise the effectiveness of cleaning.

⚠ WARNING: On completion of the cleaning cycle in the thermal-disinfector, the scaler handpiece remains at the wash temperature for a long time. Use appropriate precautions when extracting the scaler handpiece from the thermal-disinfector to prevent injury to the operator.

⚠ CAUTION: Owing to its conformation, the scaler handpiece can rotate. When not in use, the scaler handpiece must always be placed on its support.

» CLEANING PROCEDURE

Position the accessories in a metal basket. Connect the relevant adapter (supplied as an optional) to the scaler handpiece connector and then to the water jet cleaning connections of the thermal-disinfector.


Repeat the same operation for the inserts by connecting them to the appropriate adapters supplied as optionals.

Sequence and parameters applicable to the cycle:

- 1 min, rinse with cold water;
- 5 min, wash with alkaline detergent^{b)} at 55°C ±2°C;
- 1 min, neutralisation with neutralising liquid^{b)}, 1/3 cold water and 2/3 hot water;
- 1 min, rinse with water (1/3 cold water, 2/3 hot water);
- 5 min, thermal disinfection at 93°C with demineralised water.

Thermal disinfection has not been tested experimentally. In compliance with ISO 15883-1, Table B.1 [4] thermal disinfection at a temperature of 90°C for 5 min determines a value of A0 3000.

1

2

8.6 Cleaning Control

» REQUIRED MATERIALS


- Light source;
- 2.5x Magnifier.

» PROCEDURE



Once the cleaning operations have been completed, check the scaler handpiece and scaler cone under an adequate source of light, if necessary using a magnifying glass 2.5X, paying attention to any parts that may conceal soil residue (threading, cavities, grooves) and, if need be, repeat the cleaning cycle if soil is still visible. Finally, check the integrity of those parts and those elements that could have deteriorated during use;

Repeat the control operations for the other accessories (inserts, tightening wrenches), repeating if necessary the cleaning cycle.

1

2

EN

8.7 Drying and Lubrication

» REQUIRED MATERIALS

- Compressed air;
- Soft, low-lint cloth;
- Medical grade lubricant.



» PROCEDURE

Thoroughly dry all parts of the scaler handpiece, scaler cone, and light pipe by blowing them with compressed air;


⚠ CAUTION: The scaler handpiece electrical contacts must be dry at the start and end of the sterilisation cycle, before the cord is connected to the device. Always make sure that the electrical contacts of the connector are perfectly dry, if necessary dry them blowing with compressed air.

⚠ CAUTION: Before starting the sterilisation cycle, make sure that the insert is thoroughly dry both internally and externally. To do this, blow compressed air both externally and through the internal passage hole. This will prevent the appearance of stains, streaks on the surface or oxidation inside the insert.

1


2


Dry the inserts tightening wrench using a soft, low-lint cloth;

3

Lubricate the inserts torque wrench with medical-grade lubricants in the point indicated;

4

CAUTION: Do not use oil or silicone-based lubricants. Use only medical grade lubricants.


EN

8.8 Sterilisation

» PREPARATION


Seal the scaler handpiece (without inserts), scaler cone, and light pipe individually and separately in disposable sterilisation bags.

1


Seal the inserts individually inside a disposable bag for sterilisation.

2

Seal the wrench individually inside a disposable bag for sterilisation.

3

» STERILISATION METHOD

The scaler handpiece and other sterilisable accessories are manufactured with materials that resist a maximum temperature of 135°C for a maximum time of 20 minutes.

Once the scaler handpiece and the other sterilisable accessories have been put into bags individually, perform the steam sterilisation process in the autoclave.

The sterilisation process validated by MECTRON S.p.A., in an autoclave with steam, guarantees a SAL 10^{-6} by setting the parameters indicated below:

- Type of cycle: 3 times Pre-vacuum (pressure min. 60 mBar).
- Minimum sterilisation temperature: 132°C (interval 0°C ÷ +3°C).
- Minimum sterilisation time: 4 minutes.
- Minimum drying time: 20 minutes.

All stages of sterilisation must be carried out by the operator in compliance with the most current version of standards: UNI EN ISO 17665-1, UNI EN ISO 556-1 and ANSI/AAMI ST:46.

⚠ CAUTION: Do not sterilise the scaler handpiece with the insert screwed in.

⚠ WARNING: Infection control - Sterilisable parts - Carefully remove any residual organic matter before sterilisation.

⚠ CAUTION: Carry out the sterilisation by using only an autoclave with steam from water. Do not use any other sterilisation procedure (dry heat, irradiation, ethylene oxide, gas, low temperature plasma, etc.).

⚠ CAUTION: Do not exceed the permissible load of the steam-steriliser.

⚠ WARNING: On completion of the sterilisation cycle in autoclave, the scaler handpiece remains at the sterilisation temperature for a long time. Use appropriate precautions when extracting the scaler handpiece from the autoclave to prevent injury to the operator.

⚠ CAUTION: Wait for the scaler handpiece to cool down completely before use.

8.8.1 Special Information

Sterilisation parameters, in steam autoclave, used in Great Britain:

- Temperature: 134° C;
- Time: 3 minutes.

9 MAINTENANCE

If the device is not used for more than one week, observe the following recommendations:

1. Run a complete cleaning cycle on the irrigation circuit using the **LONG FLUSH** function (see *Chapter 6 on page 29*);
2. Drain the circuits of any residual water by removing the tank and starting a **SHORT FLUSH** cycle (see *Chapter 5.2.3.1 on page 19*) until all the water in the irrigation circuit has been completely drained;
3. Disconnect the device from the mains;
4. Remove the peristaltic pump by following steps from 1 to 3 described in *Chapter 9.1 on page 57*;
5. Disassemble the device as described in *Chapter 7 on page 31*;
6. Store the device in its original packaging, in a suitable place (see "Conditions of transport and storage" in *Chapter 11 on page 59*);
7. Before to use the device again, reconnect all the accessories as described in *Chapter 4.3 on page 13*;
8. Perform a cycle of cleaning and sterilisation following the instructions in *Chapter 8 on page 33*;
9. Check that the inserts are not worn, deformed or broken, with particular attention to the integrity of the tip.

⚠ WARNING: Periodically check the integrity of the electric power supply cord; when it is damaged, replace it with MECTRON original spare part.


9.1 Replacing the Peristaltic Pump

On the back of the device there is the plastic cover which covers the housing of the peristaltic pump. Remove this protection by lifting it upwards.


⚠ CAUTION: Before carrying out operations on the peristaltic pump, make sure the device is unplugged and that the liquids container is not connected.

Remove the two tubes of the pump from the respective clutches positioned below it;

1


2


Slightly rotate the peristaltic pump clockwise and remove it from its housing;

3


Insert the new peristaltic pump and make sure that the pin of the impeller and the clip (upper and lower) coincide with the respective seats. Press up to the complete insertion of the clips. Insert the water conduits in the relative sleeves for the whole length of the latter;

4

Refit the plastic cover which covers the housing of the peristaltic pump by inserting it from the top and slide it onto the side guides until complete insertion.

5

10 METHODS AND PRECAUTIONS FOR DISPOSAL

WARNING: Hospital waste.

Treat the following objects as hospital waste:

- Inserts, when worn or broken;
- Inserts tightening wrench, when worn or broken.

Disposable materials and materials that carry biological risk must be disposed of according to local regulations in force regarding hospital waste.

 **WARNING:** When handling the inserts, pay particular attention to the sharp, pointed, irregular parts to avoid any wounds or injuries.

MICROPIEZO must be disposed of and treated as waste subject to separate collection.

NOTE: MICROPIEZO packaging is made from 100% recyclable materials that are easy to sort. Specifically: Paper (PAP) and Polyethylene (PE).

Failure to comply with the previous points can result in a penalty pursuant to the directive on waste electrical and electronic equipment (WEEE).

The buyer has the faculty to deliver the device to dispose of to the retailer who supplies the new equipment; instructions for proper disposal are available at MECTRON.

11 TECHNICAL DATA

Device compliant with Regulation (EU) 2017/745	Class IIa
Classification under the IEC/EN 60601-1	II Parts applied: type B (insert) IP 20 (device) IP 22 (pedal model FS-11 or FS-05)
Essential performance	According to the standard IEC 80601-2-60 the device has no essential performance
Device for intermittent operation	55 sec. ON - 30 sec. OFF with irrigation 30 sec. ON - 120 sec. OFF without irrigation
Power Supply	100-240 V~ 50/60 Hz
Max. current	1,5A (---24V)
Battery	The device is equipped with a backup battery, which can only be replaced by an Authorised MECTRON Service Centre. Model: BR/CR 1220
Working Frequency	Automatic Scan From 24 KHz to 36 KHz
Powers type	endo perio-scaling prosthetics soft-mode

Irrigation Circuit	Peristaltic pump flow rate: Adjustable using the touch screen 7 levels of flow rate: from 0 (0 ml/min) to 6 (33 ml/min $\pm 10\%$) (see <i>Chapter 5.2.2 on page 18</i>) Tank capacity: 500 ml. Tank lighting system: LED light power risk free according to standard IEC/EN 62471.
LED system of the Handpiece	Function Light "on": the LED light is turned on when pressing the pedal and stays on for the next 3 seconds after releasing the pedal. White LED light power risk free according to standard IEC/EN 62471.
Protections of the APC Circuit	Missing handpiece; Breaking wire cord; Insert not correctly tightened or broken; Problems to the engine of the peristaltic pump; Keyboard failure; Unforeseen value of the supply voltage.
Wi-Fi / BLE operating frequency range	2400 ÷ 2483.5 MHz , maximum power limit 20 dBm
Operating Conditions	from 10°C to 35°C Relative humidity from 30% to 75% Air pressure P: 700hPa/1060hPa
Transport and Storage Conditions	from -10 °C to 45 °C Relative humidity from 10% to 90% Air pressure P: 500hPa/1060hPa
Altitude	lower than or equal to 3000 metres
Weights and dimensions	Maximum weight: 940 g 161 x 214 x 171 mm (W x L x H)^{a)}

a) W = width ; L = length ; H = height

11.1 Software

Following is the list of third-party software installed and used in the MICROPIEZO device:

- Espressif IoT Development Framework (esp-idf), Copyright© 2015-2019 Espressif Systems;
- FreeRTOS, MIT open source license.

11.2 Electromagnetic Compatibility IEC/EN 60601-1-2

⚠ WARNING: Contraindications. Interference with other equipment
Even if compliant with standard IEC 60601-1-2, MICROPIEZO may interfere with other devices in the vicinity.
MICROPIEZO must not be used near or stacked with other devices. However, if this should become necessary, it is necessary to check and monitor the correct functioning of the device in that configuration.

⚠ WARNING: Portable and mobile radio communication equipment may influence the correct operation of the device.

⚠ WARNING: Contraindications. Interference from other equipment
An electrosurgical scalpel or other electrosurgical units positioned near the MICROPIEZO device may interfere with the correct operation of the device itself.

⚠ WARNING: The device requires particular EMC precautions and must be installed and put into service according to the EMC information provided in this chapter.

⚠ WARNING: The use of cables and accessories not supplied by MECTRON may adversely affect the EMC performances.

EN

11.3 Guide and Manufacturer’s Declaration - Electromagnetic Emissions

MICROPIEZO is designed to operate in the electromagnetic environment specified below. The client or user of MICROPIEZO should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment Guidance
RF Emissions CISPR 11	Group 1	MICROPIEZO uses RF energy only for its internal operation. Its RF emissions are therefore very low and unlikely to cause any interference with nearby electronic devices.
RF Emissions CISPR 11	Class B	MICROPIEZO is suitable for use in a professional healthcare environment including hospitals, medical clinics, surgical centres, and specific healthcare facilities where the equipment and systems are managed by qualified and specialised healthcare personnel.
Harmonic emissions IEC 61000-3-2	Class A	
Emissions of fluctuations voltage/flicker IEC 61000-3-3	Compliant	

11.4 Accessible Parts of the Casing

MICROPIEZO is designed to operate in the electromagnetic environment specified below. The client or user of MICROPIEZO should ensure that it is used in such an environment.

Phenomenon	Basic EMC standard or test method	Immunity test levels	Electromagnetic Environment Guidance
Electrostatic discharge (ESD)	IEC 61000-4-2	±8 kV on contact ±2 kV, ±4 kV, ±8 kV, ±15 kV in air	The floor must be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic material, the relative humidity should be at least equal to 30%.
Radiated RF EM fields ^{a)}	IEC 61000-4-3	3 V/m ^{f)} 80 MHz - 2,7 GHz ^{b)} 80% AM a 1 kHz ^{c)}	Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended distances, calculated from the equation applicable at the frequency of the transmitter.
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See Chapter 11.6 on page 67	
RATED power frequency magnetic fields ^{d)}	IEC 61000-4-8	30 A/m 50 Hz o 60 Hz	The magnetic fields at the mains frequency should have levels characteristic of a typical location in a commercial or hospital environment.
Proximity magnetic fields	IEC 61000-4-39	See Chapter 11.7 on page 68	Portable and mobile RF communication devices shall be used with a separation distance of at least 0,15 m from the field sources.

- a) The interface between the PATIENT physiological signal simulation, if used, and the MICROPIEZO shall be located within 0,1 m of the vertical plane of the uniform field area in one orientation of the MICROPIEZO.

b) MICROPIEZO that intentionally receives RF electromagnetic energy for the purpose of its operation shall be tested at the frequency of reception. Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS. This test assesses the BASIC SAFETY and ESSENTIAL PERFORMANCE of an intentional receiver when anambient signal is in the passband. It is understood that the receiver might not achieve normal reception during the test.
- c) Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.

d) Applies only to ME EQUIPMENT and ME SYSTEMS with magnetically sensitive components or circuitry.

e) Void.

f) Before modulation is applied.

11.5 Guide and the Manufacturer's Declaration - Electromagnetic Immunity

11.5.1 Power Connection BC Input

MICROPIEZO is designed to operate in the electromagnetic environment specified below. The client or user of MICROPIEZO should ensure that it is used in such an environment.

Phenomenon	Essential EMC standard or test method	Immunity test values	Electromagnetic Environment Guidance
Electrical fast transient/burst l) o)	IEC 61000-4-4	± 2 kV on contact 100 KHz repetition frequency	The quality of the network voltage should be that of a typical commercial or hospital environment.
Pulses Differential mode b) j) o)	IEC 61000-4-5	± 0.5 kV, ± 1 kV	The quality of the network voltage should be that of a typical commercial or hospital environment.
Pulses common mode b) j) k) o)	IEC 61000-4-5	± 0.5 kV, ± 1 kV, ± 2 kV	The quality of the network voltage should be that of a typical commercial or hospital environment.
Conductive disturbances induced by RF fields c) d) o)	IEC 61000-4-6	3 V ^{m)} 0.15 MHz - 80 MHz 6 V ^{m)} in the ISM bands between 0.15 MHz and 80 MHz ⁿ⁾ 80 % AM at 1 KHz ^{e)}	Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended and calculated distances from the equation applicable at the frequency of the transmitter.
Voltage dips f) p) r)	IEC 61000-4-11	0% UT; 0,5 cycle ^{g)} At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° ^{q)}	The quality of the network voltage should be that of a typical commercial or hospital environment.
		0 % UT; 1 cycle ^e 70 % UT; 25/30 cycle ^{h)} Single phase: at 0°	
Voltage interruptions f) i) o)	IEC 61000-4-11	0% UT; 250/300 cycle ^{h)}	The quality of the network voltage should be that of a typical commercial or hospital environment.

- a) Void.
- b) All MICROPIEZO cables are attached during the test.
- c) Calibration for current injection clamps shall be performed in a 150 Ω system.
- d) If the frequency stepping skips over an ISM or amateur band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.
- e) Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- f) ME EQUIPMENT and ME SYSTEMS with a d.c. power input intended for use with a.c.-to-d.c. converters shall be tested using a converter that meets the specifications of the MANUFACTURER of the ME EQUIPMENT or ME SYSTEMS. The IMMUNITY TEST LEVELS are applied to the a.c. power input of the converter.
- g) Applicable only to MICROPIEZO connected to single-phase a.c. mains.
- h) E.g. 10/12 means 10 periods at 50 Hz or 12 periods at 60 Hz.
- i) ME EQUIPMENT and ME SYSTEMS with RATED input current greater than 16 A / phase shall be interrupted once for 250/300 cycles at any angle and at all phases at the same time (if applicable). MICROPIEZO with battery backup shall resume line power operation after the test. For ME EQUIPMENT and ME SYSTEMS with RATED input current not exceeding 16 A, all phases shall be interrupted simultaneously.
- j) ME EQUIPMENT and ME SYSTEMS that does not have a surge protection device in the primary power circuit may be tested only at ± 2 kV line(s) to earth and ± 1 kV line(s) to line(s).
- k) Not applicable to CLASS II MICROPIEZO.
- l) Direct coupling shall be used.
- m) r.m.s. , before modulation is applied.
- n) The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.
- o) Applicable to ME EQUIPMENT and ME SYSTEMS with RATED input current less than or equal to 16 A / phase and ME EQUIPMENT and ME SYSTEMS with RATED input current greater than 16 A / phase.
- p) Applicable to ME EQUIPMENT and ME SYSTEMS with RATED input current less than or equal to 16 A / phase.
- q) At some phase angles, applying this test to ME EQUIPMENT with transformer mains power input might cause an overcurrent protection device to open. This can occur due to magnetic flux saturation of the transformer core after the voltage dip. If this occurs, the ME EQUIPMENT shall provide BASIC SAFETY during and after the test.
- r) For ME EQUIPMENT and ME SYSTEMS that have multiple voltage settings or auto ranging voltage capability, the test shall be performed at the power input voltage specified in Table 1 - "Power input voltages and frequencies during the tests" of the IEC 60601-1-2:2014/AMD1:2020.

11.5.2 Points of Contact with the Patient

MICROPIEZO is designed to operate in the electromagnetic environment specified below. The client or user of MICROPIEZO should ensure that it is used in such an environment.

Phenomenon	Essential EMC standard or test method	Immunity test values	Electromagnetic Environment Guidance
Electrostatic discharges (ESD) ^{c)}	IEC 61000-4-2	±8 kV on contact ±2 kV, ±4 kV, ±8 kV, ±15 kV in air	The floor must be made of wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Conductive disturbances induced by RF fields ^{a)}	IEC 61000-4-6	3 V ^{b)} 0.15 MHz - 80 MHz 6 V ^{b)} in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 KHz	Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended distances, calculated from the equation applicable at the frequency of the transmitter.

a) The following applies:

- All connection cables with the patient must be tested, either individually or grouped together.
- The connection cables with the patient must be tested using a current clamp unless the current clamp is not suitable. If a current clamp is not suitable, an EM clamp must be used.
- In any case, no intentional decoupling device should be used between the injection site and POINT OF CONNECTION TO THE PATIENT.
- The tests can be performed with other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- The tubes that are intentionally filled with conductive liquids and intended to be placed in contact with the PATIENT must be considered connection cables with the patient.
- If an ISM or amateur radio band is not present among the frequency samples, as appropriate, an additional test frequency has to be used in the ISM

band or in the amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.

- The ISM bands (industrial, scientific and medical) between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz. The amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz and 50.0 MHz to 54.0 MHz.

b) R.M.S., applied before modulation.

c) The discharges must be applied without connection to an artificial hand and without connection to the PATIENT simulation. The PATIENT simulation can be connected after the test, if necessary, to verify BASIC SAFETY and ESSENTIAL PERFORMANCE.

11.5.3 Parts Accessible to the Input / Output Signals

MICROPIEZO is designed to operate in the electromagnetic environment specified below. The client or user of MICROPIEZO should ensure that it is used in such an environment.

Phenomenon	Essential EMC standard or test method	Immunity test values	Electromagnetic Environment Guidance
Electrostatic discharges (ESD) ^{e)}	IEC 61000-4-2	±8 kV on contact ±2 kV, ±4 kV, ±8 kV, ±15 kV in air	The floor must be made of wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst ^{b) f)}	IEC 61000-4-4	±1 kV on contact 100 KHz repetition frequency	The quality of the network voltage should be that of a typical commercial or hospital environment.
Pulses Common mode ^{a)}	IEC 61000-4-5	± 2kV	The quality of the network voltage should be that of a typical commercial or hospital environment.
Conductive disturbances induced by RF fields ^{d) g) j) k)}	IEC 61000-4-6	3 V ^{h)} 0.15 MHz - 80 MHz 6 V ^{h)} in the ISM bands between 0.15 MHz and 80 MHz ⁱ⁾ 80 % AM at 1 KHz ^{c)}	Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended distances, calculated from the equation applicable at the frequency of the transmitter.

- a) This test applies only to output lines intended to connect directly to outdoor cables.
- b) SIP/SOPS whose maximum cable length is less than 3 m in length are excluded.
- c) Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- d) Calibration for current injection clamps shall be performed in a 150 Ω system.
- e) Connectors shall be tested per 8.3.2 and Table 4 of IEC 61000-4-2:2008. For insulated connector shells, perform air discharge testing to the connector shell and the pins using the rounded tip finger of the ESD generator, with the exception that the only connector pins that are tested are those that can be contacted or touched, under conditions of INTENDED USE, by the standard test finger shown in Figure 6 of the general standard, applied in a bent or straight position.
- f) Capacitive coupling shall be used.

- g) If the frequency stepping skips over an ISM or amateur radio band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.
- h) R.M.S., before modulation as applied.
- i) The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0, 15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18, 17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29, 7 MHz and 50,0 MHz to 54,0 MHz.
- j) See IEC 61000-4-6:2013, Annex B, for modified start frequency versus cable length and equipment size.
- k) SIP/SOPS whose maximum cable length is less than 1 m are excluded.

11.6 Specifications of the Tests for the Immunity of the Accessible Parts of the Casing to the Wireless RF Communications Device

MICROPIEZO is designed to operate in an electromagnetic environment in which radiated RF disturbances are under control. The client or operator of MICROPIEZO can help prevent electromagnetic interferences by guaranteeing a minimum distance between the mobile and portable RF communication devices (transmitters) and MICROPIEZO, as recommended below, in relation to the maximum output power of the radio communication devices.

Test Freq. (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Max power (W)	Distance (m)	Immunity test level (V/m)
385	380 - 390	TETRA 400	Pulse modulation ^{b)} 18 Hz	1.8	0.3	27
450	430 - 470	GMRS 460 FRS 460	FM ^{c)} ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704 - 787	LTE band 13, 17	Pulse modulation ^{b)} 217 Hz	0.2	0.3	9
745						
780						
810	800 - 960	GSM 800/900 TETRA 800 iDEN 820 CDMA 850 Band LTE 5	Pulse modulation ^{b)} 18 Hz	2	0.3	28
870						
930						
1720	1700 - 1990	GSM 1800 CDMA 1900 GSM 1900 DECT LTE Band 1, 3, 4, 25 UMTS	Pulse modulation ^{b)} 217 Hz	2	0.3	28
1845						
1970						
2450	2400 - 2570	Bluetooth WLAN 802.11 b/g/n RFID 2450 Band LTE 7	Pulse modulation ^{b)} 217 Hz	2	0.3	28
5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0.2	0.3	9
5500						
5785						

a) For some services, only the uplink frequencies are included.

b) The carrier shall be modulated using a 50% duty cycle square wave signal.

c) As an alternative to FM modulation, the carrier may be pulse modulated using a 50% duty cycle square wave signal at 18 Hz. While it does not represent actual modulation, it would be worst case.

NOTE: If necessary to reach the immunity test level, the distance between the transmitter antenna and MICROPIEZO can be reduced to 1 m. The test distance of 1 m is allowed by IEC 61000-4-3.

⚠ WARNING: Portable RF communication equipment (including peripheral devices such as antenna cables and external antennas) must not be used closer than 30 cm to any part of the MICROPIEZO device, including the cables specified by the manufacturer. Otherwise, there may be a performance degradation of these devices.

11.7 Immunity to Proximity Magnetic Fields in the Frequency Range 9 kHz to 13,56 MHz

The following table reports the test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields in the frequency range 9 kHz to 13,56 MHz.

Test Frequency	Modulation	Immunity test level (A/m)
30 kHz ^{a)}	CW	8
134,2 kHz	Pulse modulation ^{b)} 2,1 kHz	65 ^{c)}
13,56 MHz	Pulse modulation ^{b)} 50 kHz	7,5 ^{c)}

- a) This test is applicable only to devices intended for use in the HOME HEALTHCARE ENVIRONMENT.

b) The carrier shall be modulated using a 50% duty cycle square wave signal.

c) r.m.s., before modulation is applied.

12 TROUBLESHOOTING

12.1 Diagnostic System and Symbols on the Keyboard

MICROPIEZO is provided with a diagnostic circuit that allows the detection of malfunctions and viewing of their type on the touch keyboard by means of a symbol. Users, by using the following table, are guided to identifying and the possible resolution of the malfunction detected.

Symptom	Possible Cause	Solution
 Yellow "handpiece" symbol active and two short acoustic signals when the pedal is pressed.	Handpiece/cord electrical contacts wet.	Dry contacts with compressed air. ^{b)}
	MICROPIEZO handpiece not connected to the device.	Connect the handpiece. ^{b)}
	Handpiece not detected. Faulty Handpiece.	Replace the handpiece. ^{b)}
 Yellow "insert" symbol active and two short acoustic signals when the pedal is pressed.	Insert not present or not properly tightened on the handpiece.	Unscrew the insert and screw it back correctly using the torque wrench (see <i>Chapter 5.4 on page 25</i>). ^{b)}
	Insert broken, worn or deformed.	Replace the insert. ^{b)}
	Electrical contacts wet cord/handpiece.	Dry contacts with compressed air. ^{b)}
	Faulty Handpiece.	If the previous solutions do not solve the problem, replace the handpiece. ^{b)}
	Tuning circuit malfunction.	Contact an Authorised MECTRON Service Centre.

EN

Symptom	Possible Cause	Solution
 Yellow "wrench" symbol active and prolonged acoustic signal.	Handpiece back in place during treatment. The handpiece has been placed in the appropriate housing available on the system with the pedal pressed and a power and/or irrigation level higher than 0.	Release the pedal. Remove the handpiece from housing available on the system and press the pedal again to continue with the treatment.
	FLUSH cycle started with handpiece placed in the appropriate housing available on the system.	With handpiece put back, the FLUSH cycle will not start. If you DO NOT want to perform a FLUSH cycle: • Press twice the button FLUSH to exit short FLUSH mode. • Press once the button FLUSH to exit long FLUSH mode. If you want to perform a flush cycle: • Remove the handpiece from its seat. • Press the pedal and release it to start the procedure.
	Malfunction peristaltic pump.	Check that there are no impediments for the rotation of the pump. Check that the pump and the two pipes are correctly installed.
	Faults on the electrical network or excessive electrostatic discharge or internal faults.	Turn off and wait 5 seconds before switching on the device again. If the warning persists, contact an Authorised MECTRON Service Centre.
	No communication with the keyboard.	Contact an Authorised MECTRON Service Centre.
	Unforeseen value of the supply voltage.	Replace the handpiece. If the problem persists, contact an Authorised MECTRON Service Centre.

Symptom	Possible Cause	Solution
 Yellow “wrench” symbol active and no acoustic signal.	Corrupt or invalid firmware.	In case of corrupt or invalid firmware, the system will automatically attempt a reset. Wait a few minutes and check that the warning disappears. If the warning persists, contact an Authorised Mectron Service Centre.
 The button FLUSH and the indication of the “wrench” flash. ^{c)} The device emits two prolonged and consecutive beeps. ^{a)}	No LONG FLUSH cycle has been performed during the previous 8 hours of operation (device turned on) or 24 hours.	Perform a “long” FLUSH cycle.
 The value 0 of the irrigation scale and the indication of the wrench flash. ^{c)} The device emits two prolonged and consecutive beeps. ^{a)}	The peristaltic pump has worked for more than 400 hours or for more than 1 calendar year.	1. Turn the machine off. 2. Replace the peristaltic pump (see <i>Chapter 9.1 on page 57</i>). 3. Turn on the machine and press 2 consecutive times the value 0 of the irrigation scale to reset the error signal. The device emits a confirmation tone.

Symptom	Possible Cause	Solution
 The value 0 of the power scale and the indication of the wrench flash. The device emits two prolonged and consecutive beeps.^{a)}	The system was started with the pedal pressed.	Check that during system start-up the pedal is not pressed. If the pedal is not pressed, check that the pedal returns correctly to the rest position and there is nothing that prevents the movement. If the warning persists, contact an Authorised MECTRON Service Centre
 Wi-Fi button flashing.	The Wi-Fi interface is active, but the device was unable to establish a connection with the MECTRON server	Check that it is possible to access the internet via the Wi-Fi network to which MICROPIEZO is connected. Perform the procedure for connecting to the wi-fi network.

- a) For 10 seconds when the machine is switched on.
- b) The error condition is manifested when the pedal is pressed and with a “Power Level” different to 0. If there is no error condition is present any more, the warning signal will be updated at the next pedal pressure.
- c) When the pedal is pressed, the lighted diagnostic indications will be deactivated. When the pedal is released, if the error conditions persist, the relative indications will be reactivated.

12.2 Quick Troubleshooting

If the device does not seem to work properly, read the instructions again and then check the following table.

Symptom	Possible Cause	Solution
The device does not switch on after pressing the Power button for 2 consecutive seconds.	The terminal of the power supply cable is incorrectly inserted in the rear socket of the device.	Check that the power supply cable is securely connected.
	The electric power cable is defective.	Check that the supply socket is working. Replace the power supply cable.
The device is on but is not working. The display reports no errors.	The pin of the pedal is not properly inserted in the socket of the device.	Properly insert the plug of the pedal in the socket on the back of the device.
	The pedal does not work.	Contact an Authorised MECTRON Service Centre.
	The power level on the keyboard is set to "0".	Adjust the power level to a value greater than 0.
The device is on but is not working. On the keyboard there are symbols lit.	See <i>Chapter 12.1 on page 69</i> for the possible cause depending on the symbol.	See <i>Chapter 12.1 on page 69</i> for the action to be taken depending on the symbol.
During operation a faint whistling noise can be heard coming from the MICROPIEZO handpiece.	The insert is not correctly tightened on the handpiece.	Unscrew and screw the insert correctly using the torque wrench MECTRON (see <i>Chapter 5.4 on page 25</i>).
	The irrigation circuit has not been completely filled.	Fill the irrigation circuit using the FLUSH function (See <i>Chapter 5.4 on page 25</i>).

Symptom	Possible Cause	Solution
No liquid flows out from the insert during operation.	The insert is of the type which does not provide the passage of liquids.	Use an insert of the type with the passage of liquids.
	The the insert is blocked.	Unscrew the insert from the handpiece and release the water passage of the insert by blowing compressed air through it. If the problem persists, replace the insert with a new one
	The handpiece is blocked.	Contact an Authorised MECTRON Service Centre.
	The irrigation level on the keyboard is set to "0".	Adjust the irrigation level to a value greater than 0.
	The liquid tank is empty.	Fill the tank.
	The tank is not properly installed.	Properly connect the tank to the machine body.
	The silicone tubes of the pump are not correctly installed.	Check the connections of the tubes.
	The peristaltic pump is worn.	Replace the peristaltic pump (See <i>Chapter 9.1 on page 57</i>).
Poor performance.	The insert is not correctly tightened on the handpiece.	Unscrew and screw the insert correctly using the torque wrench MECTRON (see <i>Chapter 5.4 on page 25</i>)
	Insert broken, worn or deformed.	Replace the insert with a new one.

12.3 Shipping to an Authorised MECTRON Service Centre

If technical assistance is required on the machine, contact an Authorised MECTRON Service Centre or your Dealer. Do not try to repair or modify the device and its accessories. Clean and sterilise all parts that need to be sent to an Authorised MECTRON Service Centre by following the instructions in *Chapter 8 on page 33*.

Leave the sterilised parts in the bag which certifies the sterilisation process.

The cleaning and sterilisation demands comply with the mandatory requirements for workplace health and safety protection pursuant to Legislative Decree 81/08 and Italian laws.

If the client fails to comply with the requirements, MECTRON reserves the right to charge them the cost of cleaning and sterilisation or to refuse goods that arrive in unsuitable conditions, returning them at the client's own expense for proper cleaning and sterilisation.

The device must be returned suitably packaged and accompanied by all the accessories and a form including:

- Owner's data with telephone number
- Product Name
- Serial number and/or batch number
- Reason for return / description of failure
- Photocopy of the bill or invoice for the purchase of the device

CAUTION: Packaging

Pack the device in its original packaging to prevent damage during transport. If the device is returned in non-original or otherwise inadequate packaging with consequent damage caused by transport, MECTRON reserves the right to refuse assistance and to return the device at the expense of the end client.

Once the material is received at the Authorised MECTRON Service Centre, the qualified technical personnel will evaluate the problem. The repair will be made only upon acceptance by the end client. For further details please contact your nearest Authorised MECTRON Service Centre or dealer.

Unauthorised repairs can damage the system and void the warranty, releasing MECTRON of all liability for direct or indirect damages to persons or property.

13 WARRANTY

Before being marketed, all MECTRON devices, are subjected to a thorough final check that verifies the full functionality.

MECTRON warrants its products, purchased new from a MECTRON dealer or importer, against defects in material and workmanship for a period of 2 (two) years for the device and 1 (ONE) YEAR for the handpiece, from the date of purchase.

During the warranty period, MECTRON undertakes to repair (or, at its discretion, to replace) the parts of products free of charge, which, according to its judgement, are proven to be defective. The assessments of MECTRON for the validity or otherwise of the warranty, can also take place on the basis of data collected through the use of IoT technology. Complete replacement of Mectron products is excluded.

MECTRON disclaims all liability for direct or indirect damage to persons or property, in the following cases:

- the device is not used according to the intended use for which it is provided;
- the device is not used in accordance with all the instructions and requirements described in this manual;
- the electrical system of the places where the equipment is used do not comply with the laws in force and the related regulations;
- Assembly operations, extensions, re-adjustments, modifications, replacements and repairs are carried out by personnel not authorised by Mectron or in breach of what is provided in this manual also in regard to the origin of the authorised material;
- The environmental conditions for preservation and storage of the device are not complying with the requirements indicated in *Chapter 11 on page 59*;

- The installation or transport of the device is not performed as specified in this manual or in other documentation provided by MECTRON, or in any case available on the website of the latter;
- The device or component thereof is purchased from a subject not authorised by MECTRON;
- The device, including its subcomponents, parts or assemblies, are altered or modified with respect to the indications in this manual;
- Accidents, misuse, abuse, abnormal use, negligent use, intentional misconduct or use exceeding the allowed and recommended limits of the device, or in case of its normal wear or deterioration;
- If the defect or non-conformity has not been promptly and readily communicated in writing to MECTRON as set out in this manual;
- If the damage, costs, or expenses are caused by force majeure;
- The connection of the device has been carried out at a voltage different from the one envisaged, including warning lights, knobs and all the accessories.

In any case, MECTRON will not recognise indemnity or compensation for loss of use, inconvenience, loss of profits, loss of business, loss of business opportunities, damage to reputation, and any incidental or consequential damages arising out of or relating to the device.

The expected service life of the device is minimum 5 years.

The service life / duration does not establish a limit of use; the service life of the device defines the period of time after installation and/or commissioning, during which the original performance, or in any case performance suited to the intended use, is guaranteed without there being any degradation such as to compromise functionality and reliability.

The service life is a minimum qualitative

objective of the design, therefore, individual parts or components may guarantee superior performance and reliability with respect to that declared by the manufacturer.

The service life assumes compliance with the maintenance schedules set out in this manual, does not include components normally subject to "wear", and is not linked to the warranty period: the service life does not establish any implicit or explicit extension of the warranty period.

CAUTION

The warranty starts from the date of purchase of the device, which evidence is given by the delivery note/purchase invoice issued by the Dealer / Importer or, in case of device with activation code, from the date of activation of the same.

In order to avail of the warranty service, the client must return the device, at their own expense, to the MECTRON dealer/importer from which they purchased the product.

The device must be returned together with the original packaging, accompanied by all

the accessories and by a form containing:

- The data of the owner and telephone number;
- The data of the Dealer / Importer;
- Photocopy of the delivery note/purchase invoice of the device by the owner where are reported the date, the name of the device and the serial number;
- Description of the failure.

The transport and the damage caused by transport are not covered by the warranty.

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© Mectron S.p.A. 2024. Tutti i diritti riservati. Nessuna parte di questo documento può essere riprodotta, in qualsiasi forma, senza il consenso scritto del proprietario del copyright.