mectron

medical technology

USE AND MAINTENANCE MANUAL

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EN

ES

turbodent touch

CE
0051

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turbodent touch

Copyright

© Mectron S.p.A. 2025. All rights reserved. No part of this document may be reproduced, in any form whatsoever, without the prior written consent of the copyright holder.
The images are for demonstration purposes only.

EN

CONTENTS

1	Introduction	1
1.1	Intended Use	2
1.2	Description of the Device	2
1.2.1	Eligible Patient Group	3
1.2.2	Patient Selection Criteria	3
1.2.3	Indications for Use	3
1.2.4	Users	3
1.2.5	Operating environment	3
1.3	Disclaimer	4
1.4	Safety Requirements	4
1.5	Symbols	6
2	Identification Data	8
2.1	Device Nameplate	8
2.2	Accessories Identification Data	8
3	Delivery	9
3.1	List of components	9
4	Installation	12
4.1	First Installation	12
4.2	Safety Requirements During Installation	12
4.3	Connecting the components	13
5	Use	16
5.1	Power On/Off	16
5.2	Description of the Touch Surface	17
5.3	Safety Requirements Before and During Use	19
5.4	Using the Handpiece	21
5.5	“refill” function	24
5.6	“flush” function	25
5.7	Stand-by Mode (Energy Saving)	27
6	Cleaning	27
6.1	Disassembling Parts for Cleaning and Sterilisation	28
6.2	Preparation	30
6.3	Cleaning the Non-Sterilisable Parts	31
6.4	Cleaning Sterilisable components	32
6.4.1	Pre-Cleaning	32
6.4.2	Manual Cleaning	33
6.4.3	Automatic Cleaning	39
7	Cleaning Check	40
8	Drying	41
9	Sterilisation	42
10	Disinfecting the Irrigation Circuit	44
11	Decontamination of the Irrigation Circuit	50
12	Maintenance	50
12.1	Maintenance After Every Treatment	50
12.2	Daily Maintenance	50
12.3	Transport or Long Periods of Non-use	51
12.4	Unclogging the Air-polishing Circuit	52
12.5	Unclogging the AIR-POLISHING handpiece	53
12.6	Cleaning and/or Replacing the Water Filter and the Water Inlet O-ring	55
12.7	Cleaning the Powder Containers and Caps	57
12.8	Replacing the Air-polishing Cord O-Rings	58

12.9	Eliminating the Condensate	59
12.10	Maintenance table	59
13	Methods and Precautions for Disposal	60
14	Technical Data	61
14.1	Electromagnetic Compatibility IEC 60601-1-2	62
14.2	Guide and Manufacturer's Declaration - Electromagnetic Emissions	62
14.3	Accessible Parts of the Casing	63
14.4	Guide and the Manufacturer's Declaration - Electromagnetic Immunity	64
14.4.1	Input A.C. Power Connection	64
14.4.2	Points of Contact with the Patient	66
14.4.3	Parts Accessible to the Input / Output Signals	67
14.5	Specifications of the Tests for the Immunity of the Accessible Parts of the Casing to the Wireless RF Communications Device	68
14.6	Proximity Magnetic Field Immunity in the Frequency Range from 9 kHz to 13,56 MHz	69
15	Troubleshooting	70
15.1	Diagnostic Reports and Solutions	70
15.2	Quick Troubleshooting	71
15.3	Replacing the Fuses	74
15.4	Shipping to an Authorised MECTRON Service Centre	75
16	Warranty	76

PAGE LEFT INTENTIONALLY BLANK

1 INTRODUCTION

This manual applies to the following medical devices:

- turbodent touch (referred to as "device" in the text)
- AIR-POLISHING PERIO handpiece (referred to as "accessory/ies" in the text)
- AIR-POLISHING 90° handpiece (referred to as "accessory/ies" in the text)
- AIR-POLISHING 120° handpiece (referred to as "accessory/ies" in the text)

Read this manual carefully before installing, operating, performing maintenance or otherwise working on the device and its accessories.

This manual must always be available to the operator.

Important: To avoid damages to persons or property, carefully read all the "Safety Requirements" in the manual.

Depending on the level of severity, the safety information is classified with the following indications:

 **WARNING:** always refers to personal injury.

 **CAUTION:** refers to possible damage to property.

The aim of this manual is to make operators aware of the safety regulations, installation procedures, instructions for proper use and maintenance of the device and its components. Do not use this manual for purposes other than those strictly related to the installation, operation and maintenance of the device.

The information and illustrations in this manual were updated on the edition date shown on the last page.

MECTRON is committed to continuously updating its products, which may result in modifications to the device and its components.

If discrepancies are found between the content of this manual and the devices in your possession, you may:

- check for any available updates in the *MANUALS* section of the MECTRON website¹;
- ask clarifications to Your Dealer;
- contact the MECTRON After-Sales Service.

¹ <https://manuals.mectron.com/>

1.1 Intended Use

It is a device equipped with a dedicated handpiece that utilizes air-polishing technology for dental prophylaxis. Its purpose is to remove supragingival stains and both supra- and subgingival plaque from natural and prosthetic structures as well as orthodontic appliances in the oral cavity.

AIR-POLISHING 90° / AIR-POLISHING 120°

It is a device used in conjunction with air-polishing technology for supra- and subgingival dental prophylaxis.

AIR-POLISHING PERIO

It is a device used in conjunction with air-polishing technology for subgingival dental prophylaxis to be used in combination with SUBGINGIVAL PERIO TIPS.

 **WARNING:** Only use the device and its accessories for their intended purpose. Failure to observe this regulation can cause serious injury to the patient, the operator and damage/faults on the device.

1.2 Description of the Device

turbodent touch is a table-top device, equipped with a polisher handpiece, which uses air-polishing technology to offer complete treatment in dental prophylaxis. Its aim is to remove stains at supragingival level on natural teeth and restorations, as well as eliminate plaque at supra- and subgingival level on teeth, braces, restorations and implants.

The operating principle of the polisher is based on the mechanical action obtained from a jet of various crystal types accelerated by a flow of compressed air. The kinetic energy thus imprinted at the particles, dissipates almost completely due to impact against the surface of the enamel, producing a gentle but effective cleaning action. The action is completed by a jet of water which, using the vacuum created around the nozzle, has a bell-shape around the main flow, thus producing a double effect: to prevent much of the rebound and the leakage of the cloud of powder and perform continuous washing of the treated area, dissolving the powder.

 **WARNING:** The device and its accessories must be used in a dental practice or clinic, or in professional oral hygiene and preventive care centres. Do not use the device and its accessories in environments where the atmosphere is saturated with flammable gases (anaesthetic mixtures, oxygen, etc.).

 **WARNING: Qualified and specialist personnel.** The device and its accessories must be used exclusively by specialised personnel such as doctors/dentists or dental hygienists. No specific training is required for use of the device and of its accessories. The device and its accessories produce no side effects when used correctly.  **WARNING:** For composite resin surfaces, the jet should be directed exclusively at the part to be treated, for an average time of 2-3 seconds per tooth.

NOTE: For detailed information on SUBGINGIVAL PERIO TIPS, always refer to the relevant instructions for use.

1.2.1 Eligible Patient Group

This medical device is designed to be used with the following patient groups:

- Children;
- Adolescents;
- Adults;
- Seniors.

This medical device can be used on patients of any weight, height, sex or nationality, and who are over 2 years of age.

1.2.2 Patient Selection Criteria

It is not recommended to use the device and its relative handpiece in the following cases:

- | | |
|---|---|
| 1. Upper respiratory tract infections, chronic bronchitis/asthma; | 3. Patients undergoing treatment (radiotherapy, chemotherapy, antibiotics); |
| 2. Pregnant and breastfeeding women; | 4. Acute infectious oral lesions. |

It is not recommended to use the powders in the following cases:

1. Allergy to the aroma of the powder;
2. Patients on a low-sodium diet or suffering from severe respiratory problems, such as chronic bronchitis, asthma, emphysema, etc. unless specifically indicated by the surgeon (for bicarbonate powders only).

All MECTRON devices are intended for professional use only. Therefore, it is up to the operator to decide, based on his or her own experience, how and whether to treat patients.

NOTE: It is recommended to consult the Instructions for Use of the powders before operation to check for any contraindications or specific indications.

1.2.3 Indications for Use

The use of turbodent touch, equipped with a polisher handpiece, is indicated for all intended patients (see *Chapter 1.2.1 on page 3*) who require treatment within the scope of the intended use of the device (see *Chapter 1.1 on page 2*).

1.2.4 Users

turbodent touch and the polisher handpiece must be used exclusively by specialised and properly trained personnel, specifically the dental surgeon/dentist or dental hygienist, who must be able-bodied, adult, of any weight, age, height, gender and nationality.

1.2.5 Operating environment

The device is portable. It is intended for use in outpatient, private or hospital settings; in the absence of flammable mixtures, liquids, powders; far from other devices and/or electro-medical equipment.

1.3 Disclaimer

The manufacturer MECTRON disclaims all liability, express or implied, and cannot be held liable for direct or indirect personal injury and/or property damage, occurring as a result of incorrect procedures related to use of the device and its components/accessories.

The manufacturer MECTRON cannot be held liable, expressly or implicitly, for any type of personal injury and/or damage to property, inflicted by the user of the product and its components/accessories and which occurs in the following cases, by way of example and not of limitation:

- Misuse or use during procedures other than those specified in the destination of use of the product;
- The environmental conditions for preservation and storage of the device are not complying with the requirements indicated in *Chapter 14 on page 61*;
- The device is not used in accordance with all the instructions and requirements outlined in this manual;
- The electrical system of the facilities where the device is used does not comply with the electrical code compliance standards in force and the relative electrical safety precautions;

- Assembly operations, extensions, re-adjustments, upgrades and repairs of the device are carried out by personnel not authorised by MECTRON;
- Misuse, abuse, abnormal use, negligent use, intentional misconduct or use exceeding the limits of the device indicated and allowed and/or normal wear or deterioration, ill-treatment and/or incorrect interventions;
- Any attempt to tamper with or modify the devices under any circumstances;
- Shortage of spare materials (handpiece, wrenches) to be used in case of faults or problems;
- Use of non-original MECTRON subgingival nozzles.
- Incorrect/omitted maintenance compared to what is stated in *Chapter 12 on page 50*;
- Breach of the requirements and the information contained in *Chapter 5 on page 16*;
- Unauthorised repairs in accordance with the indications contained in *Chapter 15.4 on page 75*.

1.4 Safety Requirements

WARNING: Interference from other equipment.

An electrosurgical scalpel or other electrosurgical units positioned near the device turbodent touch may interfere with the correct operation of the device.

 **WARNING: Interference with other equipment.** Though compliant with standard IEC 60601-1-2, turbodent touch and its accessories may nonetheless interfere with other devices nearby. turbodent touch and its accessories must not be used near to or stacked on other devices. However, if this were to prove necessary, you must check and monitor correct operation of the device in that configuration.

 **WARNING: Risk of explosion.** The device cannot operate in environments where the atmosphere is saturated with flammable gases (anaesthetic mixtures, oxygen, etc.).

 **CAUTION:** If the end user, operating in their dental practice or clinic, is required to conduct periodic checks to comply with regulatory requirements, the testing procedures for evaluating the safety of medical devices and systems must be carried out in accordance with the IEC 62353 standard: "Medical electrical equipment – periodic inspections and tests to be carried out after repair of medical electrical equipment". The frequency of periodic inspections in the intended conditions of use described in this manual

is once per year or every 2000 hours of use, whichever condition is satisfied first.

⚠ CAUTION: Maintenance and/or repair work must only be carried out by specialised personnel.

⚠ WARNING: Checking the condition of the device and its accessories before treatment.

Always check that there is no water under the device. Before each treatment, always check the perfect working order of the device and the efficiency of the components/accessories. If operation anomalies are observed do not carry out the treatment. If the faults concern the device, contact a MECTRON Authorised Service Centre.

⚠ CAUTION: The electrical system located on the premises in which the device is installed and used must be compliant to the electrical code compliance standards in force and the relative electrical safety precautions.

⚠ CAUTION: To avoid the risk of electrical shocks, this device must be strictly connected to a grounded power supply.

⚠ WARNING: Sterile single-use material.

Before using a sterile disposable object, check the integrity of the packaging to ensure that it is sterile. The object will no longer be sterile if the packaging is broken or damaged.

⚠ WARNING: Cleaning and sterilising new or repaired instruments. All components of the new or repaired devices, with the exception of SUBGINGIVAL PERIO TIPS, are not sterile. When first used and after each treatment, they must be cleaned and sterilised in strict accordance with the instructions in *Chapters 6, 7, 8, 9*.

⚠ WARNING: Infection control. For maximum patient and operator safety, before using all reusable parts and components/accessories, make sure they are previously cleaned and sterilised according to the instructions (see *Chapters 6, 7, 8, 9*).

⚠ CAUTION: After having autoclaved the handpiece, or any other sterilisable component/accessory, wait for them to completely cool down before reusing them.

⚠ WARNING: Use only original MECTRON accessories and replacement parts.

⚠ CAUTION: No modification of this device and its components/accessories is permitted.

⚠ WARNING: Jet polisher. Patients wearing contact lenses must remove them prior to receiving treatment with the jet polisher.

⚠ WARNING: Contraindications - Supragingival powder jet polisher. Do not aim the air, water, or supragingival powder jet onto the soft tissues or inside the gingival sulcus. Failure to comply with this prescription can cause a gingival tissue emphysema (emphysema of the mucous and/or subcutaneous). For this type of application, use only subgingival powder.

⚠ WARNING: Contraindications - Jet polisher. Do not use the device and its accessories near areas subject to recent dental extraction and in traumatised/damaged areas (or areas nearby) due to the risk of emphysema.

⚠ WARNING: Temperature of the water spray - Jet polisher. The device is equipped with a double safety mechanism that controls the temperature of the water spray. Before treatment, it is however recommended to instruct the patient to inform the operator if he perceives an excessive increase in the temperature of the water.

⚠ WARNING: In case of an adverse event and/or accident attributable to the device and/or related accessory during correct use and in accordance with the intended use, a report must be made to the Competent Authority and manufacturer indicated on the product label.

1.5 Symbols

Symbol	Description	Symbol	Description
	Class IIa device compliant with Regulation (EU) 2017/745		Nemko brand Compliance with UL - CSA regulations
	Importer		Distributor
	Manufacturer		Date of manufacture
	Type B applied part		Medical device
	Unique Device Identifier		Health Industry Bar Code
	Batch number		Serial Number
	Model number		Catalogue number
	Caution		Consult instructions for use or consult electronic instruction for use
	Sterilisable up to a max. temperature of 135 °C		Non-sterile
	Alternating current		Quantity in package: 1
	Power switch set to "on"		Power switch set to "off"
	Temperature limit		Foot switch connection
	Humidity limitation		Atmospheric pressure limitation
	Separate collection for waste of electrical and electronic equipment		This side up
	Fragile		Keep dry

Symbol	Description	Symbol	Description
	General warning sign ^{a)}		Warning: electricity
	Wear protective glasses ^{b)}		Wear protective mask ^{b)}
	Handpiece water flow		Handpiece air flow
IP**	International Protection Code (IP) for mechanical enclosure		Authorized Representative for the indicated Country or Jurisdiction (XX)
Rx only	For US market only CAUTION: Federal (US) law restricts the marketing of this device to prescription or licensed practitioners		

a) The symbol is represented by a yellow warning triangle and a black graphical symbol.

b) The symbol is represented by a blue warning circle and a white graphical symbol.

The following table presents the additional symbols applicable to sterile devices and/or devices with an expiry date, used in combination with the Turbodont Touch medical equipment and its associated accessories.

Symbol	Description	Symbol	Description
	Sterilised with Ethylene Oxide (EO)		Do not re-sterilise
	Do not reuse		Use by date
	Single sterile barrier system		Do not use if the packaging is damaged and consult instruction for use

NOTE: For other symbols, refer to Chapter 15.1 on page 70

2 IDENTIFICATION DATA


A proper description of the model and serial number of the device and its accessories will enable the After-Sales Service to provide

quick and effective answers. Always provide this information every time you contact the MECTRON After-Sales Service.

2.1 Device Nameplate

Each turbodent touch comes with a nameplate indicating the main technical features and traceability data, including the UDI. This nameplate is located under the device. The complete technical specifications are provided in *Chapter 14* on page 61.


NOTE: The complete list of symbols and their description is provided in *Chapter 1.5* on page 6.

2.2 Accessories Identification Data

Traceability data, including the UDI code, are laser-engraved on each AIR-POLISHING handpiece.

NOTE: The complete list of symbols and their description is provided in *Chapter 1.5* on page 6.

3 DELIVERY

3.1 List of components

Please refer to the figure below. turbodent touch is supplied with standard equipment and a set of components/accessories that can be ordered separately, which vary depending on the configuration and client requests (see Tables on page 11).

NOTE: Both items included in the standard supply and all components can be ordered separately by the client.

The device's packaging must be handled with care to prevent strong impacts, as it contains sensitive electronic components. Both transport and storage should be conducted with caution.

Do not stack various boxes in order to avoid pinching the underlying packaging.

All material sent by MECTRON is checked at the time of dispatch.


The device and its components are shipped properly protected and packaged.

Upon receipt, check if the device was damaged during transport and, if damage/defects are found, submit a claim to the carrier.

Keep the packaging for possible return to a MECTRON Authorised Service Centre.

WARNING: Before starting the treatment, always make sure to have spare materials (AIR-POLISHING handpieces, SUBGINGIVAL PERIO TIPS, K9 wrench) to use in the event of failures or drawbacks.

Components and Accessories for turbodent touch		
Ref	Item	Notes
A	Device Body	
B	AIR-POLISHING 120°	120° air-polishing handpiece
C	Foot pedal FS-11	
D1	Prophy Powder Container	
D2	Perio Powder Container	
E	Powder container cap	
F	Use and Maintenance Manual	
G	Power cable	Plug compatible with destination country
H	Air supply tubing with quick coupling	
I	Water supply tubing with quick coupling	
L	Water filter	
M	Water circuit disinfection kit	
N	Cleaning needle kit	
O	O-ring kit turbodent touch	
P	K9 Wrench	
Q	SUBGINGIVAL PERIO TIPS	40 pc set of SUBGINGIVAL PERIO TIPS
R	AIR-POLISHING PERIO	Perio air-polishing handpiece
S	AIR-POLISHING 90°	90° air-polishing handpiece

4 INSTALLATION

4.1 First Installation

The device and its accessories must be installed in a suitable place allowing its comfortable use.

⚠ WARNING: The place where the device is installed must meet the requirements set out in *Chapter 4.2 on page 12*.

turbodent touch can be purchased ready for use or may need to be enabled by typing in an activation key.

If your device requires an activation key, the procedures may vary from country to country. Always refer to your dealer and the leaflet included in the packaging for further details.

4.2 Safety Requirements During Installation

⚠ WARNING: Contraindications. Interference with other equipment.
Though compliant with standard IEC 60601-1-2, turbodent touch and its accessories may nonetheless interfere with other devices nearby. turbodent touch and its accessories must not be used near to or stacked on other devices. However, if it were to prove necessary, you must check and monitor correct operation of the device in that configuration.

⚠ WARNING: Contraindications. Interference from other equipment.
An electrosurgical scalpel or other electrosurgical units positioned near the device turbodent touch and/or its related accessories may interfere with the correct operation of the device.

⚠ CAUTION: The electrical system located on the premises in which the device is installed and used must be compliant to the electrical code compliance standards in force and the relative electrical safety precautions.

⚠ CAUTION: To avoid the risk of electrical shocks, this device must be strictly connected to a grounded power supply.

⚠ CAUTION: The device shall be positioned such that the power plug remains readily accessible at all times, as it serves as the means of disconnection from the mains supply.

⚠ WARNING: Risk of explosion.
The device cannot operate in environments where there are atmospheres saturated with flammable gases (aesthetic mixtures, oxygen, etc.).

⚠ WARNING: Install the device in a place protected against impact or accidental splashing of water or liquids.

⚠ WARNING: Do not install the device above or near heat sources. Ensure the installation allows adequate air circulation around the device.

⚠ CAUTION: Do not expose the device and its related accessories to direct sunlight or UV light sources.


⚠ CAUTION: The device and its related accessories is portable but must be handled with care when moved. Position the foot pedal on the floor so that it can only be activated intentionally by the operator.

4.3 Connecting the components

Connect the foot pedal to the back of the device in the socket marked with the symbol  by inserting the foot pedal plug until you hear a 'click'.


CAUTION: Exercise caution when positioning the foot pedal, ensuring it can only be activated intentionally by the operator.

1


Drain the condensate from the compressed air system. Connect the supply hose to the compressed air circuit in the medical practice using a suitable reduction and shut-off valve (not included in the MECTRON supply). Connect the quick coupling to the male coupling on the back of the device.

2


Connect the water supply hose to the water circuit in the medical practice using a suitable reduction and shut-off valve (not included in the MECTRON supply). Connect the quick coupling to the male coupling on the back of the device.

3

Plug the power cable into its socket at the back of the device and then to the wall socket.

4

CAUTION: The device shall be positioned such that the power plug remains readily accessible at all times, as it serves as the means of disconnection from the mains supply.

Unscrew the cap of the prophy powder container on the left.

⚠ CAUTION: The prophy powder container is marked "prophy" and is positioned on the left.

⚠ CAUTION: Before extracting the powder container or unscrewing the caps, check that 'refill' function has been performed and that the corresponding LED remains in pulsating mode or that the device is turned off (see *Chapter 5.5 on page 24*).

Pour the specific supragingival powder distributed by MECTRON for turbodent touch into the container, making sure the level stays below the diffuser located inside.


⚠ WARNING: Insert only specific supragingival powder distributed by MECTRON for turbodent touch in the prophy powder container.

⚠ CAUTION: Correct level of powder in the container.



Minimum level: The level of powder in the container must not be lower than 10 mm to prevent cleaning performance from dropping.

Maximum level: The level of the powder in the container must remain below the diffuser (at least 5 mm).

5


6

Screw the cap onto the container as far as it will go.

⚠ CAUTION: Always check that the threaded parts of the powder containers are perfectly clean.


7

Unscrew the cap of the perio powder container on the right.

⚠ CAUTION: The perio powder container is marked "perio" and is positioned on the right.


8

Pour the specific subgingival powder distributed by MECTRON for turbodent touch into the container, making sure the level stays below the diffuser located inside.

⚠ WARNING: Insert only specific subgingival powder distributed by MECTRON for turbodent touch in the perio powder container.


9

⚠ CAUTION: Correct level of powder in the container.

Minimum level: The level of powder in the container must not be lower than 10 mm to prevent cleaning performance from dropping.


Maximum level: The level of the powder in the container must remain below the diffuser (at least 5 mm).

Close the tap on the container without tightening it excessively.

CAUTION: Always check that the threaded parts of the powder containers are perfectly clean.

10

5 USE

5.1 Power On/Off

Switching on the device

The switch is positioned at the rear of the device, on the left.

Flick the switch into the "I" position, taking care not to press the foot pedal.

All the warning lights on the device will turn on and off again. The touch surface stays off for a few moments until an acoustic signal indicates the end of the diagnostic cycle. At the end of the diagnostic cycle, the system loads the preset configuration and is ready for use.

Preset configuration:


- "perio" function,
- "low" level

NOTE: The settings can only be changed when the foot pedal is released.

Switching off the device

The switch is positioned at the rear of the device, on the left.

Flick the switch into the "O" position, taking care not to press the foot pedal.


CAUTION: Place the device in such a way that the power plug is easily accessible as it is considered a disconnection device.

5.2 Description of the Touch Surface

The user can configure the device by simply touching the touch surface.

NOTE: An audible signal confirms the button selection. A short audible signal indicates that the button pressed cannot be selected in the current configuration.

NOTE: All LEDs indicating the selection of different functions and power levels are green. Error notifications are represented by a wrench symbol with a yellow LED.



"prophy" and "perio" FUNCTIONS:

Depending on the type of application, one of two treatments can be selected as follows:

- **"prophy"**: dedicated to the clinical indications of the specific supragingival powder distributed by MECTRON for turbodent touch;
- **"perio"**: dedicated to the clinical indications of the specific subgingival powder distributed by MECTRON for turbodent touch.

Different power levels can be selected for the two types of treatment using the dedicated buttons (high, medium, low).


"no powder" FUNCTIONS

The "no powder" functions can be used to disable the emission of powder from the handpiece.


- The key marked with the symbols:   disables the emission of powder from the handpiece, allowing only water and air to be dispensed.
- The key marked with the symbol:  disables the emission of powder and water from the handpiece, allowing only air to be dispensed.

NOTE: Both functions disable the emission of powder from the handpiece.


"refill" FUNCTION

The "refill" function must be used to depressurise the powder containers, in order that they can be opened or removed from the device (see *Chapter 5.5* on page 24).

IRRIGATION

The water circuit flow rate can be regulated continuously through the right hand knob;

“flush” FUNCTION



The “flush” function allows a rinsing cycle of the irrigation circuit and empties both channels simultaneously.

The “flush” function can be used to change the type of powder, cleaning the circuit of the previously used type. (See Chapter 5.6 on page 25)

NOTE: The “flush” function can be activated by lifting the handpiece and checking that the powder containers are correctly inserted into the device.

SYMBOLS

turbodent touch is equipped with a diagnostic circuit allowing the detection and display of any errors on the touch surface via a wrench-shaped symbol (see Chapter 15.1 on page 70).

EN

5.3 Safety Requirements Before and During Use

⚠ WARNING: Use only original MECTRON accessories and replacement parts.

⚠ WARNING: Before starting work, always ensure you have spare materials (handpiece, wrenches) to be used in case of faults or problems.

⚠ WARNING: Checking the device and its accessories before treatment. Always check that there is no water under the device. Before each treatment, always check the perfect working order of the device and the efficiency of the components/accessories. If operation anomalies are observed do not carry out the treatment. Contact a MECTRON Authorised Service Centre with regard to errors concerning the device and/or its related accessories.

⚠ WARNING: Infection control. First use: All reusable parts and components

(new or returned from a MECTRON Authorised Service Centre) are delivered in NON-STERILE conditions and must be treated before each use, according to the instructions provided in Chapters 6, 7, 8, 9.

Subsequent uses: After each treatment, clean and sterilise all reusable parts and components according to the instructions in Chapters 6, 7, 8, 9.

⚠ WARNING: The patient must not come into contact with the device core unit or with the foot pedal.

⚠ WARNING: During the intervention on the patient, do not perform any maintenance tasks on the system.

⚠ WARNING: Personal injury. Ensure the cables do not obstruct the free movement of personnel.

⚠ WARNING: Contraindications. Patients on a sodium restricted diet must not

receive treatment with sodium bicarbonate powders.

 WARNING: Contraindications.

Patients suffering from severe respiratory problems, such as chronic bronchitis, asthma, emphysema, etc. must not receive prophylaxis treatment, unless otherwise specified by the surgeon.

 **CAUTION:** Before using powders, always refer to the Instructions for Use specific to the product in use. Carefully follow the manufacturer's instructions to ensure safe and effective use.

 **WARNING: Contraindications.** Patients wearing contact lenses or glasses must remove them prior to receiving treatment with the jet polisher.

  **WARNING: Patient Preparation**
During treatment, patients are required to wear protective glasses.

   **WARNING: Operator Preparation**
During treatment, the operator is required to wear protective glasses and a protective mask.

 **WARNING: Contraindications - Jet polisher.** Do not aim the air, water, or supragingival powder jet onto the soft tissues or inside the gingival sulcus. Failure to comply with this prescription can cause a gingival tissue emphysema (emphysema of the mucous and/or subcutaneous). For this type of application, use only subgingival powder.

 **WARNING: Contraindications.** Do not use the device near areas subject to recent dental extraction and in traumatised/damaged areas (or areas nearby) due to the risk of emphysema.

 **WARNING: Temperature of the water spray.** The device is equipped with a double safety system that monitors the temperature of the water spray. Before the treatment, it is however recommended to instruct the patient to inform the operator if he perceives an excessive increase in the temperature of the water.

 **WARNING: Infection control and cleaning the water and air circuits.** For maximum patient and operator safety, after each treatment, follow all the instructions outlined in *Chapters 6, 7, 8, 9.*

 **WARNING: Do not use the device without water.** Check that the device is connected to the water circuit and that the water tap is open.

 **WARNING: Treatment preparation.** High vacuum suction (HVS) and saliva ejector are recommended during treatment. They evacuate the powder deviated by the treated tooth.


 **CAUTION:** Do not attempt to unscrew the cap of the powder container before having performed the "refill" cycle.

5.4 Using the Handpiece

Before using the handpiece, ensure that all components are connected as described in *Chapter 4.3 on page 13*.

Lift the handpiece cord;


1

Choose the type of AIR-POLISHING handpiece depending on the operation that needs to be performed;


2

⚠ CAUTION: The SUBGINGIVAL PERIO TIPS can only be used with the AIR-POLISHING PERIO handpiece.

On the quick coupling of the air-polishing handpiece cord, slide the ring and hold it in position, then insert the air-polishing handpiece until it comes into contact, and finally release the ring again;

3


4

Select the desired "prophy" or "perio" function.

- "prophy" function: use the specific supragingival powder distributed by MECTRON for turbodent touch;
- "perio" function: use the specific subgingival powder distributed by MECTRON for turbodent touch.

⚠ WARNING: Before performing any other operations, make sure the air-polishing handpiece is properly and completely inserted in the quick coupling (the quick coupling ring and air-polishing handpiece must be in contact) and that the ring has been returned to the initial position.

⚠ WARNING: Treatment preparation. High vacuum suction (HVS) and saliva ejector are recommended during treatment. They evacuate the powder deviated by the treated tooth.


5

Insert the SUBGINGIVAL PERIO TIPS on the AIR-POLISHING PERIO handpiece, pushing it in until it comes into contact.

⚠ CAUTION: If the AIR-POLISHING PERIO handpiece has been selected, the SUBGINGIVAL PERIO TIPS must be inserted on the front part.

⚠ CAUTION: Handle the SUBGINGIVAL PERIO TIPS with care.

⚠ WARNING: Before using the SUBGINGIVAL PERIO TIPS (supplied sterile), check the integrity of the packaging and inspect the product to ensure there is no damage. Do not use the SUBGINGIVAL PERIO TIPS if the packaging is open or damaged. The SUBGINGIVAL PERIO TIPS loses sterility if the packaging is broken or compromised. If the packaging is damaged, dispose of the tip properly. Do not re-sterilize or reuse the tip.

⚠ CAUTION: With the AIR-POLISHING PERIO handpiece, only specific subgingival powder distributed by MECTRON for turbodent touch must be used.

⚠ WARNING: Make sure the SUBGINGIVAL PERIO TIPS is properly inserted on the air-polishing handpiece; the two pieces must be in contact.

⚠ CAUTION: Use only original MECTRON components.

Press the foot pedal to start treatment. The flow of water can be regulated via the right hand knob until the desired amount is reached.

⚠ ⚙ WARNING: Patient Preparation During treatment, patients are required to wear protective glasses.



⚠ ⚙ ⚙ WARNING: Operator Preparation During treatment, the operator is required to wear protective glasses and a protective mask.

If the user switches from “perio” to “prophy” or vice versa, the device automatically performs the circuit emptying process, which lasts for about 3 seconds.


NOTE: During the transition from the “perio” function to “prophy” or vice versa, do not point the handpiece toward the patient or the operator.

At the end of the treatment, stow the air-polishing handpiece in its housing.

6

7

5.5 “refill” function


The “refill” function allows the powder containers to be depressurised, enabling their opening or removal and preventing powder from escaping. This function must be used whenever a container is to be reloaded or cleaned.

NOTE: The powder containers are only pressurised when, after being selected, the foot pedal is pressed.

CAUTION: During the “refill” cycle, air and powder are expelled from the air-polishing handpiece.

Press the “refill” button.

1


2

Wait for the end of the cycle, during which the “refill” key LED will flash together with that of the relative powder container, emitting an audible signal.

NOTE: If both the powder containers are pressurised, the “refill” cycle will depressurise both.


The “refill” function will depressurise the containers one at a time, and not simultaneously.

NOTE: The duration of the “refill” cycle varies depending on whether one or both pressurised powder containers are present.

3


At the end of the cycle, the “refill” key LED will be lit and will start to pulse rather than blink, indicating that the containers can now be opened or removed.

At this point the containers can simply be closed or re-inserted to continue using the device.

NOTE: If the powder containers are neither opened nor removed, it suffices to simply press the “refill” key to return to the last used program.

4

EN

5.6 “flush” function

The “flush” function is used to perform the irrigation circuit rinse cycle and to empty both channels simultaneously.


⚠ CAUTION: Flush function.

The “flush” function must be used after each treatment, before starting the cleaning and sterilisation procedures.

Lift the AIR-POLISHING handpiece;


⚠ **CAUTION:** Make sure that both powder containers are inserted in the device before activating the “flush” function.

1

2


Place the cord and/or handpiece over a container or sink to catch the liquid and powder expelled during the “flush” cycle;

To enter "flush" mode, select "flush" on the touch surface: the corresponding LED will blink.


All the other selectable options on the touch surface will be disabled.

3

Press the foot pedal once and release it to start the cycle. The LEDs of the function and powder containers will flash. A sequence of short acoustic signals indicates that the cycle is running;

4


NOTE: Before running the "flush" cycle, check that the right-hand tap is open.

CAUTION: The execution of the "flush" cycle cannot be interrupted even by returning the handpiece to its holder.


The "flush" cycle is run on both powder containers simultaneously. Once the cycle is finished, the "flush" LED will continue to blink until any key is pressed.

5

Once the irrigation circuit has been subjected to the "flush" cycle, proceed to disassemble, clean and sterilise the individual parts (see Chapters 6, 7, 8, 9).

6

5.7 Stand-by Mode (Energy Saving)

After 20 minutes of inactivity, the device automatically enters stand-by mode (energy saving).

In this phase, the system disables water heating and stops the continuous air flow of the cleaning handpiece.

The keypad remains active and retains the configured settings. The device is immediately reactivated by pressing any key.

6 CLEANING

⚠ CAUTION: To disinfect the device and/or its components/accessories, it is recommended to use water-based disinfectant solutions with neutral pH (pH 7).

DO NOT USE the following as disinfectants:

- Alcohol-based products above 70%;
- Products containing peracetic acid, formaldehyde, glutaraldehyde or other equivalent solutions;
- Very alkaline products (pH > 9);
- Products containing sodium hypochlorite;
- Products containing hydrogen peroxide;
- Products containing abrasive substances;
- Very acid products (pH < 4);
- Products containing aldehyde, amine and/or phenols;
- Acetone;
- Methyl ethyl ketone;

as they may discolour and/or damage the materials of the device and its components.

The manufacturer disclaims all liability for any damage caused by the substances mentioned above. In case of damage caused by those substances, the Warranty will be void.

⚠ CAUTION: Methods not covered in the following sections must not be used.

⚠ WARNING: The sterile single-use SUBGINGIVAL PERIO TIPS must be used for one treatment and one patient only. The single-use SUBGINGIVAL PERIO TIPS must not be reused. Sort and dispose of each single-use SUBGINGIVAL PERIO TIPS in accordance with current hospital waste disposal laws.

NOTE: Repeated reconditioning has a minimal effect on the devices and components. The end-of-life of the devices and components is generally determined by wear or damage resulting from use. MECTRON guarantees the integrity of its air-polishing handpieces for up to 250 reconditioning cycles.

6.1 Disassembling Parts for Cleaning and Sterilisation


Before performing the cleaning procedures outlined in the following sections, disconnect all components of turbodent touch.

» PROCEDURE

⚠ WARNING: Switch the device off.

Always switch the device off using the switch and disconnect the power cable from the wall socket and device core unit before carrying out cleaning and sterilisation tasks.

1


2

Disconnect the foot pedal from the device:

grasp the foot pedal connector, press the release tab and pull the connector back;


⚠ CAUTION: Do not try to unscrew or turn the connector during disconnection: this might damage the connector.

⚠ CAUTION: To disconnect the cable from the foot pedal, only ever grab the cable connector. Never pull the cable itself.


3

On the quick coupling of the air-polishing handpiece cord, slide the ring and hold it in place, pull out the air-polishing handpiece, then release the handpiece ring;

If the AIR-POLISHING PERIO handpiece was used, remove the SUBGINGIVAL PERIO TIPS using the K9 wrench and proceed with disposal (see *Chapter 13* on page 60);


NOTE: Insert the K9 wrench exactly in the position shown in the figure.

⚠ CAUTION: Before extracting the powder container or unscrewing the caps, check that the device is switched off and disconnected from the power mains, or start the "refill" function (see *Chapter 5.5* on page 24).


Extract the powder container from the device, remove the cap and empty it. Proceed with polishing. (see *Chapter 12.7* on page 58)

Repeat the procedure on the device's second powder container if both have been used.

Once cleaning is complete, place the containers back in the device.

⚠ CAUTION: Insert the powder containers so that the groove in the bottom of the device coincides with the curvature of the powder container.

6.2 Preparation

1. Run the "flush" function (see Chapter 5. 6 on page 26);
2. Disinfect the irrigation circuit as described in Chapter 10 on page 44.
3. Check that all of the following components have been removed/disconnected from the device core unit (see Chapter 6.1 on page 28):
 - Power supply cable;
 - Foot pedal;
 - AIR-POLISHING handpiece;
 - Water and air hose.

⚠ CAUTION: The cleaning and sterilising operations described in the following sections are to be performed at first use and at all subsequent uses.

⚠ WARNING: Always switch the device off using the O/I switch and disconnect it from the power mains before carrying out the post-preparation cleaning tasks.

⚠ CAUTION: Do not immerse the handpiece in non-recommended solutions as it may be damaged.

⚠ CAUTION: In case of excessive powder humidity, remove/disconnect the AIR-POLISHING handpiece from the quick coupling on the cord, free the AIR-POLISHING handpiece channel of any residual powder using the supplied Ø 0,4 mm cleaning needle, then run the "flush" function without the air-polishing handpiece.

Machine core unit and foot pedal:

The only applicable cleaning method is manual cleaning.

Handpiece, K9 wrench and irrigation circuit cleaning kit:

The operator can choose between the following cleaning procedures:

- Manual cleaning;
- Automatic cleaning.

The procedures applicable to the turbodent touch components in this manual have been independently validated.

6.3 Cleaning the Non-Sterilisable Parts

The following procedure must be performed on all non-sterilisable parts of the device. The parts concerned are:

- Device core unit;
- Foot pedal and relative connection cable to the machine
- Cleaning cord and its quick coupling for the handpiece

⚠ WARNING: Always switch the device off using the O/I switch and disconnect it from the power mains before carrying out cleaning tasks.


» NECESSARY MATERIAL

- Clean, soft, low-lint cloths;
- Cleaning solution (pH 4-9).

» PROCEDURE


Clean the surface of the parts in question using a clean, dry, low-lint cloth, dampened with a detergent solution (pH 4-9);

1

Dry the parts using a dry, non-abrasive, low-lint cloth.

2

⚠ CAUTION: Do not sterilise the non-sterilisable parts. They may stop working and cause damage to people and/or property.

⚠ WARNING: The device and its non-sterilisable parts are not protected against the penetration of liquids. Do not spray

liquids directly onto the surface of the device and its parts sold as non-sterilisable.

⚠ CAUTION: Do not use running water to clean the non-sterilisable parts.

⚠ CAUTION: Do not use solutions other than those recommended as they may damage the non-sterilisable parts.

6.4 Cleaning Sterilisable components

The Operator can choose between the following cleaning procedures:

- Manual cleaning;
- Automatic cleaning.

The procedures applicable to the turbodent touch components have been validated by an independent body.

The sterilisable parts of the device are:

- AIR-POLISHING handpiece;
- K9 Wrench;
- Irrigation circuit cleaning kit.

NOTE: The cleaning procedures must be performed immediately after each use. Immerse the instrument in demineralised water or in an enzymatic detergent solution immediately after use. Do not leave residue or blood deposits on the instruments; eliminate larger impurities with a disposable cloth or paper towel.

CAUTION: The instructions provided below have been validated by the manufacturer to allow the reuse of the medical device.

The process manager is responsible for ensuring that the processes repeated are effectively performed using the equipment, materials and staff in the reprocessing structure in order to obtain the desired result. This generally requires the validation and systematic monitoring of the process.


Similarly, all deviations from the instructions provided by the processes manager must be adequately assessed to judge their efficiency and potential undesired consequences.

6.4.1 Pre-Cleaning

The “pre-wash” cycle is optional and preparatory to the chosen cleaning method. After following the instructions provided in *Chapter 6.2 on page 30*, proceed as follows:

Thoroughly clean the outer surfaces with ready-to-use aldehyde-free wipes (with less than 35% alcohol) until visually clean. Ensure the surfaces are sufficiently damp. Note the action time of the detergent as indicated by the manufacturer.

1

6.4.2 Manual Cleaning

Manual cleaning can be carried out as an alternative to the automatic cleaning described in Chapter 6.4.3 on page 39.

» NECESSARY MATERIAL

- Enzymatic detergent at pH 4-9;
- Water;
- Container for immersion in the enzymatic liquid;
- Ultrasonic tank;
- Clean, soft, low-lint cloths;
- Brush with soft nylon bristles;
- Syringe;
- Demineralised water


» PROCEDURE - AIR-POLISHING HANDPIECES

Prepare a solution of enzymatic detergent^{a)} with pH 4-9, according to the manufacturer's instructions;

⚠ CAUTION: Once used, dispose of the enzymatic detergent correctly, do not recycle it.

a) Process validated by an independent body with enzymatic ENZYMEC, 0,8% v/v.


1

2


Blow compressed air into the central hole of the handpiece, from the front side as shown in the figure.

⚠ CAUTION: If the single-use SUBGINGIVAL PERIO TIPS is on the AIR-POLISHING handpiece, remove and dispose of it.


3

Place the air-polishing handpiece in a clean container, in a horizontal position, and add enough enzymatic detergent solution to completely cover the handpiece;


Leave to soak for 10 minutes at $22^{\circ}\text{C} \pm 2^{\circ}\text{C}$ ($71,6^{\circ}\text{F} \pm 3,6^{\circ}\text{F}$). This procedure reduces the amount of blood, protein and mucous present on the AIR-POLISHING handpiece;

4


After 10 minutes of immersion in the enzymatic solution, gently brush the surface of the AIR-POLISHING handpiece using the brush with soft nylon bristles for at least 20 seconds, and eliminate all visible residue. Thoroughly clean the hard-to-reach areas such as the edges, recesses, and joints;

5


Rinse the AIR-POLISHING handpiece under hot running water ($40^{\circ}\text{C} \pm 5^{\circ}\text{C}$ ($104^{\circ}\text{F} \pm 9^{\circ}\text{F}$)) for about 2 minutes, moving it slightly to allow the water to reach the entire surface;

6

Place the AIR-POLISHING handpiece in the ultrasonic tank immersed in the enzymatic detergent solution for at least 10 minutes at $22^{\circ}\text{C} \pm 2^{\circ}\text{C}$ ($71,6^{\circ}\text{F} \pm 3,6^{\circ}\text{F}$);


7

8


Gently brush the surface of the AIR-POLISHING handpiece using the brush with soft nylon bristles;

Rinse the AIR-POLISHING handpiece under hot running water ($40^{\circ}\text{C} \pm 5^{\circ}\text{C}$ ($104^{\circ}\text{F} \pm 9^{\circ}\text{F}$)) for about 2 minutes, moving it slightly to allow the water to reach the entire surface;


9

Place the AIR-POLISHING handpiece in the ultrasonic tank immersed in demineralised water for at least 10 minutes at $40^{\circ}\text{C} \pm 2^{\circ}\text{C}$ ($104^{\circ}\text{F} \pm 3,6^{\circ}\text{F}$);

10

After 10 minutes of immersion in the ultrasonic tank with demineralised water, rinse the internal channels of the AIR-POLISHING handpiece by injecting 20 ml of demineralised water at ambient temperature ($20^{\circ}\text{C} - 25^{\circ}\text{C}$ ($68^{\circ}\text{F} - 77^{\circ}\text{F}$)) using a syringe with needle, five times in a row.

» PROCEDURE - K9 WRENCH AND IRRIGATION CIRCUIT CLEANING KIT

Prepare a solution of enzymatic detergent^{a)} with pH 4-9, according to the manufacturer's instructions;

CAUTION: Once used, dispose of the enzymatic detergent correctly, do not recycle it.

a) Process validated by an independent body with enzymatic ENZYMEC, 0,8% v/v.

Place the wrench in a clean container. Add enough enzymatic detergent solution to completely cover it.

Leave the wrench to soak in the enzymatic detergent solution for 10 minutes at ambient temperature $22\text{ °C} \pm 2\text{ °C}$ ($71,6\text{ °F} \pm 3,6\text{ °F}$);

After soaking in the enzymatic solution, gently brush all the surfaces until all visible soil has been removed.


Use a clean brush with soft nylon bristles for the external surfaces, and a clean pipe cleaner with soft nylon bristles for the internal cavities and grooves.

NOTE: Thoroughly brush all parts of the wrench for about 20 seconds.


Remove the wrench from the enzymatic detergent solution.

Thoroughly rinse and brush all the surfaces of the wrench (see previous point) under running water for at least 2 minutes;


11


12

13


14


Place the wrench in the ultrasonic tank immersed in the enzymatic detergent solution at $24^{\circ}\text{C} \pm 2^{\circ}\text{C}$ ($75,2^{\circ}\text{F} \pm 3,6^{\circ}\text{F}$) and run a cycle for at least 20 minutes;

15


Remove the wrench from the ultrasonic tank and rinse it under running water so as to eliminate all detergent residue. Brush the internal and external surfaces of the wrench with a clean brush with soft nylon bristles, under running water.

16


Place the cleaning kit in a clean container. Add enough enzymatic detergent solution to completely cover it. Leave it to soak in the enzymatic detergent solution for 10 minutes at ambient temperature ($22^{\circ}\text{C} \pm 2^{\circ}\text{C}$ ($71,6^{\circ}\text{F} \pm 3,6^{\circ}\text{F}$)).

17


Remove the tube and fitting from the enzyme detergent solution. Thoroughly rinse all the surfaces under running water for at least 2 minutes.

18


Place the cleaning kit in the ultrasonic tank immersed in the enzymatic detergent solution at $24\text{ }^{\circ}\text{C} \pm 2\text{ }^{\circ}\text{C}$ ($75,2\text{ }^{\circ}\text{F} \pm 3,6\text{ }^{\circ}\text{F}$) and run a cycle for at least 10 minutes.

19


Remove the cleaning kit from the enzymatic solution and rinse it under warm running water ($40\text{ }^{\circ}\text{C} \pm 5\text{ }^{\circ}\text{C}$ ($104\text{ }^{\circ}\text{F} \pm 9\text{ }^{\circ}\text{F}$)) for at least 2 minutes.

20

Rinse the inner channel with a 20 ml syringe (without needle) previously filled with demineralised water. Repeat three times.

21

6.4.3 Automatic Cleaning

Automatic cleaning can be carried out as an alternative to the manual cleaning described in Chapter 6.4.2 on page 33.

NOTE: Procedure validated with:

- Miele PG8536 washer/disinfector;
- Miele DES-VAR-TD program;

- Alkaline detergent: neodisher® FA (0.2% v/v);
- Neutralising liquid: neodisher® Z (0.1% v/v)

» NECESSARY MATERIAL

- Alkaline detergent: neodisher® FA (0.2% v/v);
- Neutralising liquid: neodisher® Z (0.1% v/v);
- Water;
- Metal basket;
- Adapters;
- Thermal-disinfector.

NOTE: Ensure that the components/accessories are properly secured in the basket and cannot move during washing. Any blows could damage them. Position the instruments in a way that the water can flow through all the surfaces, even internal.

⚠ WARNING: Avoid overloading

the thermal-disinfector as this could compromise the effectiveness of cleaning.

⚠ WARNING: On completion of the cleaning cycle in the washer-disinfector, the handpiece and the K9 key remain at the wash temperature for a long time. Use appropriate precautions when extracting the handpiece and the K9 key from the washer-disinfector to prevent injury to the operator.

⚠ CAUTION: Due to its form, the handpiece can roll. When not in use, the handpiece must always be placed on its support.

» PROCEDURE


Place the components/accessories in a metal basket.

Sequence and parameters applicable to the cycle:

- 1 min, Rinse with cold water;
- 5 min, Wash with alkaline detergent at 55 °C $\pm$ 2 °C (131 °F $\pm$ 3,6);
- 1 min, Neutralisation with neutralising liquid and water at 32 °C $\pm$ 2 °C (89,6 °F $\pm$ 3,6 °F);
- 1 min, Rinse with water at 32 °C $\pm$ 2 °C (89,6 °F $\pm$ 3,6 °F);
- 5 min, Thermal disinfection at 93 °C (199,4 °F) with demineralised water.

Automatic thermal disinfection is not experimentally tested. In compliance with ISO 15883-1, Table B. 1[4] thermal disinfection at a temperature of 90°C (194°F) for 5 min determines a value of A0 3000.

1

7 CLEANING CHECK

» NECESSARY MATERIAL


- Light source;
- 2, 5X Magnifier.

» PROCEDURE


Once the cleaning operations have been completed, check the handpiece under an adequate source of light, if necessary using a magnifying glass 2.5X, paying attention to any parts that may conceal soil residue (threading, cavities, grooves) and, if need be, repeat the cleaning cycle if soil is still visible. Finally, check the integrity of those parts and those elements that could have deteriorated during use;

Repeat the checks for the other components, repeating the cleaning cycle if necessary.

1

2

8 DRYING

» NECESSARY MATERIAL

- Compressed air;
- Soft, low-lint cloth;

» PROCEDURE


Dry all parts of the AIR-POLISHING handpiece thoroughly by blowing compressed air from the tip of the handpiece;

⚠ CAUTION: Before starting the sterilisation cycle, make sure that the AIR-POLISHING handpiece is thoroughly dry both inside and out. To do so, blow compressed air both externally and through the internal through-holes. This will prevent the appearance of stains, streaks on the surface, or oxidation inside the AIR-POLISHING handpiece.


Dry the K9 wrench using a soft, low-lint cloth.

Dry the irrigation circuit cleaning kit.


1

2

3

4

9 STERILISATION

The sterilisable parts of the device are:

- AIR-POLISHING handpiece;
- K9 Wrench;
- Irrigation circuit cleaning kit.

The sterilisable parts are manufactured with materials able to withstand a maximum temperature of 135 °C (275 °F) for a maximum time of 20 minutes.

After individually bagging the AIR-POLISHING handpiece and the other sterilisable

components, proceed with the sterilisation process in a steam autoclave.


CAUTION: Use sterilisation pouches compliant with standard UNI EN ISO 11607-1.

CAUTION: Ensure that the pouch is large enough to hold the instrument without needing to force its closure/sealing or tearing the packaging.

» PREPARATION


Seal the AIR-POLISHING handpieces individually inside a single-use sterilisation pouch.

1


Seal the K9 wrench separately inside a single-use sterilisation pouch.

2

Individually seal the tube and the fitting separately in a single-use pouch legally marketed for steam sterilisation.

3

» **PROCEDURE**

The sterilisation process validated by MECTRON, in a steam autoclave, guarantees a SAL 10^{-6} by setting the parameters indicated below:

Procedure A:

- Cycle type: 3 times Pre-vacuum (minimum pressure 60 mBar).
- Minimum sterilization temperature: 132 °C (0 °C ÷ +3 °C) / 269.6 °F (0 °F ÷ +5.4 °F).
- Minimum sterilization time: 4 min.
- Drying time: 20 min.

Or Procedure B:

- Cycle type: 3 times Pre-vacuum (minimum pressure 60 mBar).
- Minimum sterilization temperature: 134 °C (-1 °C ÷ +1 °C) / 273.2 °F (-1.8 °F ÷ +1.8 °F).
- Minimum sterilization time: 3 min.
- Drying time: 20 min.

For the sterilization of products with moist heat, follow the guidelines specified in ISO 17665 and comply with the legal requirements in your country.


⚠ CAUTION: The device core unit and foot pedal cannot be sterilised in an autoclave.

⚠ WARNING: Infection control – Sterilisable parts – Carefully remove any residual organic matter before sterilisation.

⚠ CAUTION: Carry out the sterilisation by using only an autoclave with steam from water. Do not use any other sterilisation procedure (dry heat, irradiation, ethylene oxide, gas, low temperature plasma, etc.).

⚠ CAUTION: Do not exceed the permissible load of the steam steriliser as this could cause excess residual moisture that could impair sterilisation.

1

⚠ WARNING: On completion of the autoclave sterilisation cycle, the sterilised parts remain at the sterilisation temperature for a long time. During the removal of the sterilised parts from the autoclave adopt suitable precautionary measures in order to prevent injury to the operator.

⚠ CAUTION: Wait until the sterilised parts have completely cooled down before use.

10 DISINFECTING THE IRRIGATION CIRCUIT

Before disinfecting the irrigation circuit, check that the “flush” function has been performed (see Chapter 5.6 on page 25)

NOTE: Before starting the disinfection procedure, check that the AIR-POLISHING handpiece has been disconnected from the device (Chapter 6.1 on page 28).

To disinfect the irrigation circuit, proceed as follows:

» NECESSARY MATERIAL


- CIDEX® OPA disinfectant solution;
- Distilled water;
- Irrigation circuit disinfection kit
- 60ml syringe with eccentric Luer Slip cone, without needle.

» PROCEDURE

⚠ WARNING: Switch the device off.


Switch the device off using the switch (Ref. 1) and disconnect the power cable from the wall socket (Ref. 2) and machine core unit (Ref. 3) before carrying out cleaning and sterilisation tasks.

1


2

Disconnect the quick connector of the external water circuit from the device (see Chapter 4.3 on page 13);


Unscrew the male coupling for the external water circuit connection (Ref. A); remove the water filter (Ref. B) - (see Chapter 12.6 on page 55);

3


Check the integrity of the filter, rinse it with compressed air, making sure to remove any residual impurities;

4


Completely immerse the male coupling for the external water circuit connection in the disinfectant solution (CIDEX® OPA). Leave to soak for 12 minutes at 20 °C $\pm$ 2 °C (68 °F $\pm$ 3,6 °F).

5


Remove the male coupling for the external water circuit connection from the disinfectant solution (CIDEX® OPA) and rinse it under running water; Place the male coupling for the external water circuit connection in a clean and safe place for later use.

6

Disinfect the internal surfaces of the seat where the male coupling for the external water circuit connection will be inserted using a clean, soft, low-lint cloth dampened with the disinfectant solution (CIDEX® OPA).


7

Screw in the male coupling supplied with the "Irrigation circuit disinfection kit" in its seat until it comes into contact;


8

⚠ CAUTION: Check the integrity of the "irrigation circuit disinfection kit". Replace the "irrigation circuit disinfection kit" in case of wear or damage.


If previously disconnected (for example, for cleaning), connect the hose supplied with the kit to the male water coupling;

9

Connect the power supply cable to the device (Ref. 1) and to a wall socket (Ref. 2). Switch on the device using the switch located at the back (Ref. 3);


10

DISINFECTING THE IRRIGATION CIRCUIT


Lift the air-polishing cord and place it over a container or sink able to hold the liquid that will spill out during the procedure;

11


Fully open the knob at the front right of the device;

12


Using the syringe, aspirate 60ml of disinfectant (CIDEX® OPA), making sure there is no air inside;

13

Connect the syringe to the end of the hose previously connected to the male water coupling (see point 9 of the procedure);

14


EN

Select "prophy" or "perio" on the touch surface.

NOTE: Make sure the correspondent powder container is correctly inserted into the device and the cap is screwed on.


15

Press the foot pedal until 50 ml of disinfectant has been injected, leaving the remaining 10 ml inside the syringe.


CAUTION: Do not press too hard on the piston of the syringe. The irrigation circuit hoses are quite small and it is normal for injection to occur slowly.

16

Wait 12 minutes;


17

Press the foot pedal and inject the remaining 10 ml into the syringe.


CAUTION: Do not press too hard on the piston of the syringe. The irrigation circuit hoses are quite small and it is normal for injection to occur slowly.

18


Disconnect the syringe and repeat points 13 to 18.

19


At the end of the disinfection procedure, remove the syringe (Ref. 1), the hose (Ref. 2) and unscrew the male coupling used (Ref. 3).

20


Re-insert the previously cleaned water filter (see point 4 of the procedure) in its seat - Ref. 1 -, screw in the original and previously disinfected male coupling (see points 5 and 6 of the procedure) - Ref. 2 - and reconnect the quick coupling of the external water circuit - Ref. 3.

21

Perform the "flush" function (see *Chapter 5.6 on page 25*) three times in a row.

22

11 DECONTAMINATION OF THE IRRIGATION CIRCUIT

The device's irrigation line should be kept clean to prevent any microbial contamination. A regular cleaning and maintenance protocol should be adopted to clean and protect the device's irrigation lines. The manufacturer recommends the use of EPA-registered dental device irrigation line detergents (e.g. MicroCLEAR).

Chlorine dioxide decontamination must be carried out daily at the end of the working day, according to the following protocol:

- Perform "flush" cycles with water to remove any residual traces before disconnecting the handpieces (*Chapter 5.6 on page 25*);

- Disconnect the device from the water supply (see *Chapter 4.3 on page 13*);
- Follow *Chapter 10 on page 44* and inject irrigation line detergent into the water mains connection using a 60 cc syringe;
- Leave the germicidal agent to act in the water hoses overnight and over the weekend;
- Before the first use on the next working day, reconnect the irrigation lines and run three "flush" cycles to flush the water lines.

12 MAINTENANCE

12.1 Maintenance After Every Treatment

At the end of each treatment, proceed as follows:

Run a complete cleaning cycle on the irrigation circuits using the "flush" function (see *Chapter 5.6 on page 26*).

Immediately disassemble the various parts (see *Chapter 6.1 on page 28*) and proceed with their cleaning and sterilisation (*Chapter 6.4 on page 32*).

Clean the non-sterilisable parts (*Chapter 6.3 on page 31*).

12.2 Daily Maintenance

Regardless of the time elapsed since the last treatment and use of the device, at the end of the day proceed as follows:

1. Start the refill function if the powder containers are still pressurised (see *Chapter 5.5 on page 24*).
2. Remove and empty both powder containers (see *Chapter 6.1 on page 28*).
3. Blow compressed air into the powder container housings each time they are removed, to remove any powder residue. Do not use water or lubricants.
4. Clean the powder containers with compressed air (see *Chapter 12.7 on page 58*). Place the empty powder containers back in their position in the device.
5. Start "flush" (see *Chapter 5.6 on page 25*).
6. Clean the parts of the device that cannot be sterilised (see *Chapter 6.3 on page 31*).

12.3 Transport or Long Periods of Non-use

If the device is not used for more than 72 hours, proceed as follows:

1. Empty the powder containers;
2. Run a complete cleaning cycle on the irrigation circuit using the "flush" function (see *Chapter 5.6 on page 25*).
3. Completely disinfect the irrigation circuit (see *Chapter 10 on page 44*).
4. Drain the irrigation circuit of any residual water.
5. Remove the condensate from the air filter (see *Chapter 12.9 on page 59*).
6. Disconnect the device from the power supply and the air and water circuits.
7. Clean and dry the water filter (see *Chapter 12.6 on page 55*).
8. Follow the steps described in *Chapter 10*, filling the syringe with air instead of disinfectant this time.

9. Repeat all the steps again; it is not necessary to wait 12 minutes after the injection of air.

10. In case of long periods of non-use, store the device in its original packaging in a safe place.

Before using the device and its accessories again after a long period of non-use:


11. Run a complete cleaning cycle on the irrigation circuit using the "flush" function (see *Chapter 5. 6 on page 26*).
12. Completely disinfect the irrigation circuit again (see *Chapter 10 on page 44*).
13. Clean and sterilise the handpiece and components according to the instructions in *Chapters 6, 7, 8, 9*.

 **WARNING:** periodically check the integrity of the power supply cable; when it appears to be damaged replace it with an original MECTRON spare part.

12.4 Unclogging the Air-polishing Circuit


⚠ WARNING: Switch the device off. Switch the device off using the switch (Ref. 1) and disconnect the power cable from the wall socket (Ref. 2) and core unit of the device (Ref. 3) before carrying out the cleaning and sterilisation tasks.

1


2

Disconnect the air supply and irrigation hoses from the device;


3

Remove both powder containers from their housing;


4

On the quick coupling of the AIR-POLISHING handpiece cord, slide the ring and hold it in place, pull out the AIR-POLISHING handpiece, then release the handpiece ring;


Take the air-polishing cord connected to the device and blow compressed air into the front of the cord connector;

5

The air expelled from the powder container housings will clear the powder channel inside the device, in the event it is clogged.


6

12.5 Unclogging the AIR-POLISHING handpiece


Disconnect the AIR-POLISHING handpiece if used (see Chapter 5.4 on page 20);

1

Remove the obstruction by inserting the 0,4 mm diameter needle from the front of the AIR-POLISHING handpiece;


2

CAUTION: Clean the nozzle channel only using the Ø 0.4 mm cleaning needle supplied with the device.


Remove the obstruction by inserting the 0,4 mm diameter needle from the back of the AIR-POLISHING handpiece;

3


Now remove any residual powder by blowing compressed air into the AIR-POLISHING handpiece from the head side, as shown in the figure.

4

If the mechanical actions are not effective, leave the handpiece to soak in a 2% acetic acid solution;

5


12.6 Cleaning and/or Replacing the Water Filter and the Water Inlet O-ring

Check and clean the water filter monthly, performing the following operations:

Disconnect the water supply hose from the male coupling.
Unscrew the knurled bush of the male coupling;


⚠ CAUTION: Before cleaning and/or replacing the water filter, make sure the device is disconnected from the mains power and that the liquid container is not connected.

1

2


Extract the filter and wash it under running water to eliminate any impurities obstructing it;

3

Re-introduce the filter into its seat and screw the knurled bush back into its housing tightly until it is fully home.

NOTE: Replace the filter with a new one if it is damaged or washing is not effective.


⚠ CAUTION: Regularly check the wear condition of the O-ring and replace it when necessary. Replacement must be carried out at least once every 2 years.

Replace the worn O-ring carefully, ensuring that the surface where it is positioned is not damaged or scratched. Place the new O-ring in its seat, taking care not to damage it.

⚠ CAUTION: Do not use sharp tools or instruments to insert the new O-ring into its seat.

4

12.7 Cleaning the Powder Containers and Caps

Check cleanliness of the powder container and, in particular, the cap as powder residues in the presence of moisture could solidify and make opening and closing operations difficult.

⚠ CAUTION: Do not clean the containers with cleaning solutions, use only compressed air.

⚠ CAUTION: Always switch the device off using the O/I switch and disconnect it from the power mains before cleaning the powder containers and caps.


⚠ WARNING: Before blowing compressed air into the powder containers, make sure that they have been emptied.

Blow compressed air into the container and onto the threads of both the powder containers and the caps;


Once cleaning is complete, place the containers back on the device.

⚠ CAUTION: Insert the powder containers so that the groove in the bottom of the device coincides with the curvature of the powder container.

1

2

12.8 Replacing the Air-polishing Cord O-Rings

⚠ CAUTION: Regularly check the wear condition of the O-rings and replace them when necessary. Replacement must be carried out at least once every 6 months.

EN

Remove worn O-rings, taking care not to damage and/or scratch the surface where they are located;


NOTE: A tool is required to extract and insert the O-rings. This tool is not included in the standard supply.

Thoroughly clean and, if necessary, dry the O-ring seats by blowing compressed air before inserting the new O-rings.


Place the new O-ring in its seat, taking care not to damage it.

⚠ CAUTION: Do not use sharp tools or instruments to place the new O-ring into its seat.


1

2

3

12.9 Eliminating the Condensate


The device is equipped with an air filter that captures any impurities and condensate in the pneumatic circuit.

To prevent condensate from circulating in the device, check and empty the air filter weekly, proceeding as follows:

Place an absorbent cloth under the device to collect the condensate;

⚠ CAUTION: This task must be carried out with the device switched on in order to ensure proper pressurization of the air circuit.


1

2

With the device switched on, in the “no powder” configuration (by selecting  or ), and in a perfectly horizontal position: press the air filter vent valve at the bottom of the device until only air is expelled.

NOTE: However, it is recommended to use dry compressors and introduce a dehumidifier into the practice's pneumatic circuit.

12.10 Maintenance table

Description	Qty. per unit	Replace every:
Air-polishing cord O-ring 2,5x1,0	1	6 months
Air-polishing cord O-ring 6,0x1,0	1	6 months
Water filter	1	1 year
Powder containers	2	18 months
O-ring in blue silicone (External irrigation line)	1	2 years

13 METHODS AND PRECAUTIONS FOR DISPOSAL

 **WARNING: Hospital waste.** Treat the following objects as hospital waste:

- AIR-POLISHING handpiece, when worn or broken;
- Cleaning needles, when worn or broken;
- K9 wrench, when worn or broken;
- SUBGINGIVAL PERIO TIPS, after each use.

Disposable materials and materials that carry biological risk must be disposed of according to local regulations in force regarding hospital waste.

turbodent touch and related accessories must be disposed of and treated as waste subject to separate collection.

Failure to comply with the previous points can result in a penalty pursuant to the directive on waste electrical and electronic equipment (WEEE).

The buyer has the option to return the device to the retailer when receiving a new one for disposal. Disposal instructions are available at MECTRON for proper waste management.

14 TECHNICAL DATA

EN

Device compliant with Regulation (EU) 2017/745	Class IIa
Classification under the IEC 60601-1	I Applied part: type B
Classification under the IEC 60529	IP 20 (device) IP 22 (foot pedal model FS-11)
Essential performance	According to the standard IEC 80601-2-60 the device has no essential performance
Power Supply	100-240 V~ 50/60 Hz
Max. Power Consumption	90 VA
Fuses	Type 5 x 20 mm, T 2AL, 250V
Air-polishing functions	Can be selected via touch screen: "prophy" function - "perio" function
Irrigation	Continuous regulation via knob. Water heating via heating element.
Water supply:	Operating pressure from 1.0 to 5.0 bar. Water circuit cleaning function, air-polishing part - see <i>Chapter 5. 6 on page 26</i> . Connection via the supplied hose with quick coupling through an incorporated and removable filter.
Air supply:	Inlet pressure between 4.5 and 8.0 bar. Air circuit cleaning function - See <i>Chapter 12.9 on page 59</i> Connection via the supplied hose with quick coupling through an incorporated filter and pressure reducer.
Operating Conditions	from 10 °C to 30 °C (50 °F to 86 °F) Relative humidity from 30% to 75% Air pressure P: 700hPa/1060hPa
Transport and Storage Conditions (Excluding Powders)	from -10 °C to 60 °C (14 °F to 140 °F) Relative humidity from 10% to 90% Air pressure P: 500hPa/1060hPa
Altitude	lower than or equal to 3000 metres
Weights and dimensions	4.2 Kg 310 x 260 x 145 mm (L x W x H) ^{a)}

a) L = length ; W = width; H = height

14.1 Electromagnetic Compatibility IEC 60601-1-2

EN

⚠ WARNING: Contraindications. Interference with other equipment
Though compliant with standard IEC 60601-1-2, turbodent touch and its accessories may nonetheless interfere with other devices nearby. turbodent touch and its accessories must not be used near to or stacked on other devices. However, if this were to prove necessary, you must check and monitor correct operation of the device -and of all equipments- in that configuration.

⚠ WARNING: Portable and mobile radiocommunications equipment can influence the correct operation of the device and its related accessories.

⚠ WARNING: Contraindications. Interference by other equipment
An electrosurgical scalpel or other electrosurgical units positioned near the turbodent touch device and its related accessories may interfere with the correct operation of the device itself.

⚠ WARNING: turbodent touch with its related accessories, requires special EMC precautions and must be installed and commissioned in compliance with the EMC information in this chapter.

⚠ WARNING: The use of cables and components not supplied by MECTRON may adversely affect EMC performance.

14.2 Guide and Manufacturer’s Declaration - Electromagnetic Emissions

turbodent touch, with related accessories, is designed to operate in the electromagnetic environment specified below. The purchaser or user of turbodent touch should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment Guidance
RF Emissions CISPR 11	Group 1	turbodent touch uses RF energy only for its internal operation. Therefore, its RF emissions are very low and probably do not cause any interference with nearby electronic devices.
RF Emissions CISPR 11	Class B	turbodent touch is suitable for use in all buildings, including domestic buildings, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Emissions of fluctuations voltage/flicker IEC 61000-3-3	Compliant	

14.3 Accessible Parts of the Casing

turbodent touch, with related accessories, is designed to operate in the electromagnetic environment specified below.

The purchaser or user of turbodent touch should ensure that it is used in such an environment.

Phenomenon	Essential EMC standard or test method	Immunity test values	Electromagnetic Environment Guidance
Electrostatic discharge (ESD)	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV in air	The floor must be made of wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated EM RF fields ^{a)}	IEC 61000-4-3	3 V/m ^{f)} 80 MHz - 2,7 GHz ^{b)} 80% AM at 1 kHz ^{c)}	Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended and calculated distances from the equation applicable at the frequency of the transmitter.
RF proximity fields from wireless communication equipment	IEC 61000-4-3	See Chapter 14.5 on page 68	
Magnetic field at network frequency ^{d)}	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	Power supply frequency magnetic fields should have levels characteristic of a typical location in a commercial or hospital environment.
Proximity magnetic fields	IEC 61000-4-39	See Chapter 14.6 on page 69	Portable and mobile RF communication equipment must be used at a separation distance of at least 0,15 m from the field sources.

a) If used, the interface between the PATIENT'S physiological signal simulation and the device must be positioned within a 0,1 m radius of the vertical plane of the uniform field area in the same direction as the device .

b) The device, which intentionally receives RF electromagnetic energy for its operation, must be tested at the receiving frequency. The test can be performed with other modulation frequencies identified by the RISK MANAGEMENT PROCESS. This test evaluates the BASIC SAFETY and ESSENTIAL PERFORMANCE of an intentional receiver when

an environmental signal is in the passband. It is understood that the receiver may not receive normally during the test.

c) The test can also be conducted at other modulation frequencies identified through the RISK MANAGEMENT PROCESS.

d) Only applicable to equipment and systems with magnetically sensitive components or circuits.

e) Empty.

f) Before applying the modulation.

14.4 Guide and the Manufacturer's Declaration - Electromagnetic Immunity

14.4.1 Input A.C. Power Connection

turbodent touch, with related accessories, is designed to operate in the electromagnetic environment specified below.

The purchaser or user of turbodent touch should ensure that it is used in such an environment.

Phenomenon	Essential EMC standard or test method	Immunity test values	Electromagnetic Environment Guidance
Electrical fast transient/burst ^{l) o)}	IEC 61000-4-4	± 2 kV contact 100 KHz repetition frequency	The quality of the network voltage should be that of a typical commercial or hospital environment.
Pulses differential mode ^{b) j) o)}	IEC 61000-4-5	± 0,5 kV, ± 1 kV	The quality of the network voltage should be that of a typical commercial or hospital environment.
Pulses common mode ^{b) j) k) o)}	IEC 61000-4-5	± 0,5 kV, ± 1 kV, ± 2kV	The quality of the network voltage should be that of a typical commercial or hospital environment.
Conductive disturbances induced by RF fields ^{c) d) o)}	IEC 61000-4-6	3 V ^{m)} 0,15 MHz - 80 MHz 6 V ^{m)} in ISM bands between 0,15 MHz and 80 MHz ⁿ⁾ 80% AM at 1 KHz ^{e)}	Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended and calculated distances from the equation applicable at the frequency of the transmitter.
Voltage dips ^{f) p) r)}	IEC 61000-4-11	0% UT; 0,5 cycle ^{g)} A 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° ^{q)}	The quality of the network voltage should be that of a typical commercial or hospital environment.
		0% UT; 1 cycle ^{e)} 70% UT; 25/30 cycle ^{h)} Single phase: at 0°	
Voltage interruptions ^{f) i) o)}	IEC 61000-4-11	0% UT; 250/300 cycle ^{h)}	The quality of the network voltage should be that of a typical commercial or hospital environment.

- a) Empty.
- b) During the test, all the device cables must be connected.
- c) The calibration of the current injection terminals must be performed in a system at 150 Ω .
- d) If an ISM or amateur radio band is not present among the frequency samples, as appropriate, an additional test frequency has to be used in the ISM band or in the amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.
- e) The test can be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- f) Devices and systems with a Direct Current (DC) input power supply using AC to DC converters must be tested with a converter that complies with the MANUFACTURER'S specifications. Immunity test levels are applied to the AC power input of the converter.
- g) Only applicable to devices and systems connected to a single-phase Alternating Current (AC) power supply.
- h) For example, 10/12 means 10 periods at 50 Hz or 12 periods at 60 Hz.
- i) Devices and systems with nominal input current above 16 A / phase must be disconnected from the power supply once every 250/300 cycles at any angle and from all phases simultaneously (if applicable). Devices and systems with battery backup, after the test, must resume operation using the power supply line. For devices and systems with a nominal input current lower than 16 A, all phases must be disconnected simultaneously.
- j) Devices and systems that do not have a surge protection device in the primary power circuit can only be tested at ± 2 kV between the line(s) and the earth (common mode) and at ± 1 kV between line(s) and line(s) (differential mode).
- k) Not applicable to CLASS II devices and systems.
- l) Direct coupling must be used.
- m) R. M. S., applied before modulation.
- n) The ISM bands (industrial, scientific and medical) between 0,15 MHz and 80 MHz are 6,765 MHz and 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.
- o) Applicable to devices and systems with NOMINAL input current less than or equal to 16 A / phase and devices and systems with NOMINAL input current greater than 16 A / phase.
- p) Applicable to devices and systems with NOMINAL input current less than or equal to 16 A / phase.
- q) At some phase angles, the application of this test to devices with a transformer on the input power supply may cause the opening of an overcurrent protection device. This can occur due to the saturation of the magnetic flow in the transformer core after the voltage drop. If this happens, the device must ensure BASIC SAFETY during and after the test.
- r) For equipment and systems that have multiple voltage settings and self-regulating voltage capacity, the test must be carried out at the input supply voltage specified in Table 1 - "Power input voltages and frequencies during the tests" of standard IEC 60601-1-2:2014/AMD1:2020.

14.4.2 Points of Contact with the Patient

turbodent touch, with related accessories, is designed to operate in the electromagnetic environment specified below. The purchaser or user of turbodent touch should ensure that it is used in such an environment.

Phenomenon	Essential EMC standard or test method	Immunity test values	Electromagnetic Environment Guidance
Electrostatic discharges (ESD) ^{c)}	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV in air	The floor must be made of wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Conductive disturbances induced by RF fields ^{a)}	IEC 61000-4-6	3 V ^{b)} 0,15 MHz - 80 MHz 6 V ^{b)} in ISM bands between 0,15 MHz and 80 MHz 80% AM at 1 KHz	Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended distances, calculated from the equation applicable at the frequency of the transmitter.

a) The following applies:

- All connection cables with the patient must be tested, either individually or grouped together.
- The connection cables with the patient must be tested using a current clamp unless the current clamp is not suitable. If a current clamp is not suitable, an EM clamp must be used.
- In any case, no intentional decoupling device should be used between the injection site and POINT OF CONNECTION TO THE PATIENT.
- The tests can be performed with other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- The tubes that are intentionally filled with conductive liquids and intended to be placed in contact with the PATIENT must be considered connection cables with the patient.

• If an ISM or amateur radio band is not present among the frequency samples, as appropriate, an additional test frequency has to be used in the ISM band or in the amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.

- The ISM bands (industrial, scientific and medical) between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The non-professional bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

b) R. M. S., applied before modulation.

c) The discharges must be applied without connection to an artificial hand and without connection to the PATIENT simulation. The PATIENT simulation can be connected after the test, if necessary, to verify BASIC SAFETY and ESSENTIAL PERFORMANCE.

14.4.3 Parts Accessible to the Input / Output Signals

turbodent touch, with related accessories, is designed to operate in the electromagnetic environment specified below. The purchaser or user of turbodent touch should ensure that it is used in such an environment.

Phenomenon	Essential EMC standard or test method	Immunity test values	Electromagnetic Environment Guidance
Electrostatic discharges (ESD) ^{e)}	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV in air	The floor must be made of wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst ^{b) f)}	IEC 61000-4-4	± 1 kV contact 100 KHz repetition frequency	The quality of the network voltage should be that of a typical commercial or hospital environment.
Pulses Common mode ^{a)}	IEC 61000-4-5	± 2kV	The quality of the network voltage should be that of a typical commercial or hospital environment.
Conductive disturbances induced by RF fields ^{d) g) j) k)}	IEC 61000-4-6	3 V ^{h)} 0,15 MHz - 80 MHz 6 V ^{h)} in ISM bands between 0,15 MHz and 80 MHz ⁱ⁾ 80% AM at 1 KHz ^{c)}	Portable and mobile RF communication devices should not be used near any part of the product, including cables, except when they respect the recommended distances, calculated from the equation applicable at the frequency of the transmitter.

- a) This test applies only to output lines connected directly to the external cables.
- b) SIP/SOPS with maximum cable length less than 3 m are excluded.
- c) The tests can be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- d) The calibration of the current injection terminals must be performed in a system at 150 Ω.
- e) The connectors must be tested in accordance with Section 8. 3. 2 and Table 4 of IEC 61000-4-2:2008. For insulated connector housings, perform the air discharge test on the connector housing and pins using the probe with the rounded tip of the ESD generator, with the exception that only the connector pins that are tested are those that can be reached or touched, under the conditions of INTENDED USE, by the standard probe shown in Figure 6 of the general standard, applied in a bent or straight position.
- f) The capacitive coupling must be used.
- g) If an ISM or amateur radio band is not present among

the frequency samples, as appropriate, an additional test frequency has to be used in the ISM band or in the amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.

- h) R. M. S., applied before modulation.
- i) The ISM bands (industrial, scientific and medical) between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The non-professional bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.
- j) See IEC 61000-4-6:2013, Annex B, for the modified start frequency with respect to cable length and device size.
- k) SIP/SOPS with maximum cable length less than 1 m are excluded.

14.5 Specifications of the Tests for the Immunity of the Accessible Parts of the Casing to the Wireless RF Communications Device

turbodent touch, with related accessories, is designed to operate in an electromagnetic environment in which radiated RF disturbances are under control. The purchaser or operator of turbodent touch can help prevent electromagnetic interferences by guaranteeing a minimum distance between the mobile and portable RF communication devices (transmitters) and turbodent touch , as recommended below, in relation to the maximum output power of the radio communication devices.

Test Freq. (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation	Max power (W)	Distance (m)	Immunity test value (V/m)
385	380 - 390	TETRA 400	Pulse modulation ^{b)} 18 Hz	1,8	0,3	27
450	430 - 470	GMRS 460 FRS 460	FM ^{c)} ± 5 kHz deviation 1 kHz sine	2	0,3	28
710	704 - 787	LTE band 13,17	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9
745						
780						
810	800 - 960	GSM 800/900 TETRA 800 iDEN 820 CDMA 850 Band LTE 5	Pulse modulation ^{b)} 18 Hz	2	0,3	28
870						
930						
1720	1700 - 1990	GSM 1800 CDMA 1900 GSM 1900 DECT LTE Band 1, 3, 4,25 UMTS	Pulse modulation ^{b)} 217 Hz	2	0,3	28
1845						
1970						
2450	2400 - 2570	Bluetooth WLAN 802. 11 b/g/n RFID 2450 Band LTE 7	Pulse modulation ^{b)} 217 Hz	2	0,3	28
5240	5100 - 5800	WLAN 802. 11 a/n	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9
5500						
5785						

- a) For some services, only uplink frequencies are included.
- b) The carrier must be modulated using a square-wave signal with a duty cycle of 50%.
- c) As an alternative to FM modulation, the carrier can be pulse-modulated using an 18 Hz square-wave signal with 50% duty cycle. Although this does not represent actual modulation, it would be worst case.

NOTE: If necessary to reach the immunity test level, the distance between the transmitter antenna turbodent touch can be reduced to 1 m. The test distance of 1 m is allowed by IEC 61000-4-3.

 **WARNING:** Portable RF communication equipment (including peripheral devices such as antenna cables and external antennas) must not be used closer than 30 cm to any part of the turbodent touch device, including the cables specified by the manufacturer. Otherwise, there may be a performance degradation of these appliances.

EN

14.6 Proximity Magnetic Field Immunity in the Frequency Range from 9 kHz to 13,56 MHz

The following table outlines the specifications of the test for IMMUNITY of the CASING to proximity magnetic fields in the frequency range from 9 kHz to 13,56 MHz.

Test frequency	Modulation	Immunity test level (A/m)
30 kHz ^{a)}	CW	8
134,2 kHz	Pulse modulation ^{b)} 2,1 kHz	65 ^{c)}
13,56 MHz	Pulse modulation ^{b)} 50 kHz	7,5 ^{c)}

- a) This test applies only to devices intended for use in the HOME HEALTHCARE ENVIRONMENT.
- b) The carrier must be modulated using a square-wave signal with 50% duty cycle.
- c) r.m.s., before modulation is applied.

15 TROUBLESHOOTING

15.1 Diagnostic Reports and Solutions

turbodent touch is provided with a diagnostic circuit that allows the detection of malfunctions and viewing of their type on the touch surface by means of a symbol. Users, by using the following table, are guided to identifying and the possible resolution of the malfunction detected.

EN

Symbol on touch surface	Possible cause	Solution
	Incorrect start-up procedure: the device was started with the foot pedal pressed.	Check that the foot pedal is not pressed down. If the problem persists, disconnect the foot pedal and potentially contact a MECTRON Authorised Service Centre.
	Incorrect power-on procedure: the device was operated during the calibration phase.	Switch off and wait 5 seconds before switching on the device again. Do not operate the foot pedal during the calibration phase. If the warning persists, contact a MECTRON Authorised Service Centre.
	The self-diagnostic functions detected a fault or malfunction.	Switch off and wait 5 seconds before switching on the device again. If the warning persists, contact a MECTRON Authorised Service Centre.
	Mains faults, electrostatic discharges.	Turn off and wait 5 seconds before switching the device on again. If the warning persists, contact a MECTRON Authorised Service Centre.
	A pressurised powder container was opened without running the "refill" cycle.	Before opening one of the powder containers, the "refill" cycle must be run (See Chapter 5.5 on page 24).
	The device fails to regulate the air pressure correctly.	Check for compressed air in the system and check that the connector is correctly connected. Switch off and wait 5 seconds before switching on the device again. Do not operate the foot pedal during the calibration phase. If the warning persists, contact a MECTRON Authorised Service Centre.
	Condensate drainage was attempted without activating one of the two "no powder" functions.	Switch the device off and on again, set the "No Powder" function (by selecting   or ) and attempt to drain the condensate again.

NOTE: For diagnostic warnings not included in this list, contact the technical service desk.

15.2 Quick Troubleshooting

Problem	Possible Cause	Solution
The device does not start after the switch is set to "I".	The terminal of the electrical power cable has not been properly inserted in the rear socket of the device.	Check that power supply cable is firmly connected.
	The power supply cable is faulty.	Check that the supply socket is working. Replace the power supply cable.
	The fuses are out of order.	Replace the fuses (See <i>Chapter 15.3 on page 74</i>).
The device is on but not working. No errors are reported on the touch surface.	The foot pedal plug is not properly inserted in the device socket.	Properly insert the foot pedal plug into the connector at the back of the device (see <i>Chapter 4.3 on page 13</i>).
	The foot pedal does not work properly.	Contact an Authorised MECTRON Service Centre.
The device is on but not working. The following symbol appears on the touch surface: 	See <i>Chapter 15.1 on page 70</i> for the possible cause	See <i>Chapter 15.1 on page 70</i> for details of the necessary action.
No liquid flows out from the AIR-POLISHING handpiece during operation.	Device not connected to the water circuit.	Check the connection to the water circuit (See <i>Chapter 4.3 on page 13</i>).
	The AIR-POLISHING handpiece is clogged.	Disconnect the AIR-POLISHING handpiece and clear the water passage by blowing compressed air through it. See <i>Chapter 12 on page 50</i> See <i>Chapter 6.4 on page 32</i> If the problem persists, replace the AIR-POLISHING handpiece with a new one.
	The water tap on the device is closed.	Adjust the water flow using the dedicated knob See <i>Chapter 5.2 on page 17</i> .
	The foot pedal connector is not properly inserted in the device socket.	Properly insert the foot pedal connector at the back of the device (see <i>Chapter 4.3 on page 13</i>).

Problem	Possible Cause	Solution
No powder flows out from the AIR-POLISHING handpiece during operation.	The AIR-POLISHING handpiece is clogged due to an excessive moisture content in the powder or insufficient cleaning/ maintenance.	Remove the powder from the container and clean it using a dry cloth. Disconnect the AIR-POLISHING handpiece and clear the water passage by blowing compressed air through it. See <i>Chapter 12 on page 50</i> See <i>Chapter 6.4 on page 32</i> If the problem persists, replace the AIR-POLISHING handpiece with a new one.
	The amount of powder in the container exceeds the maximum level.	Remove the powder from the container and clean it using a dry cloth. Restore the correct level of powder in the container (See <i>Chapter 4 on page 12</i>).
	Unsuitable powder.	Use the correct powder to ensure proper functioning of the device.
Powder leak through the powder container cap.	Powder residue in the thread.	Clean the powder container thread (See <i>Chapter 12.7 on page 58</i>).
Poor cleaning performance.	Insufficient pressure in the air supply circuit.	Check air supply circuit pressure (4,5 - 8,0 bar).
	Insufficient or excessive powder level in the container.	Restore the correct level of powder in the container See <i>Chapter 4.3 on page 13</i> .
	Unsuitable powder.	Use the correct powder to ensure proper functioning of the device.
	The AIR-POLISHING handpiece is clogged due to an excessive moisture content in the powder or insufficient cleaning/ maintenance.	Remove the powder from the container and clean it using a dry cloth. Disconnect the AIR-POLISHING handpiece and clear the water passage by blowing compressed air through it. See <i>Chapter 12 on page 50</i> See <i>Chapter 6.4 on page 32</i> If the problem persists, replace the AIR-POLISHING handpiece with a new one.

Problem	Possible Cause	Solution
One of the powder container caps does not unscrew.	The device is on and the powder container is pressurised.	Run the "refill" cycle before opening one of the powder containers (See <i>Chapter 5.5 on page 24</i>).
	The "refill" cycle has been performed but the powder containers have remained pressurised because the AIR-POLISHING handpiece is clogged.	Read the section on cleaning the AIR-POLISHING handpiece (See <i>Chapter 6.4 on page 32</i>).
	The "refill" cycle has been performed but the powder containers have remained pressurised because the air-polishing cord is clogged.	Contact an Authorised MECTRON Service Centre.

15.3 Replacing the Fuses


⚠ WARNING: Switch the device off.
Always switch the device off using the main switch and disconnect it from the

power socket before performing the next procedure.

EN

Use a flat tool, if necessary, to open the fuse box located below the power socket;


1

Extract the fuse holder compartment;


⚠ WARNING: Replace the fuses in respect of the characteristics indicated in Chapter 14 on page 61.

2

Reinsert the compartment in its place.

3

15.4 Shipping to an Authorised MECTRON Service Centre

If technical assistance is required on the device, contact an Authorised MECTRON Service Centre or your Dealer. Do not try to repair or modify the device and its components/accessories.

Clean and sterilise all parts to be sent to a MECTRON Authorised Service Centre following the instructions in *Chapters 6, 7, 8, 9*. Leave the sterilised parts in the bag which certifies the sterilisation process.

The cleaning and sterilisation demands comply with the mandatory requirements for workplace health and safety protection pursuant to Legislative Decree 81/08 and Italian laws.

If the client fails to comply with the requirements, MECTRON reserves the right to charge them the cost of cleaning and sterilisation or to refuse goods that arrive in unsuitable conditions, returning them at the client's own expense for proper cleaning and sterilisation.

The device and its components/accessories must be returned properly packaged with all components and a slip stating:

- Details of owner and contact number;
- Product name;
- Serial number and/or batch number;
- Reason for return / description of failure;
- Photocopy of purchase bill or invoice for the device.

CAUTION: Packaging

Pack the device and/or its accessories in the original packaging to avoid damage during transport.

Once the material is received at the Authorised MECTRON Service Centre, the qualified technical personnel will evaluate the problem. The repair will be made only upon acceptance by the end client. For further details contact your nearest Authorised MECTRON Service Centre or your dealer.

Unauthorised repairs can damage the system and void the warranty, releasing MECTRON of all liability for direct or indirect damages to persons or property

16 WARRANTY

All MECTRON devices undergo a thorough final inspection to ensure full functionality before being marketed.

MECTRON provides a warranty for turbodent touch and related accessories, purchased new from a MECTRON dealer or importer, which covers defects in material and workmanship for a period of:

- 2 (TWO) YEARS for the device from the date of purchase;
- 1 (ONE) YEAR for the handpiece from the date of purchase.

The other components/accessories are not covered by the warranty.

During the warranty period, MECTRON undertakes to repair (or, at its discretion, to replace) the parts of products free of charge, which, according to its judgement, are proven to be defective.

Full replacement of MECTRON products is not covered by the warranty.

The manufacturer's warranty and device approval are not valid in the following cases:

- The device is not used in accordance with the intended use.
- The device, with compatible accessories, is not used in accordance with all the instructions and requirements outlined in this manual.
- The electrical system of the facilities where the device is used does not comply with the electrical code compliance standards in force and the relative electrical safety precautions.
- Assembly operations, extensions, adjustments, updates and repairs are carried out by personnel not authorised by MECTRON.

- The environmental conditions for preservation and storage of the device and its accessories do not comply with the requirements indicated in *Chapter 14 on page 61*.
- Use of non-original MECTRON components/accessories and replacement parts that may compromise the correct operation of the device and cause injury to the patient.
- Accidental breakage during transport.
- Damage due to incorrect use or carelessness, or due to connection to a voltage other than that envisaged.
- Expired warranty.

The expected service life of the device is minimum 5 years.

The service life / duration does not establish a limit of use; the service life of the device defines the period of time after installation and/or commissioning, during which the original performance, or in any case performance suited to the intended use, is guaranteed without there being any degradation such as to compromise functionality and reliability.

The service life is a minimum qualitative objective of the design, therefore, individual parts or components may guarantee superior performance and reliability with respect to that declared by the manufacturer.

The service life assumes compliance with the maintenance schedules set out in this manual, does not include components normally subject to "wear", and is not linked to the warranty period: the service life does not establish any implicit or explicit extension of the warranty period.

CAUTION

The warranty starts from the date of purchase of the device, which evidence is given by the delivery note/purchase invoice issued by the Dealer / Importer or, in case of device with activation code, from the date of activation of the same.

In order to benefit from the warranty service, the client must return the device, at their own expense, to the MECTRON dealer/importer from which they purchased the product.

The device must be returned in its original packaging, accompanied by all components/accessories and a slip stating:

- The data of the owner and telephone number;
- The data of the Dealer / Importer;
- Photocopy of the purchase bill/invoice for the device issued to the owner, stating the purchase date, name of the device and serial number;
- Description of the failure.

The transport and the damage caused by transport are not covered by the warranty.

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MECTRON S.p.A.